ECA MEDICAL INSTRUMENTS.

Trusted Surgical Solutions

ECA Medical Instruments
Instrument Use and Care

Manufactured by:

ECA Medical Instruments

1107 Tourmaline Drive, Newbury Park, CA 91320

805-376-2509

MDSS GmbH

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30175 Hannover, Germany

Distributed by:

stryker[®]

Spine

0459

Single Procedure Instrumentation

The CE mark is only valid if it is also printed on the product label.

Warning

ECA Torque Limiting Devices should only be used by qualified personnel trained in the use of surgical instruments and the relevant surgical procedures. The surgeon should thoroughly understand all aspects of the surgical procedures and the limitations of the instrumentation. It is strictly prohibited to carry out any modifications on any ECA Medical Instruments manufactured device. ECA Medical Instruments disclaims any and all liability resulting from use of such modified devices.

- **DO NOT REUSE or REPROCESS. Device is a Single Procedure Torque Limiting Instrument.** A second PROCEDURE use is prohibited, and may reduce the performance of the instrument and puts the patient and the surgeon at risk. The instrument must be properly disposed of after each PROCEDURE.

General guidelines:

- No more than 60 actuations of the torque limit set point. Exceeding this number of actuations may void torque setting.
- Do not use if there are any indications of damage, excessive wear or performance failure.

IMPORTANT: Anyone of the above conditions is an absolute indication to dispose of the instrument and to discontinue use.

Precautions

All health care workers should apply Universal Precautions when handling sharp instruments and instruments using torque limiting set points of any force.

Instruction of use

Any Torque Limiting Instrument **MUST** be restricted to use by a physician or qualified health care professional and used in accordance with manufacturer's instruction. The surgeon or qualified health care professional is responsible for the

prescribed use of the Torque Limiting Instrument. The Torque Limiting function may drift out of calibration if used beyond the prescribed useful life. Improper use may void the calibration and warranty. Do not use the Torque Limiting device beyond manufacturers' recommendations. Before using device in situ, actuate the device one to two times to confirm audible actuation. The Torque Limiting Instrument should not be used for any other function than the intended manufacturers' recommendations and design. Exceeding, disregarding and NOT following these parameters may cause the device to break and or not perform correctly. Efficient screw tightening requires a delicate calibration.

Be careful not to allow the device to be mishandled. If the torque increases over the maximum set limit, the coupling will slip to protect against risk of breakage. At this point, an audible Torque Limiting actuation has occurred.

Sterility Information

- This instrument is supplied by the manufacturer in a sterile condition; it has been sterilized by Gamma radiation at a dosage between 25-40kGy (noted as STERILE-R on the label).
- For shelf life expiration, refer to date printed on the product label. Meticulously inspect the package before use. Do not use the instrument if the package has been opened, is torn or punctured, or if the seal has been broken. Use the device as soon as package has been opened; otherwise the sterility has been compromised.

Storage and Handling

- Devices should be handled with care at all times. Storage zones for surgical instruments should be away from areas of humidity and must be out of contact with UV rays and sources of electro-magnetic radiation.
- Do not use the sterile surgical instrument after the expiration date indicated on the package and label. After this date, the sterility of the instruments can no longer be assured.

Sterilization

- Steam sterilization, including flash sterilization, of this instrument, is NOT recommended.

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© 2013, ECA Medical Instruments. All rights reserved. ECA Medical Instruments is a trademark or registered trademark of ECA Medical Instruments. Rev A (03/13)