


Suggestions for skills assessment with the **LIFEPAK® 35** Monitor/Defibrillator

Automated external defibrillation (AED)

Name: _____ Unit: _____

Reviewer: _____ Date: _____

This skills checklist assesses a trained user's ability to perform defibrillation with the LIFEPAK 35 monitor/defibrillator. This limited checklist does not cover all information and skills required to operate the device safely and effectively. This checklist is designed to be completed after receiving or viewing in-service training by a qualified instructor. Refer to the LIFEPAK 35 monitor/defibrillator Operating Instructions for complete directions for use, indications, contraindications, warnings, cautions, and potential adverse events. CAPITAL LETTERS indicate labels, screen messages, and voice prompts. "Push" is used when referring to keypad buttons, and "Press" is used when referring to softkeys on the touchscreen indicated by italics

Assessment criteria	Complete	Incomplete	Comments
1. verifies the patient is unconscious, pulseless, and not breathing normally.			
2. Pushes the POWER button.			
3. Describes skin prep for the pads according to local protocol or Operating Instructions.			
4. Connects the therapy electrode pads to the therapy cable. Confirms therapy cable is connected to the device.			
5. Applies the therapy electrode pads to the patient's chest in the anterior-lateral position			If the defibrillator turns on in Manual mode, pushes the ANALYZE button to change to AED mode.
6. Stops all motion, including CPR as needed. 7. Pushes the ANALYZE button to initiate the analysis. 8. Confirms PEDIATRIC AED or ADULT AED mode.			PEDIATRIC setting is for patients weighing less than 25 kg (55 lbs) or who are less than 8 years old. The pediatric energy setting is 50J, 75J, and 90J.
9. Continue to follow the screen messages and voice prompts provided by the AED.			CPR timer is in real or elapsed time.
10. SHOCK ADVISED: <ul style="list-style-type: none"> Provides a verbal "Clear" and confirms all clear Pushes  button to deliver energy to the patient. Follows the START CPR prompt after shock delivery. 			Pressing SPEED DIAL cancels the charge or charge with disarm after 60 seconds.
11. Acknowledges CPR prompt and timer will appear and continue per duration specified per medical direction local protocols or default settings in device.			
12. NO SHOCK ADVISED: <ul style="list-style-type: none"> Begins CPR per START CPR voice prompts, follows screen messages, and metronome. Acknowledges when a No Shock Advised decision follows a shock, the energy level does not increase for the next shock if needed. Follows metronome sounds for compression and ventilation rate and acknowledges the CPR timer shows the remaining time. 			The metronome automatically provides compression tocks and ventilation prompts for adults at defaults of either 30/2 or 10/1 and for pediatrics at 15/2 or 10/1 accessed through the C:V tab at the bottom of the screen. To silence or unsilenced the metronome, pushes the MUTE button note.

Assessment criteria	Complete	Incomplete	Comments
<p>13. cprINSIGHT enabled:</p> <ul style="list-style-type: none"> Acknowledges cprINSIGHT technology analyzes the heart rhythm during chest compressions after the first initial analysis and may give three decisions: <ul style="list-style-type: none"> NO SHOCK ADVISED message is displayed and CPR is continued. SHOCK ADVISED message is displayed, and the device precharges during the last part of the CPR time. Once fully charged clears everyone away from the patient and pushes the SHOCK button Inconclusive in which case the user is instructed to push the ANALYZE button at the end of the CPR time and instructed to stand clear for analysis. Acknowledges the device cannot be moved or used in a moving vehicle with cprINSIGHT enabled. <p>14. Acknowledges cprINSIGHT patients with internal pacemaker the device will revert back to it's traditional algorithm the device will be instructed to push the ANALYZE button at the end of each CPR time.</p>			<p>cprINSIGHT analysis technology enables the defibrillator to analyze the patient's heart rhythm while CPR is being performed. When enabled, the user must stand clear during the first analysis.</p> <p>With an inconclusive decision, the defibrillator reverts back to traditional analysis, in which case the user is instructed to clear for analysis.</p>
15. Describes troubleshooting for: Motion Detected			
<p>16. Switching from AED to Manual Mode:</p> <ul style="list-style-type: none"> Presses the MANUAL DEFIB button OR Pushes the CHARGE button 			Access to manual mode may be restricted depending on the setup, and confirmation or passcode may be required.
<p>17. Switching back to AED Mode:</p> <ul style="list-style-type: none"> Presses the ANALYZE button OR Presses RESUME AED button 			
Optional: Preshock CPR time: When the preshock CPR time is set for 15 seconds or more, you are prompted to start CPR immediately after a shockable rhythm is detected before the shock is delivered. The device will silently charge in preparation for the shock.			The shock button is disabled during the preshock CPR interval to avoid accidental shock delivery while the defibrillator is charging, and the responder is performing CPR.

Notes: _____



Scan here for additional resources

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate Stryker's product offerings. A healthcare professional must always refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: cprINSIGHT, LIFEPAK, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.



Manufactured by:
Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
stryker.com/emergencycare