


# CHARLOTTE<sup>®</sup>

3.0 and 4.3 Multi-Use Compression Screws

SURGICAL TECHNIQUE



---

# CHARLOTTE®

## 3.0 and 4.3 Multi-Use Compression Screws

### SURGICAL TECHNIQUE

---

#### Surgical Technique as described by:

Robert Anderson, MD  
Bruce Cohen, MD  
W. Hodges Davis, MD

---



# Contents

Chapter 1	<b>6</b>	<b>Introduction</b>
		Product Information
		Design Features
Chapter 2	<b>7</b>	<b>Intended Use</b>
		Indications and Contraindications
Chapter 3	<b>8</b>	<b>Surgical Technique</b>
Appendix	<b>10</b>	<b>Ordering Information</b>

Proper surgical procedures and techniques are the responsibility of the medical professional. The following guidelines are furnished for information purposes only. Each surgeon must evaluate the appropriateness of the procedures based on his or her personal medical training and experience. Prior to use of the system, the surgeon should refer to the product package insert for complete warnings, precautions, indications, contraindications and adverse effects. Package inserts are also available by contacting Wright Medical Technology, Inc. Contact information can be found on the back of this surgical technique and the package insert is available on the website listed.

Please contact your local Wright representative for product availability.

The CHARLOTTE® Multi-Use Compression (MUC) Screw is a cannulated, headless design that is appropriate for certain osteotomies and fusions of the forefoot and midfoot. The 3.0mm screws are particularly suited to fixation of 1st metatarsal osteotomies and for small joint periarticular fixation. The 4.3mm screws are excellent for talonavicular fusions, midfoot fusions, and Hallux interphalangeal fusions.

## DESIGN FEATURES

The CHARLOTTE® MUC Screw obtains compression between two bony fragments via a differing thread pitch at the leading and trailing ends of the screw. The screw distinguishes itself from previous designs with its double-lead helical thread pattern at the leading end of the screw, which doubles the thread engagement in the distal bone fragment. In addition, its spiral-fluted cutting design allows the screw to be self-tapping and self-drilling in most bone.

## SURGICAL GOALS

- Provide maximum compression across the fusion site of two adjacent bones.
- Obtain maximum thread engagement in the distal fragment for maximum compression.
- Ensure that the head of the screw is completely countersunk, so profile-related issues may be avoided.

## SYSTEM BASICS

- All CHARLOTTE® MUC Screw implant components are manufactured from surgical grade stainless steel.
- 3.0mm diameter screw comes in 10-34mm lengths in 2mm increments.
- 4.3mm diameter screw comes in 14-50mm lengths in 2mm increments, with additional 55mm and 60mm screws.
- 4.3mm diameter screw from 36-60mm have a “long thread” version, which is threaded over half of its length, and a “short thread” version, which is threaded over one-third of its length.
- 3.0mm screw is cannulated to work over a 1mm single-tip K-Wire, which is included in the set.
- 4.3mm screw is cannulated to work over a 1.6mm single-tip K-Wire, which is included in the set.
- Cannulated Drills for the 3.0 and 4.3mm screws are included for use in hard cortical bone, an oblique approach, or when bicortical fixation is desired.
- Cannulated countersinks for the 3.0 and 4.3mm screws are included for countersinking the screw heads in hard cortical bone. These should also be used if the head will be placed in a thin, delicate section of bone.

## INDICATIONS

The CHARLOTTE® MUC Screw is indicated for fixation of bone fractures or for bone reconstruction. Examples include:

- Mono or Bi-Cortical osteotomies in the foot or hand
- Distal or Proximal metatarsal or metacarpal osteotomies
- Weil osteotomy
- Fusion of the first metatarsophalangeal joint and interphalangeal joint
- Fixation of osteotomies for Hallux Valgus treatment (such as Scarf, Chevron, etc.)
- Akin type osteotomy
- Arthrodesis base first metatarsal cuneiform joint to reposition and stabilize metatarsus varus primus
- Calcaneus/ cuboid arthrodesis
- Talar/ navicular arthrodesis

## CONTRAINDICATIONS

- Infection
- Physiologically or psychologically inadequate patient
- Inadequate skin, bone, or neurovascular status
- rreparable tendon system
- Possibility for conservative treatment
- Growing patients with open epiphyses
- Patients with high levels of activity

Prior to use of the system, the surgeon should refer to the product package insert for complete warnings, precautions, indications, contraindications and adverse effects. Package inserts are also available by contacting the manufacturer. Contact information can be found on the back of this surgical technique and the package insert is available on the website listed."

## TALONAVICULAR FUSION

### EXPOSURE/JOINT PREPARATION

One or two of the 4.3mm CHARLOTTE® MUC Screws will be used for this procedure. Expose the talonavicular joint using a standard medial approach. Distract the joint with a lamina spreader, and sharply debride the articular cartilage to expose bleeding subchondral bone. A powered drill can also be used to further penetrate the subchondral bone to ensure that bleeding bony surfaces are in apposition prior to screw insertion.

### K-WIRE PLACEMENT

Use a powered driver to place a 1.6mm K-Wire (P/N 44112008) through the medial cortex of the navicular. The K-Wire is directed from plantar medial to dorsal lateral through the talonavicular joint and advanced until it contacts (but does not penetrate) the lateral cortex of the talus. | **FIGURE 1** Verify the position of the wire fluoroscopically.




FIGURE 1 |



FIGURE 2 |

### SCREW LENGTH DETERMINATION

Use the CHARLOTTE® Cannulated Depth Gauge (P/N 44112002) over the K-Wire to measure the correct length for the 4.3mm screw. | **FIGURE 2** Long or short thread length is most easily determined radiographically. If more than half of the K-Wire length is on the far side of the joint (in the talus), choose a CHARLOTTE® Long Thread Screw. If not, choose a CHARLOTTE® Short Thread Screw.

## HEAD PREPARATION


The CHARLOTTE® MUC Screw is designed to be self-tapping and self-drilling. However, in extremely dense cortical bone, it may be difficult to install the screw. Included are: 2.0mm drill (P/N 44112004) and 3.0mm head drill (P/N 44112012) for the 3.0mm screws; 3.0mm drill (P/N 44112003); and 4.3mm head drill (P/N 44112011) for the 4.3mm screws. Use the drill and then the head drill over the K-Wire to prepare the bone.

## SCREW PLACEMENT

Load the 3mm CHARLOTTE® Cannulated Hex Driver (P/N 44112007) into the CHARLOTTE® Cannulated AO Driver Handle (P/N 44112009). Use the driver to advance the chosen 4.3mm screw over the K-Wire. | **FIGURE 3** Advance the screw until the head is completely countersunk within the bone. | **FIGURE 4** Depending on the stability of the first screw and patient related factors (obesity, post-operative compliance issues), a second screw may be used for additional fixation. Surgical closure is then performed in the normal fashion.



| FIGURE 3



| FIGURE 4



# Ordering Information

## CHARLOTTE®

PART NUMBER	DESCRIPTION
44110014	4.3MM X 14MM SHORT
44110015	4.3MM X 16MM SHORT
44110016	4.3MM X 18MM SHORT
44110017	4.3MM X 20MM SHORT
44110018	4.3MM X 22MM SHORT
44110019	4.3MM X 24MM SHORT
44110020	4.3MM X 26MM SHORT
44110021	4.3MM X 28MM SHORT
44110022	4.3MM X 30MM SHORT
44110023	4.3MM X 32MM SHORT
44110024	4.3MM X 34MM SHORT
44110025	4.3MM X 36MM SHORT
44110026	4.3MM X 38MM SHORT
44110027	4.3MM X 40MM SHORT
44110028	4.3MM X 42MM SHORT
44110029	4.3MM X 44MM SHORT
44110030	4.3MM X 46MM SHORT
44110031	4.3MM X 48MM SHORT
44110032	4.3MM X 50MM SHORT
44110033	4.3MM X 36MM LONG
44110034	4.3MM X 38MM LONG
44110035	4.3MM X 40MM LONG
44110036	4.3MM X 42MM LONG
44110037	4.3MM X 44MM LONG
44110038	4.3MM X 46MM LONG
44110039	4.3MM X 48MM LONG
44110040	4.3MM X 50MM LONG
44110041	4.3MM X 55MM LONG
44110042	4.3MM X 60MM LONG
44110001	3.0MM X 10MM
44110002	3.0MM X 12MM
44110003	3.0MM X 14MM
44110004	3.0MM X 16MM
44110005	3.0MM X 18MM
44110006	3.0MM X 20MM
44110007	3.0MM X 22MM
44110008	3.0MM X 24MM
44110009	3.0MM X 26MM
44110010	3.0MM X 28MM
44110011	3.0MM X 30MM
44110012	3.0MM X 32MM
44110013	3.0MM X 34MM
44112008	SINGLE K-WIRE, 1.6 X 150MM
44112000	SINGLE K-WIRE, 1.0 X 150MM



INSTRUMENTS	
PART NUMBER	DESCRIPTION
44112007	3MM CANN HEX DRIVER
44112011	CANN. HEAD DRILL FOR 4.3MM SCREW
44112012	CANN. HEAD DRILL FOR 3.0MM SCREW
44112004	CANN. DRILL, 3.0MM SCREW
44112003	CANN. DRILL, 4.3MM SCREW
44112001	2MM CANN HEX DRIVER
44112009	AO DRIVER HANDLE
41112017	AO QUICK CONNECT, CANNULATED
44112002	DEPTH GAUGE (CANN.)

SURGICAL TRAY	
PART NUMBER	DESCRIPTION
44112006	SURGICAL TRAY

# Notes

# Notes



**Wright Medical Technology, Inc.**  
1023 Cherry Road  
Memphis, TN 38117  
800 238 7117  
901 867 9971  
[www.wmt.com](http://www.wmt.com)

**Wright Medical EMEA**  
Atlas Arena, Australia Building  
Hoogoorddreef 7  
1101 BA Amsterdam  
the Netherlands  
011 31 20 565 9060

™Trademarks and \*Registered marks of Wright Medical Technology, Inc.  
©2014 Wright Medical Technology, Inc. All Rights Reserved.

S0028-105