

stryker

Omega

Knotless Anchor System



Eyelet



3.9mm




4.75mm



6.5mm


The **Omega Knotless Anchor System** features a decoupled eyelet and screw construct which offers intraoperative and key advantages. The eyelet is independent of the screw and has a proprietary design to reposition during suture tensioning and customise its final position based on bone density. The decoupled eyelet also allows for an intraoperative change in screw size, or sizing up, without the removal of the eyelet or rethreading of the sutures.



Self-punching all-PEEK eyelet


The eyelet is all-PEEK and is preloaded on a self-punching inserter which is designed to eliminate the need to drill.

(Drills, awls and taps are available)




Eyelet **repositioning**

As the suture limbs are tensioned, the eyelet is designed to move and customise its final position based on the patient's specific bone type.




Soft bone

In soft bone, the eyelet is designed to flip and move proximally until the eyelet barb engages with bone and reaches a secure position.



Hard bone

In hard bone, due to the bone's greater density, the eyelet is designed to engage more distally as sutures are tensioned.



Unobstructed complete cannulation

Omega screws have been designed to offer unobstructed complete cannulation due to the decoupled eyelet and screw construct and the eyelet's ability to reposition itself during deployment. Anchor cannulation may facilitate the exchange of blood and other autologous cells.

Guiding suture path

Sutures from the previously inserted eyelet will create a suture path and facilitate the insertion point for the screw.

Ease of sizing up

The design of the system allows for the surgeon to intraoperatively size up to a 6.5mm screw without removing the eyelet or rethreading the sutures.



Double-Double

The Omega system offers the ability to intraoperatively reload disposable instrumentation with a second set of implants.



Implants

Part number	Description
3910-500-391	Omega, 3.9mm PEEK Knotless Anchor System, Single (Box of 5)
3910-500-392	Omega, 3.9mm PEEK Knotless Anchor System, Double-Double (Box of 5)
3910-500-471	Omega, 4.75mm PEEK Knotless Anchor System, Single (Box of 5)
3910-500-472	Omega, 4.75mm PEEK Knotless Anchor System, Double-Double (Box of 5)
3910-500-652	Omega, 6.5mm PEEK Knotless Anchor + Eyelet (Box of 5)

Reusable instrumentation

Part number	Description
3910-500-393	Omega, 3.9mm Omega Drill
3910-500-394	Omega, 3.9mm Omega Tap
3910-500-395	Omega, 3.9mm Omega Awl
3910-500-473	Omega, 4.75mm Omega Drill
3910-500-474	Omega, 4.75mm Omega Tap
3910-500-475	Omega, 4.75mm Omega Awl
234-020-158	Multi Purpose T-Handle
234-020-117	Ratcheting Driver Handle

A second eyelet and screw are provided in the Double-Double configuration for two anchor cases.



Sports Medicine

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Omega and Stryker. All other trademarks are trademarks of their respective owners or holders.

This material is not intended for distribution outside the EU and EFTA.