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Instructions for Use Cables & Adapters

Product Name: Cables, adapters
Model/ Specification: See label
Manufacturer: APK Technology Co., Ltd.
Manufacturer (Customer Service): APK Technology Co., Ltd.
Manufacturer Address: 6 Floor, Building B2, Industry of Hengfeng,
No. 739, Zhoushi Road, Hezhou, Hangcheng Street, Bao'an
District, Shenzhen, China
Tel: +86-755-27325581
Website: www.apk-technology.com
Email: info@apk-technology.com

※ Environment Conditions

Operation temperature: +5°C ~ 40°C
Relative humidity(RH): 10% ~ 95%(No condensation)
Atmospheric pressure: 70kPa ~ 106kPa

※ Instructions for Use

1. Cable: Connect the cable to the monitor, and connect the other end to the patient end of the lead wire, SpO2 sensor, temperature probe, and IBP transducer;
2. Adapter: Connect the adapter to the lead wire or sensor, and connect the other end to the trunk cable or device;
3. Please refer to the operation manuals of related products such as lead wires, SpO2 sensors, temperature probes, and monitors for details.

※ Intended Use

1. The cables/ adapters are intended to be used with ECG, EKG, SpO2 and IBP monitoring devices etc. by connecting electrodes or sensors on the patient to a trunk cable or directly to the monitoring device for overall monitoring and/ or diagnostic evaluation by a healthcare professional. Please read and follow all the instructions that come with the device.

2. Adapter reference

PHILIPS


Disposable leadwires

Reusable leadwires

Adapters

※ **Cautions**

1. Carefully place cable/ adapter to reduce the possibility of patient entanglement or strangulation.
2. Hold the cable/ adapter (not the sensing element) while wiping the product.
3. Don't immerse or contact with liquid/ water.
4. Please hold the connector part when unplug the cable/ adapter out from the lead wire/ sensor/ monitor.
5. Please check the cable/ adapter connection status at least every 4 hours.
6. Please keep the plug dry. Moisture in plug may cause failure.
7. Please keep the metal parts clean to avoid incorrect measurement.
8. Don't pull or circle the cable with pull force $\geq 50N$, or the cable will be broken or mangled.
9. To ensure that the cable/ adapter is working properly, it can be tested with a matching monitor. Once the open circuit, short circuit, incorrect waveform, no waveform or no data, inaccurate data occurs, the cable/ adapter must have been broken and the user can't use it continuously. Any damage of the cable/ adapter shall be dealt with under relevant regulations of the hospital.
10. Please strictly follow the instructions for use of cables and adapters.

※ **Cleaning and Disinfection**

1. Please unplug the cable/ adapter from the monitor before cleaning or disinfection.
2. Clean or disinfect the cable/ adapter before attaching to a new patient.
3. Clean or disinfecting the product periodically.
4. Wipe the cable/ adapter with a soft cloth saturated with mild soap and water to clean, wipe it with alcohol to disinfect, then wipe it dry and remove all residue.
5. We recommend using 75% medical alcohol as the disinfection solution to disinfect the cable/ adapter and patient contact surface.

※ **Warning**

1. Never immerse cables/ adapters in any liquid.
2. When disconnecting, pull the plug itself, not the cable.
3. Be sure that the leadset/ sensor/ trunk cable/ device connector is properly and fully inserted into the adapter connector. All connections must be secure.

※ **Transportation and Storage**

1. Transportation: In order to prevent severe vibration, rain and snow during transportation, transport requirements shall be in accordance with the order

contract.

2. Packaged products shall be stored in a well-ventilated and non-corrosive room environment with temperature $-20 \sim +55^{\circ}C$ and relative humidity $\leq 95\%$ (No condensation).

※ **Icons Explanation**

Icon	Definition	Icon	Definition
	Instructions For Use		Don't Throw It Away
	Warning		Non-Sterile Packaging
	Batch Code		Must Be Purchased And Used Under Physician's Supervision
	Catalogue Number		Latex Free
	Manufacturer		No Immersing
	BF Type		Do Not Autoclave
	European Authorized Representative		Water Proof Level
	Medical Device		CE Mark
	Manufacture Date		Unique Identification Of medical Devices
	Quantity per pouch		

IFU Filed Date: 2021-09-13

If you have any questions, please contact APK's After-Sales Service Team or your local APK's representative.

Ruiyi Healthcare GmbH
Neuer Zollhof 3, 40221 Dusseldorf, Germany