IntelliVue Reusable ECG Lead Sets and Trunk Cables

Instructions for Use

3 ECG		4 ECG		5 ECG		5 (3+2) ECG		6 (4+2) ECG		Trunk Cable	
AAMI	M1671A, M1675A*, M1672A, M1678A*	AAMI	M1532A, M1533A,	AAMI	M1976A, M1979A*, M1978A,	AAMI	M1968A, M1973A*, M1971A, M1974A*	AAMI	M1680A, M1681A		
IEC		IEC		IEC		IEC		IEC			
AAMI	M1673A M1674A			AAMI	M1602A						
IEC				IEC							
AAMI	M1624A M1626A										
IEC											

* Products intended for use in Operating Room (OR) environments

- M1663A (5+5)
- M1665A (6+4)
- M1667A (6)
- M1668A (5)
- M1669A (3)
- M1949A (5+5)
- 989803170171 (3)*
- 989803170181 (5)*

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










EN - IntelliVue Reusable ECG Lead Sets and Trunk Cables Instructions for Use

ASSOCIATED DOCUMENT - Refer to the Addendum that accompanies this instructions for Use (IFU) for validated cleaning agents and procedures (ECG Lead Sets, Cables and Adapter Cables Care, Cleaning and Disinfection). **INTENDED PURPOSE** - The ECG lead set and trunk cable is an accessory device to be used together with vital parameter monitoring device - as the main device intended to be used for monitoring and recording of, and to generate alarms for, multiple physiological parameters of adults, pediatrics, and neonates. **Intended Use** - Philips IntelliVue reusable ECG lead sets and trunk cables are limited by the indications for use of the connected monitoring and diagnostic equipment, and are for use only by healthcare professionals. They are designed for multiple patient use and indicated for monitoring cardiac signals for both diagnostic and monitoring purposes in adult, pediatric, and neonatal patients. Philips operating room (OR) ECG trunk cables are indicated for use to monitor patient ECG in an electrosurgery (ESU) environment when used in combination with a compatible ECG patient lead set on adult, pediatric, infant, and neonatal ECG applications. **Intended Users** - Healthcare/Medical professionals only. **Indications for Use** - Philips IntelliVue reusable ECG lead sets / trunk cables and Philips OR ECG trunk cables are indicated for continuous monitoring of cardiac rhythm for both diagnostic and monitoring purposes. These devices are limited by the indications for use of the connected monitoring and diagnostic equipment in healthcare facilities. These devices are intended to interact with patient intact skin only. **Clinical Benefit** - ECG trunk cables and lead sets are intended to be used with and support the function of the connected monitor. Beyond the indirect clinical benefits of supporting a diagnosis to be made through the intended use of the connected monitor, no additional clinical benefits can be claimed. **Contraindications** - There are no known contraindications. **REACH CHEMICAL MANAGEMENT** - Visit www.philips.com/REACH for information on REACH, Registration, Evaluation, Authorization and Restriction of Chemicals substances that may be contained in some Philips products. **Product Descriptions** - Tables on page 1 define physical characteristics and applications for all ECG lead set and trunk cable products included in this IFU. Column headings within those tables are defined below.

AAMI	Check mark ✓ indicates that the lead set is AAMI color coded.	IEC	Check mark ✓ indicates that the lead set is IEC color coded.
	Check mark ✓ identifies lead sets that are appropriate for chest application.		Check mark ✓ identifies lead sets that are appropriate for arm/leg application.
3 ECG	Indicates the total number of ECG leads in each lead set, and the type of electrode connector at the end of each lead.		Check mark ✓ identifies lead sets and/or trunk cables that may be used in operating room environments.
4 ECG	= Grabber		Defines the lead set or trunk cable catalogue reference number.
5 ECG	= Snap		Identifies length(s) of lead set leads, or trunk cable.
6 ECG	= Mini Clip		
AAMI + IEC (x) or (x-x)	Defines ECG trunk cables as being both AAMI and IEC compatible, with (x) or (x-x) identifying each trunk cable's lead set connector pin configuration.		

COMPATIBILITY - IntelliVue ECG lead sets and trunk cables can be used with any ECG monitor/defibrillator for which they are listed as accessories in that product's IFU. Ensure that leads are plugged into appropriate trunk cables designed for the same number of leads. Incompatible components may result in performance issues. **OPERATING ROOM ENVIRONMENTS** - Product REF numbers shown with an asterisk * on the front cover of this IFU are intended for use in Operating Room (OR) environments. • M1675A, • M1687A, • M1557A, • M1558A, • M1979A, • M1984A, • M1973A, • M1974A, • 989803170171, • 989803170181.

WARNINGS - Modification of this equipment is not allowed. • These products are intended for use by healthcare professionals only. • Attach all applicable labels included with these products to the cables/connectors. • Always clean patient cables when they are removed from their packaging for first-time use. • Before use, read and understand all warnings listed in this IFU. Also refer to the IFU for the monitor/defibrillator for additional warnings and cautions. • Refer to the monitor/defibrillator IFU for information regarding proper lead/electrode placement that complies with standard AAMI or IEC practices. • Do not use any lead set or cable when visual inspection reveals liquid/lint contaminants in cable connectors, or other visual cable damage. • Ensure that leads are securely connected to the trunk cable, and the trunk cable is properly plugged into the monitor/defibrillator. Otherwise incorrect patient physiological data may result. • Ensure that the patient is properly grounded during Electro-Surgery (ESU) procedures to prevent patient/user injury (i.e. electrical shock). • These cables are not suitable for use in an MRI environment as there is the potential for patient burns. • To avoid patient injury, ensure that cables are carefully positioned to avoid entanglement, choking, and strangulation. • When used to monitor ECG in the OR, ensure that you use accessories that provide appropriate ESU protection (orange OR lead sets, or ICU lead sets paired with Philips TRUNK Cables 989803170171 or 989803170181) to prevent patient injury (i.e. burns). • OR lead sets and/or OR trunk cables cannot be used to measure respiration.

CAUTIONS - Federal (US) law restricts this device to sale by or on the order of a medical practitioner. • Do not use in excessively wet environments or under massive influence of fluids (e.g. rain). • Do not immerse the cable connectors in any fluid. • Do not use alcohol or use ultrasonic cleaners on the cable. • Do not clean cable connector electrical contacts or connectors with bleach. • Expected service life of these products is 18 months of typical clinical use. • Do not disconnect a cable by pulling on the lead wires. • Be sure to visually inspect the cable before attaching it to any patient or monitor/defibrillator. See *Visual Inspection Prior to Use*.

LABEL APPLICATION TO TRUNK CABLE, BEFORE FIRST TIME USE - AAMI or IEC label(s) that match your ECG cable color-codes must be selected from a "sticky-back" sheet supplied with your trunk cable. Your selected label(s) must be applied to the trunk cable connector that plugs into the patient lead set. **Label Selection** - Reference the graphics at the front pages of this IFU that show appropriate label(s) for each trunk cable REF number, and the trunk cable connector where your selected label(s) need to be applied. Select only the AAMI or IEC label(s) from the supplied label sheet that match your ECG cable color-coding and your trunk cable REF number. • **M1669A and 989803170171 Trunk Cables** - Select the adult or neonate label.

• **M1668A and 989803170181 Trunk Cables** - Select the All, (Massor) Licker or EASI label. • **M1663A and M1949A Trunk Cables** - Select one label for the top, and the other label for the bottom of your trunk cable connector. **Label Application:** 1- Peel your selected label(s) from the supplied label sheet. 2- Align the angled-cutout on the selected label(s) over the angled-cutout in the recessed area on the trunk cable connector. 3- Seal and fully-press the label(s) in place. **VISUAL INSPECTION PRIOR TO EACH USE** - Before attaching any cable to a patient or monitor, perform a visual inspection to determine if the cable has reached end-of-life. Check for cracks, blisters, peeling, exposed wires, damaged connectors, peeling labels, and similar wear or damage that could compromise accurate readings or cause patient/user injury (e.g. cuts). Whenever visual inspection reveals a cable is no longer suitable for continued use, follow appropriate product disposal procedures (see *Product Disposal*). **PATIENT APPLICATION, ECG LEADS** - Refer to the compatible patient monitor IFU for information regarding proper ECG lead/electrode placement that complies with standard AAMI or IEC practices. **CABLE CLEANING AND DISINFECTION** - Refer to the Addendum that accompanies this IFU titled *ECG Lead Sets, Cables and Adapter Cables Care, Cleaning and Disinfection* for validated cleaning and disinfecting agents and procedures. To reduce risk of cross-contamination: • Clean and disinfect cables before first-time use and before use for a different patient. • Additionally, clean and disinfect reusable cables that are applied to any patient, if they become visibly soiled.

STERILIZATION - All lead sets and cable products described in this IFU are not intended to be sterilized. **PRODUCT DISPOSAL** - Follow approved medical waste disposal methods as specified by your patient care facility or local regulations. **INCIDENT REPORTING** - Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the European Economic Area (EEA) countries, including Switzerland and Turkey in which the user and/or patient is established. **ENVIRONMENTAL SPECIFICATIONS** - Temperature, humidity and atmospheric pressure ratings identified below must be maintained for all of the lead set and cable products described in this IFU. Otherwise product damage could result. **Temperature:** Operating Temperature Range: -10°C to +55°C (14 to 131°F) at 20% RH. To 95% RH. Storage Temperature Range: -20°C to +70°C (-4 to 158°F) at 15% RH. To 90% RH. **Humidity:** Minimum Operating Humidity at +55°C (131°F): 15% RH. 24hrs; Maximum Operating Humidity at +40°C (104°F): 95% RH. 24hrs. **Atmospheric Pressure:** Operating: 650 hPa - 1080 hPa. Storage: 100 hPa - 1080 hPa. **REORDERING INFORMATION** - The following products are CE marked and available in the European Union: **AAMI Color-Coded ECG Lead Sets** • M1671A (989803145091) 3 LEADSET, GRABBER, AAMI, ICU; • M1675A (989803145131) 3 LEADSET, GRABBER, AAMI, OR; • M1673A (989803145111) 3 LEADSET, SNAP, AAMI, ICU; • M1624A (989803144941) UNSHIELDED 3 LEAD MIMIC/ICP AAMI 0.7M; • M1532A (989803144841) 4 LEAD SET GRABBER AAMI, ICU; • M1968A (989803125841) 5 LEADSET, GRABBER, AAMI, ICU; • M1978A (989803145121) 5 LEADSET, CHEST, AAMI, ICU; • M1979A (989803125811) 5 LEADSET, GRABBER, CHEST, AAMI, OR; • M1973A (989803125861) 5 LEADSET, GRABBER, AAMI, OR; • M1626A (989803144911) 5 LEAD SNAP CHEST AAMI, ICU; • M1664 (989803144991) 5 LEADSET, SNAP, AAMI, ICU; • M1680A (989803145161) 6 LEAD SET GRABBER AAMI, ICU. **IEC Color-Coded ECG Lead Sets** • M1672A (989803145101) 3 LEADSET, GRABBER, IEC, ICU; • M1678A (989803145141) 3 LEADSET, GRABBER, IEC, OR; • M1674A (989803145121) 3 LEADSET, SNAP, IEC, ICU; • M1626A (989803144951) UNSHIELDED 3 LEAD MIMIC/ICP IEC 0.7M; • M1533A (989803144851) 4 LEAD SET GRABBER IEC, ICU; • M1971A (989803125851) 5 LEADSET, GRABBER, IEC, ICU; • M1974A (989803125871) 5 LEADSET, GRABBER, IEC, OR; • M1978A (989803125891) 5 LEADSET, GRABBER, CHEST, IEC, ICU; • M1645A (989803145001) 5 LEAD SET SNAP IEC, ICU; • M1681A (989803145171) 6 LEAD SET GRABBER, IEC, ICU; **ECG TRUNK Cables** • M1669A (989803145071) 3 LEAD ECG TRUNK, AAMI/IEC 2.7M; • M1668A (989803145091) 5 LEAD ECG TRUNK, AAMI/IEC 2.7M; • M1667A (989803145051) 6 LEAD ECG TRUNK, AAMI/IEC 2.7M; • M1663A (989803144791) 10 LEAD ECG TRUNK, AAMI/IEC 2.7M; • M1949A (989803125841) 4 LEAD ECG TRUNK CABLE AAMI, IEC 2.7M; • M1665A (989803145041) 6-4 LEAD ECG TRUNK AAMI/IEC 2.7M; • 989803170171 OR 3-LEAD ECG TRUNK CABLE, AAMI/IEC; • 989803170181 OR 5-LEAD ECG TRUNK CABLE, AAMI/IEC.

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Copyright © 2023
All rights reserved
Böblingen GmbH**