

# Single-Patient-Use Disposable ECG 3, 5, 6 AAMI Lead Sets

## *Instructions for Use*


1m (3.3 ft)


989803173121


1m (3.3 ft)


989803173131


1m (3.3 ft)


989803192141


1m (3.3 ft)

989803197511


0.85m (2.8 ft)

989803173151

# PHILIPS

(EN) Instructions for Use.....	1
(ES) Instrucciones de uso .....	1
(FR) Manuel d'utilisation .....	2
(KO) 사용 설명서.....	3
(ZH-CN) 使用说明.....	4
(ZH-TW) 使用說明.....	4

**EN: CAUTION! Separate leads one at a time.** Before use, carefully separate leads from the ribbon cable one lead at a time. Separating multiple leads at the same time could damage the lead sets.

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## EN - Instructions for Use

**ASSOCIATED DOCUMENTS** - Refer to your patient monitor or defibrillator Instructions for Use (IFU) to determine product compatibility, connectivity and patient application for these ECG lead sets. **INTENDED PURPOSE** - In accordance with European Medical Device Regulation 2017/745, this accessory device enables a vital signs monitoring instrument group to record and generate alarms for multiple physiological parameters when monitoring adults and pediatrics. **Intended Use** - Philips Single-Patient-Use Disposable ECG lead sets are intended for use only by trained healthcare professionals for measurement of a patient's ECG for both diagnostic and monitoring purposes. Use is limited by the indications for use of the connected monitoring or diagnostic equipment. These Philips lead sets are intended for short-term use only (an average patient stay of 5 days). **Indication for Use** - Philips ECG lead sets are indicated for use in the monitoring of cardiac signals for both diagnostic and monitoring purposes. Use is limited by the indications for use of the connected monitoring or diagnostic equipment. **Contraindications** - There are no known contraindications. **PRODUCT DESCRIPTIONS** - These single-patient-use disposable ECG lead sets support the following patient applications: **Limb Applications** • 3-lead limb = 989803173121; • 5-lead limb = 989803173131; • 6-lead limb = 989803197511. **Chest Application** • 5-lead chest = 989803192141. **Telemetry Application** • 5-lead telemetry with shower shield = 989803173151. **Electrode Connectors** - All lead set leads terminate with an AAMI color-coded and labeled grabber-style electrode connector. **Ribbon-Style Cable** - Each individual lead peels apart from a ribbon-style cable (see lead separation Caution inside the front cover of this document and below). **COMPATIBILITY** - These ECG lead sets can be used with any ECG monitor/defibrillator for which they are listed as accessories in that product's IFU. Ensure that lead sets are plugged into appropriate trunk cables designed for the same number of leads. Incompatible components may result in performance issues.

**WARNINGS** • Modification of this equipment is not allowed. • Before use, read and understand all warnings listed in this IFU and in your patient monitor/defibrillator IFU, otherwise patient injury could result. • Refer to the monitor/defibrillator IFU and monitor configuration for information regarding proper lead/electrode placement that complies with standard AAMI practices. • Do not use any lead set when visual inspection reveals liquid/lint contaminants in lead set connectors, or other visual cable damage. • Ensure that leads are securely connected to the trunk cable, and the trunk cable is properly plugged into the monitor/defibrillator; otherwise incorrect patient physiological data may result. • Ensure that the patient is properly grounded during Electro-Surgery (ESU) procedures to prevent patient injury (i.e. electrical shock). • These cables are not suitable for use in an MRI environment as there is the potential for patient burns. • To avoid patient injury, ensure that cables are carefully positioned to avoid entanglement, choking, and strangulation. • If the outer jacket insulator on any lead wire is torn or damaged, causing the inner wire to be exposed, do not use and immediately discard the lead set. • These lead sets are for use on one patient only and may be reapplied only to the same patient; do not use on multiple patients as multi-patient use may result in patient infection or cross-infection. • These lead sets may be used in an operating room environment only if used with a Philips OR trunk cable. • For telemetry applications, be sure to use a telemetry lead set that has the included protective shield over the trunk connector.

**CAUTIONS** • Federal (US) law restricts this device to sale by or on the order of a medical practitioner. • Carefully separate leads from the ribbon cable one lead at a time as shown on the inside front cover of this IFU, otherwise lead set damage could result. • These products are not intended to be disinfected or sterilized. • Do not immerse the cable connectors in any fluid. • Do not use in excessively wet environments or under massive influence of fluids (e.g. rain). • Do not autoclave or use ultrasonic cleaners on the cable. • Do not clean cable connector electrical contacts or connectors with bleach. • Expected service life of these products is 5 days of typical clinical use. • Never disconnect by holding or pulling on cable connector wires.

**VISUAL INSPECTION PRIOR TO USE** - Before attaching any lead or cable to a patient or monitor, perform a visual inspection to verify that the lead or cable is not damaged. Check for cracks, blisters, peeling, exposed wires, damaged connectors and similar wear or damage that could compromise accurate readings or cause patient/user injury (e.g. cuts). If visual inspection reveals lead set damage, follow appropriate product disposal procedures (refer to *Product Disposal*). **PATIENT APPLICATION ECG LEADS** - Refer to the monitor or defibrillator IFU for information regarding proper ECG lead/electrode placement that complies with standard AAMI color code practices. **PRODUCT DISPOSAL** - Discard these lead sets after use on one patient. Follow approved medical waste disposal methods as specified by your patient care facility or local regulations.

**INCIDENT REPORTING** - Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. **ENVIRONMENTAL SPECIFICATIONS** - The temperature, humidity and atmospheric pressure ranges identified below must be maintained for the products described in this IFU. Otherwise product damage could result. **Temperature:** Operating Temperature Range: -5°C to +60°C (23°F to 140°F) at 20% r.H. to 95% r.H.; Storage Temperature Range: -40°C to +70°C (-40°F to 158°F) at 0% r.H. to 95% r.H. **Humidity:** Minimum Operating Humidity at 40°C (104°F) at 15% r.H. 24hrs; Maximum Operating Humidity at 40°C (104°F) at 95% r.H. 24hrs. **Atmospheric Pressure:** Operating Range: 650hPa to 1080hPa; Storage and Transportation Range: 100hPa to 1080hPa. **REORDERING INFORMATION** • 989803173121, 3-lead set, Single-Patient-Use Disposable, Bedside, AAMI; • 989803173131, 5-lead set, Single-Patient-Use Disposable, Bedside, AAMI; • 989803192141, 5-lead set, Single-Patient-Use Disposable, Chest, AAMI; • 989803197511, 6-lead set, Single-Patient-Use Disposable, Bedside, AAMI; • 989803173151, 5-lead set, Single-Patient-Use Disposable, Telemetry, AAMI.

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Copyright © 2022
All rights reserved
**Böblingen GmbH**