

Pangea™ Femur Plating System

Operative Technique

Distal Lateral Femur

5.0 Distal Medial Femur

3.5 Distal Medial Femur



Contents

Indications and contraindications	3
Pangea introduction	4
Section 01: Pangea Overview	5
Trays	6
Color coding	8
Small and large fragment trays	9
Instrumentation	11
Section 02: System Overview Pangea Femur Plating	
System Distal Femur	18
Required trays	19
Femur Plate offering.....	21
Screw offering	22
Plate details	24
Section 03: Surgical Protocol Pangea Femur Plating	
System Distal Femur	27
Preoperative planning	28
Surgical approaches and instrumentation	31
Reduction and plate contouring.....	32
Plate insertion	33
Plate positioning.....	34
Preliminary plate fixation.....	37
Non-locking screw insertion.....	38
Locking screw insertion.....	39
Lag screw technique.....	40
Small fragment femur plate	41
Cables and cable plugs	42
Final constructs	43
Section 04: System Components	46

This publication sets forth detailed recommended procedures for using Stryker devices and instruments. It offers guidance that you should heed, but, as with any such technical guide, each surgeon must consider the particular needs of each patient and make appropriate adjustments when and as required.

This document is applicable to US and Canada.

The surgeon must advise patients of surgical risks and make them aware of adverse effects and alternative treatments.

WARNING

The patient should be advised that the device cannot and does not replicate a normal healthy bone, that the device can break or become damaged as a result of strenuous activity or trauma and that the device has a finite expected service life.

Removal or revision of the device may be required sometime in the future due to medical reasons.



Pangea Platform indications

The Pangea Platform is indicated for the internal fixation and stabilization of bone fractures, osteotomies, and arthrodesis in normal and osteopenic bone, including:

- Diaphyseal, metaphyseal, epiphyseal, extra- and intra-articular fractures
- Non-unions, malunions, and deformities
- Periprosthetic fractures

The Pangea Platform is also indicated for children (2-12 years) and adolescents (12 – 21 years) for the internal fixation and stabilization of bone fractures of the diaphysis and metaphysis in which growth plates have fused or in which growth plates will not be crossed by implants.

Pangea Femur Plating System indications for use

The Pangea Femur Plating System is indicated for the internal fixation and stabilization of femur bone fractures and osteotomies in normal and osteopenic bone, including:

- Diaphyseal, metaphyseal, epiphyseal, extra- and intra-articular fractures
- Non-unions, malunions, and deformities
- Periprosthetic fractures

Compatibility with other systems

Components from the Pangea Femur Plating System may be used with the following systems:

- AxSOS 3
- Stryker Plating System (SPS)
- Dall-Miles cable system
- Pangea Platform

Please remember that the compatibility of different product systems has not been tested unless specified otherwise in the product labeling. Consult instructions for use (www.ifu.stryker.com) for a complete list of potential adverse effects, contraindications, warnings and precautions.

Contraindications

The physician's education, training, and professional judgement must be relied upon to choose the most appropriate device and treatment.

Conditions presenting an increased risk of failure include:

- Any active or suspected latent infection or marked local inflammation in or about the affected area
- Compromised vascularity that would inhibit adequate blood supply to the fracture or the operative site
- Bone stock compromised by disease, infection or prior implantation that cannot provide adequate support and / or fixation of the devices
- Material sensitivity, documented or suspected
- Patients having inadequate tissue coverage over the operative site
- Implant utilization that would interfere with anatomical structures or physiological performance
- Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care
- Other medical or surgical conditions which would preclude the potential benefit of surgery

WARNING

Follow the instructions provided in our cleaning and sterilization guide (OT-RG-1). All non-sterile devices must be cleaned and sterilized before use.


Multicomponent instruments must be disassembled for cleaning. Please refer to the corresponding assembly / disassembly instructions.

WARNING

Stryker bone screws are not approved or intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

Pangea overview

Pangea systems provide small and large fragment plating solutions for fracture treatment. Pangea's plate offerings include both utility and anatomical plates to address various fracture patterns and anatomy. The plates and non-locking screws are produced from titanium alloy (Ti6Al4V ELI), whereas the locking screws are produced from cobalt-chrome alloy (CoCr).




CoCr locking screw and Ti6Al4V ELI plate hole

Variable angle locking technology

Pangea's variable angle locking technology uses a CoCr locking screw, which is harder than the Ti6Al4V ELI plate, allowing for the screwhead's threads to form a definitive locking position in the plate's locking hole by engaging the softer, Ti6Al4V ELI material.

This technology allows the user to aim and lock the screw into the plate within a true 30° cone of the predetermined hole trajectory. The variable angle drill guide provided with the system offers guidance with respect to the limit of the 30° cone. The locking mechanism remains functional for up to three attempts at locking screw insertion.




Universal holes offer 30° cone of angulation

Hybrid LC Holes (locking/compression)

Hybrid LC Holes allow for either active compression with the use of a non-locking screw in the compression section of the hole or variable angle locking with the use of a locking screw in the universal section of the hole. If locking is not desired, the universal section of the hole also accepts non-locking screws. Each Hybrid LC Hole is designed to provide up to 2mm of compression.

Note: Hybrid LC Holes are not available with every plate type. Refer to the "Plate details" page for additional information.



A: Universal: For locking or non-locking screws
B: Compression: For non-locking screws only

Hybrid LC Hole





Section 01


Pangea Overview

Pangea overview


The Pangea systems require the use of a small or large fragment core tray, which contain the necessary instruments and screws to be used for every case. Plates are contained in separate anatomic plate trays, utility plate trays, or optional tray inserts. Some plates are offered sterile packaged only.

Anatomic and utility plate trays contain plates that will correspond to each color code within the small and large fragment core trays. To complete a case, a surgeon will often need an implant tray and its corresponding core tray.

	Small fragment		Large fragment	
Color coding	Purple	Yellow	Orange	Blue
Screw diameters	2.7mm	3.5 / 4.0mm	4.0 / 4.5 / 6.0mm	5.0mm
Screwdriver type				



Pangea Small Fragment Core Tray



Pangea Large Fragment Core Tray



Auxiliary insert and optional inserts


The Pangea small and large fragment core trays are designed to offer users modularity in their set configurations. Each core tray's standard configuration includes an auxiliary insert containing a silicone mat for storage of miscellaneous instrumentation. Listed below are optional inserts that can be ordered separately and placed into the core tray to meet the user's needs.

When using an optional insert with the Pangea small or large fragment core tray, the auxiliary insert may be replaced with one of the optional inserts.


	Auxiliary insert with silicone mat	Small fragment reduction insert	Large fragment reduction insert	Small fragment standard plate insert	Small fragment ankle plate insert	Asnis III 4.0mm cannulated screw insert
Small fragment core tray	✓	✓		✓	✓	✓
Large fragment core tray	✓		✓			✓

Optional insert configurations


Small fragment utility plate insert



Small fragment ankle plate insert



Asnis III 4.0mm cannulated screw insert
















Note: Optional inserts may not be available in all markets. Check with your local Stryker sales representative for availability.







Color coding

The Pangea systems are color-coded to allow the user to easily identify the proper instrumentation for a particular plate type or screw type. Each color represents the proper drill, drill guide, or screwdriver for a particular screw diameter. The small fragment is color-coded purple for 2.7mm screws and yellow for 3.5mm and 4.0mm screws. The large fragment is color-coded orange for 4.0/4.5/6.0mm screws and blue for 5.0mm screws.

The small and large fragment core trays offer short and long drill bits to account for various depths of bone stock. These ORIF pilot drills are calibrated for the surgeon to measure depth by referencing the associated drill guide. ORIF pilot drill bits are identified by one colored stripe, while lag screw overdrills have one colored stripe and one black stripe. The drill bit's diameter can be found on the AO quick connect.


	 Small fragment		 Large fragment	
Screw type	2.7mm cortex 2.7mm locking	3.5mm cortex 3.5mm locking 4.0mm cancellous	4.0mm locking 4.5mm cortex 6.0mm cancellous	5.0mm locking
Color code	Purple	Yellow	Orange	Blue
ORIF pilot drills	 2.0	 2.5	 3.2	 4.3
Lag screw overdrills	 2.7	 3.5	 4.5	N/A
ORIF pilot drill lengths	Short: 135mm Long: 175mm	Short: 135mm Long: 215mm	Short: 145mm Long: 215mm	Short: 145mm Long: 215mm
ORIF pilot drill calibrations	Short: 0-40mm Long: 0-80mm	Short: 0-40mm Long: 0-120mm	Short: 0-50mm Long: 0-120mm	Short: 0-50mm Long: 0-120mm
Drill guides				

Small fragment core tray set content – 2.7 / 3.5 / 4.0mm screws

T8 torx interface for all 2.7mm screws	T15 torx interface for all 3.5mm screws and 4.0mm cancellous screws
	
 <p>Ø2.7mm locking Ø2.7mm cortex</p> 	 <p>Ø3.5mm locking Ø3.5mm cortex Ø4.0mm cancellous Partially threaded Ø4.0mm cancellous Fully threaded (Variable thread length)</p> 
2.7mm lag screw drill guide	3.5mm lag screw drill guide
2.0mm pilot drill for 2.7mm screws 2.7mm overdrill for 2.7mm lag screws	2.5mm pilot drills for 3.5 / 4.0mm screws 3.5mm overdrill for 3.5mm lag screws
Torque limiting handle, 1.7Nm for 2.7mm locking screws	Torque limiting handle, 4.0Nm for 3.5mm locking screws
Depth gauge for T8 2.7mm screws	Depth gauge for T15 3.5 / 4.0mm screws
Variable angle drill guides	
Compression drill guide	
Fixed angle guides with drill sleeve inserts	
Small and large Delta screwdriver handles with AO quick connect	
Calibrated soft tissue elevator with distal suture hole for alternative plate insertion techniques	

Large fragment core tray set content – 4.0 / 4.5 / 5.0 / 6.0mm screws

T20 torx interface for all large fragment screws



4.5mm lag screw drill guide

3.2mm pilot drills for 4.0 / 4.5 / 6.0mm screws
4.3mm pilot drills for 5.0mm screws
4.5mm overdrill for 4.5mm lag screws

Torque limiting handle, 6.0Nm for 4.0mm and 5.0mm locking screws

Depth gauge for T20 4.0/4.5/5.0/6.0mm screws

Variable angle drill guides

Compression drill guide

Fixed angle guides with drill sleeve inserts

Large Delta screwdriver handles with AO quick connect

Calibrated, curved soft tissue elevator with distal suture hole for alternative plate insertion techniques

Fixed angle drill sleeve

First, the appropriate drill sleeve insert is loosely inserted into the fixed angle sleeve. Next, the entire assembly is inserted into the desired plate hole, which can be confirmed for proper placement with tactile feedback upon insertion. Lastly, turn the drill sleeve insert clockwise to tighten. This fastens the assembly securely to the plate (Fig. 1).

The drill sleeve insert is meant to be hand tightened only. The hex interface on the T8, T15, and T20 screwdriver bits facilitate removal of each respective drill sleeve insert (Fig. 3).




Fig. 1: Fixed angle drill sleeve and insert assembly

Threaded guide post

The T15 and T20 threaded guide posts are intended only for the threaded monoaxial hole in some Pangea plates and serve as the attachment point for MIS targeting capabilities (Fig. 2). Additionally, the threaded guide posts may be utilized as a plate insertion handle, joystick, K-wire sleeve, and drill guide for plates with a threaded monoaxial hole. The threaded guide post is meant to be hand tightened only (Fig. 2).

Note: Threaded monoaxial holes are not available with every plate type. Refer to the “Plate details” page for additional information.

The hex interface on the T20 screwdriver bit facilitates removal of the T20 threaded guide post. The T15 threaded guide post does not have a hex interface and may be removed using the torx portion of the T15 screwdriver bit (Fig. 3).




Fig. 2: Threaded guide post placed in threaded monoaxial hole




Fig. 3: Hex interface between screwdriver bit, fixed angle drill insert, and T20 threaded guide post

Variable angle drill guide

The ball-and-cone variable angle drill guides are used in combination with their respective drill bit to determine screw trajectory. The color-coded guide restricts the degree of screw angulation to 15° in any direction resulting in a 30° cone of the predetermined hole trajectory.


When using the ball end of the guide, gently press the guide into the plate's universal or Hybrid LC Hole. The ball end of the drill guide can be gently rotated in these holes using the handle while maintaining 15° of angulation. For small fragment indications, to ensure a precise 15° angulation, use the cone end of the variable angle drill guide by engaging the cone end of the guide into the plate hole.



Variable angle drill guide

Lag screw drill guide


Lag screw drill guides are offered to accommodate 2.7mm, 3.5mm, or 4.5mm lag screws. The color-coded guide features an overdrill sleeve with serrated tip for near cortex drilling and a self-centering pilot drill sleeve to complete the lag screw pilot hole. Lag screws may be placed through a plate or independently.



Lag screw drill guide

Compression drill guide







The color-coded compression drill guide is an eccentric drill guide used to apply compression in plates equipped with Hybrid LC Holes. Arrows are inscribed on the guide and must be aimed towards the fracture site when engaged into the plate's Hybrid LC Hole. This allows the user to properly drill the Hybrid LC Hole and achieve up to 2mm of compression upon insertion of a non-locking screw.







Compression drill guide

Handles

Pangea offers multiple handle options. These options include small and large Delta handles, and three torque limiting handles for locking screws. All handles are equipped with a small AO-coupling.

Screwdriver type	Locking screw Ø	Torque	Torque limiting handles
	2.7mm	1.7Nm	
	3.5mm	4.0Nm	
	4.0mm 5.0mm	6.0Nm	

Screwdriver type	Non-locking screw Ø	Small and large Delta handles
	2.7mm 3.5mm 4.0mm 4.5mm 6.0mm	
	2.7mm 3.5mm 4.0mm 4.5mm 6.0mm	

Depth measuring instruments

Pangea systems provide two options for determining screw length.

Drill calibrations

All pilot drills are calibrated to allow the surgeon to determine the appropriate screw length when drilling through either the fixed angle sleeve assembly, threaded guide post, or the ball end of the variable angle drill guide (Fig. 1). The cone end of the variable angle drill guide is not compatible with the calibrated drill bits to determine screw length. The calibrations when read against the selected drill guide measure the distance to the tip of the drill bit.



Fig. 1: Measuring screw length using calibrated drill bit

Depth gauge

A depth gauge can be used independently or through a plate hole. Depth gauges correspond with the screw head size i.e., T8, T15, or T20 (Fig. 2).



Fig. 2: T8, T15, and T20 depth gauges

Screw selection

To verify a screw's length, use the screw measurement scale found on the screw rack. The screw tip is placed with its tip against the back stop and its length can be read off the scale (Fig. 3).



WARNING

Always check the correct position and length of the inserted screws by fluoroscopy.




Fig. 3: Screw rack with screw measurement scale

Note: The T8 and T15 depth gauges are color-coded purple and yellow, respectively. The T20 depth gauge is not color-coded and is intended to be used for all T20 screws.

Screw capture sleeve

Pangea offers T8, T15, and T20 screw capture sleeves. The screw capture sleeves are optional devices available to offer efficient screw pick-up, insertion, and removal.

First, the T8, T15, or T20 long screwdriver shaft is inserted into its respective screw capture sleeve. Next, the screw is loaded onto the screwdriver shaft and the knob on the screw capture sleeve is depressed to capture the screw head. Once the screw is captured by the inner sleeve, the knob may be released and the screw is successfully captured. The screw may now be inserted into the pilot hole. Lastly, the screw is released by depressing the button on the end of the device prior to final tightening (Fig. 1).




Fig. 1: Screw capture sleeve assembly

Countersinks and Taps

Countersinks are available for all screws sizes to reduce screw head prominence when the screw is used independently from a plate (Fig. 2).

Taps are available for all screw types and diameters. All screws are self-tapping; however, when inserting a screw in hard bone, it is recommended to use the appropriate tap prior to screw insertion (Fig. 3).



Fig. 2: Countersink



Fig. 3: Tap

Temporary plate fixator

Temporary plate fixators may be used as a provisional plate fixation device and should be limited to the shaft holes of the plate. The temporary plate fixator functions by pushing the shaft of the plate to the bone. The temporary plate fixator is designed with a self-drilling, self-tapping tip for quick insertion into cortical bone (Fig. 4).




Fig. 4: Temporary plate fixator

Small fragment reduction instruments and retractors

These instruments facilitate fracture reduction and soft tissue management.



Periosteal elevator, round edge (705294)



Periosteal elevator, flat edge (705295)



Sharp hook (700151)



Straight reduction clamp (705297) Lobster claw (702932)



Ball spike (700153)



Hohmann retractor 6mm (700664)



Hohmann retractor 15mm (700667)

Large fragment reduction instruments and retractors

These instruments facilitate fracture reduction and soft tissue management.



T-handle large, AO coupling (700367)



Repositioning forceps, 205mm (702927)

Reduction forceps with serrated jaws (702940)




Reduction pin, 5.0mm (390083/390084)




Sharp hook (700151)



Ball spike (700153)



Periosteal elevator, round edge (705294)



Periosteal elevator, flat edge (705295)

Note: Large fragment reduction instruments are contained in an optional large fragment reduction insert tray and may not be available in all markets. Check with your local Stryker sales representative.

Bending instruments and plate cutters

These instruments are available to aid in plate bending and plate cutting. Please refer to the table for the bending and cutting specifics of each instrument.

⚠ WARNING

Repetitively bending or rebending the plates may lead to early fatigue failure or inability to lock at a screw hole. Bending the plate at the screw hole risks compromising the locking mechanism.

⚠ WARNING

In most cases the precontoured plate will fit without the need for further bending. However, should additional out of plane bending of the plate be required, refer to the following table for bending devices and respective plate sizes.

Device Type	2.7 plates	3.5 plates	5.0 plates
Plate bender	✓		
Bending irons	✓	✓ *	
Plate cutter	✓		
Table plate bender		✓	✓

*Bending irons cannot be used to bend the following 3.5 plates: straight broad, extra articular distal humerus, proximal lateral tibia, or distal medial femur




Fig. 3: Table Plate Bender




Fig. 2: Plate Cutter



Fig. 1: Bending Iron



Fig. 4: Plate Bender

Section 02

System Overview


Pangea Femur Plating System

Distal Femur

Required trays

ORIF distal lateral femur

This page details the trays required to use the Pangea Distal Lateral Femur Plates.



Pangea Large Fragment Core Tray



Pangea Distal Lateral Femur Plate Tray




Distal Lateral Femur Plate

Required trays

ORIF distal medial femur

This page details the trays required to use the Pangea Distal Medial Femur Plates.


3.5 Distal Medial Femur Plates




Pangea Small Fragment Core Tray



5.0 Distal Medial Femur Plates




Pangea Large Fragment Core Tray



Pangea Distal Medial Femur Plate Tray



3.5 Distal Medial Femur Plate



5.0 Distal Medial Femur Plate



Implants - Plates

Distal Lateral Femur Plates

160mm
4 holes196mm
6 holes231mm
8 holes267mm
10 holes302mm
12 holes337mm
14 holes371mm
16 holes*406mm
18 holes*441mm
20 holes*

5.0 Distal Medial Femur Plates

151mm
4 holes189mm
6 holes226mm
8 holes*263mm
10 holes*301mm
12 holes*337mm
14 holes*

3.5 Distal Medial Femur Plates

96mm
2 holes110mm
3 holes124mm
4 holes152mm
6 holes

* Available sterile packed only

Screw offering

Implants - Screws

The following screws are suitable for use with the Pangea Femur Plating System:



Description	Ø4.5mm cortex	Ø6.0mm cancellous Partially threaded (16mm thread)	Ø6.0mm cancellous Partially threaded (32mm thread)	Ø6.0mm cancellous Fully threaded	Ø4.0mm locking	Ø5.0mm locking	Ø5.0mm locking Flat tip
Length (increments)	14 - 150mm 14 - 50mm (2mm) 50-150mm (5mm)	30-150mm (5mm)	45-150mm (5mm)	20-150mm (5mm)	14-95mm 14 - 50mm (2mm) 50 - 95mm (5mm)	14-120mm 14-50mm (2mm) 50-120mm (5mm)	10-20mm (2mm)
Material	Titanium alloy (Ti6Al4V ELI)				Cobalt-chrome alloy (CoCr)		
Locking torque	N/A				6Nm		
Washer	Washer T20 (663201)				N/A		
Screwdriver interface	 Screwdriver bit, AO, T20 short, 93mm (542066) Screwdriver bit, AO, T20 long, 180mm (542067)						
Angulation	+/-15° (30° cone)						
Pilot drill Ø	Ø3.2mm				Ø4.3mm		
Drill bits	Ø3.2mm x 145mm (542050) Ø3.2mm x 215mm (542051)				Ø4.3mm x 145mm (542052) Ø4.3mm x 215mm (542053)		
Taps	Tap, locking, Ø5mm x 145mm (542057) Tap, locking, Ø4mm x 145mm (542024) Tap, cortex, Ø4.5mm x 145mm (702808) Tap, cancellous, Ø6mm x 180mm (705054)						




The Pangea cable plug (ref 662202S for 5.0mm plates) helps ensure positioning of a cerclage cable (Dall-Miles cable system) on the plate and can prevent slipping in oblique cable applications. These are available sterile only.

Note: Orange and blue represent the color code for the T20 screws. Ensure that the drill guides, drills, screwdrivers, and depth gauges correspond to the orange or blue color code.

T15 Implants - Screws

The following screws are suitable for use with the Pangea Femur Plating System:



Description	Ø3.5mm cortex	Ø4.0mm cancellous Partially threaded (Variable thread length)	Ø4.0mm cancellous Fully threaded	Ø3.5mm locking
Length (increments)	10-120mm 10-50mm (2mm) 50-120mm (5mm)	10-100mm 10-50mm (2mm) 50-100mm (5mm)	10-100mm 10-50mm (2mm) 50-100mm (5mm)	10-120mm 10-50mm (2mm) 50-120mm (5mm)
Material	Titanium alloy (Ti6Al4V ELI)			Cobalt-chrome alloy (CoCr)
Locking torque	N/A			4Nm
Washer	Washer T15 (663001)			N/A
Screwdriver interface	 Screwdriver bit, AO, T15 short, 93mm (542031) Screwdriver bit, AO, T15 long, 180mm (542032)			
Angulation	+/-15° (30° cone)			
Pilot drill Ø	Ø2.5mm			
Drill bits	Ø2.5mm x 135mm (542020) Ø2.5mm x 215mm (542021)			
Taps	Tap, locking, Ø3.5mm x 125mm (542023) Tap, cortex, Ø3.5mm x 125mm (702802) Tap, cancellous, Ø4mm x 125mm (702803)			




The Pangea cable plug (ref 541400S for 3.5mm plates) helps ensure positioning of a cerclage cable (Dall-Miles cable system) on the plate and can prevent slipping in oblique cable applications. These are available sterile only.

Note: Yellow represents the color code for the T15 screws. Ensure the drill guides, drills, screwdrivers, and depth gauges correspond to the yellow color code.

Distal Lateral Femur Plate

- 8 distal T20 screw holes including 2 kickstand screw holes and 1 threaded monoaxial hole
- 2 variable angle kickstand holes with a predetermined trajectory targeting the medial femoral condyle
- 1 threaded monoaxial hole for use with T20 threaded guide post and serves as a parallel reference to the joint line when using a K-wire, drill, or screw
- Staggered shaft holes to accommodate screw placement around a nail or hip prosthesis
- 2.0mm K-wire holes
- Converging and diverging screw trajectories in the metaphyseal screw holes to increase screw spread in the articular segment
- Rounded and tapered ends designed for insertion under soft tissue

Screw options

 T20 – 5.0mm and 4.0mm locking screws, 4.5mm cortex screws, and 6.0mm cancellous screws.

Note: 4.0mm locking screws are not compatible with the T20 threaded monoaxial hole.


Variable angle locking technology

Universal holes that accept both locking and non-locking screws within a 30° cone.

SOMA designed

These plates are designed using Stryker Orthopedics Modeling and Analytics (SOMA) software.¹


- Evidence-based screw hole placement designed to offer a wide range of trajectory options for particular anatomy¹
- Anatomical fit, creating contoured plates with left and right specific options¹



5.0 Distal Medial Femur Plate

- 5 distal T20 screw holes
- 2.0mm K-wire holes
- Proximal end of the plate twists to lie on the anterior femur
- May be used as medial based fixation for meta-diaphyseal fractures or as a supplement to lateral plate fixation when dual plate fixation is desired
- Rounded and tapered ends designed for insertion under soft tissue

Screw options

 T20 – 5.0mm and 4.0mm locking screws, 4.5mm cortex screws, and 6.0mm cancellous screws.

Variable angle locking technology

Universal holes that accept both locking and non-locking screws within a 30° cone.

SOMA designed

These plates are designed using Stryker Orthopedics Modeling and Analytics (SOMA) software.¹

- Evidence-based screw hole placement designed to offer a wide range of trajectory options for particular anatomy¹
- Anatomically influenced plate fit, creating contoured plates with left and right specific options¹



3.5 Distal Medial Femur Plate

- 10 distal T15 screw holes
- 2.0mm K-wire holes
- May be used for isolated medial condylar fractures or as a supplement to lateral plate fixation when dual plate fixation is desired
- Rounded and tapered ends designed for insertion under soft tissue

Screw options

- **T15** T15 – 3.5mm locking and cortex screws and 4.0mm cancellous screws.


Variable angle locking technology

Universal holes that accept both locking and non-locking T15 screws within a 30° cone.

SOMA designed

These plates are designed using Stryker Orthopedics Modeling and Analytics (SOMA) software.¹

- Evidence-based screw hole placement designed to offer a wide range of trajectory options for particular anatomy¹
- Anatomically influenced plate fit, creating contoured plates with left and right specific options¹



Section 03

Surgical Protocol

Pangea Femur Plating System

Distal Femur

Plate sizes

Selecting the appropriate plate length involves considerations including bone quality, fracture configuration and location, and the type of bone healing expected.





**Distal Lateral
Femur Plate****3.5 Distal Medial
Femur Plate****5.0 Distal Medial
Femur Plate**

Plate laser markings

The plate's reference number and left/right orientation is written on the distal end of the plate. The Distal Lateral Femur Plate also has the number of holes written on the plate. Additionally, the number of plate holes, reference number, left/right orientation, and plate length are written on the implant tray.





Preoperative planning

Digital templates

The use of digital templates in association with adequate X-ray/ fluoroscopy or CT scans may assist in the selection of an appropriately sized implant. Pangea digital templates are licensed to template service providers e.g., PeekMed, Sectra AB, Lexi, and Merge Healthcare. Template service providers provide software tools and access to 3D models in order to allow surgeons to perform preoperative planning.




Preoperative planning using 3D planning software




Scale marker

Details: Scale, marker, and product description



Rotational marker

Example of 2D femur templates using ".dxf" files



3D scanned bone model

".stl" plate model

".igs" axis data

".stl" screw model

Example of 3D femur templates using ".stl" and ".igs" files

NOTICE If digital X-ray images are used, correct magnification has to be verified prior to use.

Note: Digital templates may also be accessed by contacting your local Stryker sales representative.

Surgical approaches

The approaches described below are common approaches for the insertion of the selected distal femur plate. These plates can be inserted through a variety of approaches based on surgeon preference, soft tissue, and fracture characteristics.

Distal lateral femur

Surgeons may use a standard lateral, modified lateral (anterolateral), or lateral parapatellar approach to the distal lateral femur (Fig.1).




Fig. 1: Lateral approaches to the distal femur

Distal medial femur

Surgeons may use a standard medial, medial parapatellar, or medial subvastus approach to the distal medial femur (Fig. 2).



Fig. 2: Medial approaches to the distal femur

Periosteal elevators

The soft tissue elevators can be utilized to separate or lift the soft tissues from the bone (Fig. 3). The curved soft tissue elevator has been designed to create a pathway for the plate and is calibrated allowing the surgeon to monitor depth or measure plate length (Fig. 4).

The surgeon should not use the curved soft tissue elevator to tunnel directly on the medial or posteromedial side of the femur but allow the twist of the plate to guide the implant anteriorly.




Fig. 3: Periosteal elevators




Fig. 4: Curved soft tissue elevator, T20

Part #	Instrumentation
706416	Curved soft tissue elevator, T20
705294	Periosteal elevator, round edge 6mm
705295	Periosteal elevator, straight edge 13mm

Reduction

Indirect fracture fragment reduction can be achieved through traction, manipulation, external fixation, or the surgeon's preferred technique. Direct fragment reduction can be achieved using elevators, clamps, or K-wires. Anatomic reduction of the fracture can be achieved through direct visualization and use of clamps.

The use of K-wires may be helpful in maintaining provisional reduction. Independent lag screws may definitively stabilize articular reduction prior to plate insertion (Fig. 1).




Fig. 1: Reduction of distal femur

Plate contouring

Should bending of the plate be required, the Table Top Bender (ref 702900) may be used (Fig. 2).




Fig. 2: Table Top Bender

Plate insertion

When a plate insertion handle is desired, the surgeon may use either the T20 threaded guide post or the T20 fixed angle sleeve assembly. The T20 threaded guide post is attached by threading it into the plate's threaded monoaxial hole. The T20 fixed angle sleeve assembly may be attached to any of the plate's universal holes (Fig. 1).

Using the T20 threaded guide post or T20 fixed angle sleeve assembly, the surgeon may insert the selected distal femur plate submuscularly (Fig. 2). For submuscular insertion, slide the plate proximally using the plate to feel the femur to confirm the plate is not being inserted anterior or posterior to the femoral shaft. Once the plate has been inserted, the surgeon may use either plate insertion handle as a joystick to adjust plate positioning (Fig. 3).

Note: T20 threaded guide post is only compatible with the Distal Lateral Femur Plate.

Assembly instructions of the fixed angle sleeve and drill sleeve insert can be found on page 11.




Fig. 1: T20 threaded guide post or T20 fixed angle sleeve assembly used as a plate insertion handle or joystick

CAUTION

Avoid plate insertion through the muscle to prevent intramuscular vessel disruption. Minimize periosteal disruption while inserting the plate to help preserve bone blood supply.




Fig. 2: Plate insertion using the T20 threaded guide post as a plate insertion handle




Fig. 3: T20 fixed angle sleeve assembly used as a joystick for final plate adjustment

Distal Lateral Femur Plate positioning

Once reduction is achieved, the plate is placed on the lateral surface of the distal femur.

These plates have been designed using SOMA software¹. In most cases the precontoured plate will fit without the need for additional adjustments.

Confirm that the capsule edges and iliotibial band are not trapped under the plate, as these layers will need to be available for layered wound closure. Confirm that the plate is submuscular, not intramuscular.

The plate is in the proper position when it is sitting adjacent to the anterior margin of the articular surface.

The Distal Lateral Femur Plate's threaded monoaxial hole serves as a reference hole and its predetermined trajectory should be parallel with the joint line. The T20 threaded guide post and 2.0mm K-wire sleeve may be used to provisionally fixate the plate to the femur (Fig. 1).

Due to the anterolateral positioning of the plate, it is recommended to check plate position with a true medial-lateral view of the **plate**, as a true medial-lateral view of the **femur** will show the plate appearing too anterior (Fig. 2).




Fig. 1: Distal Lateral Femur Plate reference hole




Fig. 2: Positioning of the Distal Lateral Femur Plate



Plate positioning



3.5 Distal Medial Femur Plate

Once reduction is achieved, the plate is placed on the anteromedial surface of the distal femur.

These plates have been designed using SOMA software¹. In most cases the precontoured plate will fit without the need for additional adjustments.

The plate is in the proper position when the distal end of the plate is adjacent to the medial articular margin of the trochlea and anterior to the medial epicondyle.

The surgeon must be familiar with the vascular anatomy of the thigh to ensure that proximal screw insertion will not put vascular structures at risk.




5.0 Distal Medial Femur Plate

Once reduction is achieved, the plate is placed on the anteromedial surface of the distal femur with the plates shaft twisting anteriorly along the femur.

These plates have been designed using SOMA software¹. In most cases the precontoured plate will fit without the need for additional adjustments.

The plate is in the proper position when the distal end of the plate is adjacent to the medial articular margin of the trochlea and anterior to the medial epicondyle.

The surgeon must be familiar with the vascular anatomy of the thigh to ensure that proximal screw insertion will not put vascular structures at risk.



Provisional plate fixation

Provisional plate fixation can be achieved in various ways:

- 2.0mm K-wire placed through the K-wire holes
- 2.0mm K-wire placed through the T20 fixed angle sleeve and drill sleeve insert or T20 threaded guide post utilizing the 2.0mm K-wire sleeve insert
- Temporary plate fixator

The temporary plate fixator is designed to provide provisional fixation and can be used to push the shaft of the plate to the bone (Fig. 2). There is a self-drilling, self-tapping tip for quick insertion into cortical bone. Bicortical purchase should be confirmed with fluoroscopy.

In order to protect surrounding soft tissues during pin insertion, the temporary plate fixator sleeve must be preassembled onto the temporary plate fixator pin with the self-drilling tip of the pin being flush with the tip of the sleeve (Fig. 3).

Once the device is inserted through the far cortex (Fig. 4), the threaded outer sleeve resting on the plate is turned clockwise until the desired plate position is achieved (Fig. 5).

If replacing the temporary plate fixator with definitive screw fixation is required, the surgeon must re-drill the hole using the appropriate drill and drill guide. This ensures that the pilot hole for the screw is the proper diameter and within the 30° cone.

⚠ CAUTION

Be careful when using sharp instruments such as drills, taps, K-wires and temporary plate fixators or when inserting screws to avoid damage to the soft tissue or vessels by going too far past the far cortex.




Fig. 1 Provisional fixation of Pangea Distal Lateral Femur Plate





Fig. 2: Temporary plate fixator



Note: The Ø2.0mm K-wire sleeve insert is only compatible with the T20 fixed angle sleeve assemblies and T20 threaded guide post.

Non-locking screw insertion

Insertion of a non-locking screw is started with the insertion of the appropriate drill guide for the screw hole.

To achieve the predetermined trajectory of the universal holes, use the T20 fixed angle sleeve with the 3.2mm drill sleeve insert. Assembly instructions of the fixed angle sleeve and drill sleeve insert can be found on page 11. When using the threaded monoaxial hole on the Distal Lateral Femur Plate, the T20 threaded guide post must be used (Fig. 1). If a variable angle trajectory is desired, use the T20, 3.2mm variable angle drill guide (Fig. 2). Variable angle drill guides are not compatible with the threaded monoaxial hole.

Next, using the 3.2mm drill bit, create a pilot hole by drilling through the selected drill guide (Fig. 3). The drill trajectory may be verified under fluoroscopy if required.

The depth may be measured utilizing the T20 depth gauge or the drill bit calibrations. The selected screw is then inserted into the pilot hole using the T20 screwdriver bit (Fig. 4).

A T20 screw capture sleeve may be used to aid in retention between the screw and screwdriver shaft during screw insertion.

CAUTION

If excessive resistance is felt during insertion, or if the bone is dense, it is recommended to use a tap.

CAUTION

Use bi-cortical fixation when possible.

CAUTION

It is recommended to insert the screws by hand. If power tools are used, use those at low speed to avoid improper alignment.




Fig. 1: Drill guides for predetermined screw trajectory




Fig. 2: Variable angle drill guide for variable angle drilling




Fig. 3: Drill pilot hole with 3.2mm drill bit




Fig. 4: Screw placement with T20 screwdriver

Locking screw insertion

Locking screws can be placed within a 30° cone in any universal screw hole.

Uni-cortical fixation is recommended when bi-cortical fixation cannot be achieved safely. For example, when bi-cortical fixation will cause screw penetration into the articular surface. In situations of uni-cortical screws, an increased screw count may be needed to obtain sufficient fixation.

To achieve the predetermined trajectory of the universal holes, use the T20 fixed angle sleeve with its respective drill sleeve insert for the desired diameter locking screw. Assembly instructions of the fixed angle sleeve and drill sleeve insert can be found on page 11. When using the threaded monoaxial hole on the Distal Lateral Femur Plate, the T20 threaded guide post must be used (Fig. 1). If a variable angle trajectory is desired, use the appropriate T20 variable angle drill guide (Fig. 2). Variable angle drill guides are not compatible with the threaded monoaxial hole.

Next, using the 3.2mm or 4.3mm drill bit, create a pilot hole by drilling through the selected drill guide (Fig. 3). The drill trajectory may be verified with fluoroscopy if required.

The depth may be measured utilizing either the T20 depth gauge or the drill bit calibrations. The selected screw is then inserted into the pilot hole using the T20 screwdriver bit.

The T20 screw capture sleeve may be used to aid in retention between the screw and screwdriver shaft during screw insertion.

Use the 6Nm torque limiting T-Handle and T20 screwdriver bit to ensure proper seating of the locking screw. The torque limiter will produce an audible “click” when the required torque is achieved (Fig. 4).

⚠ WARNING

Always perform final tightening by hand using the appropriate torque limiter and screwdriver bit as final tightening with the power tool can cause over-torquing or damage to the screw-plate interface, which can lead to breaking or stripping screw heads.




Fig. 1: Drill guides for predetermined screw trajectory




Fig. 2: Variable angle drill guide for variable angle drilling




Fig. 3: Drill pilot hole with appropriate drill bit



Fig. 4: Screw placement with 6Nm torque limiter

Lag screw technique

Independent lag screw

To insert an independent lag screw, select the appropriate lag screw drill guide for the desired lag screw size. Then use the serrated end of the guide and the appropriate overdrill to drill the near cortex (Step 1a).

Next, insert the self-centering end of the lag screw drill guide into the gliding hole and use the appropriate pilot drill to drill through the far cortex (Step 2a).

Measure the depth of the hole using the pilot drill or appropriate depth gauge and insert the selected cortex screw. If required, countersink the gliding hole or apply a washer to the screw. Upon screw insertion, this technique will serve to lag the far cortex towards the near cortex, thus applying compression (Step 3a).

Lag screw through plate

To insert a lag screw through a plate hole, use the serrated end of the appropriate lag screw drill guide and its respective overdrill to drill the near cortex through the plate hole (Step 1b).

Next, drill the far cortex with the appropriate size pilot drill by placing either the variable angle drill guide into the plate hole or the self-centering end of the lag screw drill guide into the gliding hole through the plate. When drilling the far cortex using a variable angle drill guide, ensure the trajectory of the pilot drill is co-linear with the gliding hole (Step 2b).

Measure the depth of the hole using the pilot drill or appropriate depth gauge and insert the selected cortex screw. Upon screw insertion, this technique will serve to lag the far cortex towards the near cortex, thus applying compression through the plate. Screw holes in the plate may be populated to complete the construct (Step 3b).




T15 3.5 Distal Medial Femur Plate

When using the 3.5 Distal Medial Femur Plate, a surgeon may follow the same steps outlined in previous pages for plate insertion, provisional fixation, screw insertion, and screw measurement utilizing the instruments in Fig. 1.

The 3.5 Distal Medial Femur Plate comprises T15 holes that accommodate 3.5mm locking screws, 3.5mm cortex screws and 4.0mm cancellous screws. The instrumentation is color-coded yellow and is contained in the Pangea Small Fragment Core Tray.

Ref #	T15 Instrumentation
542036	Ø1.6mm x 100mm Olive tip K-wire
542025 542026	Fixed angle sleeve, T15 Drill sleeve insert, T15, 60mm, Ø2.5mm
542028	Variable angle drill guide, T15, Ø2.5mm
542020 / 542021	Drill bit, AO, Ø2.5mm x 135mm, x 215mm
542017	Depth gauge, T15, 0-120mm
542031 / 542032	Screwdriver bit, AO, T15, 93mm / 180mm
542034	Torque limiting T-handle, T15, 4Nm
541400S	3.5mm cable plug, T15

Fig. 1: T15 Instrumentation for 3.5mm Distal Medial Femur Plate



Cables and cable plugs

Cable plug insertion

In any universal hole, cable plugs can be utilized in conjunction with Ø2.0mm vitallium cables from the Dall-Miles cable system.

The use of a cable plug helps to ensure positioning of a cerclage cable (Dall-Miles cable system) on the plate and can prevent slipping in oblique cable applications. These are only available in sterile packaging.

Insert a cable plug by “clicking” it into the desired screw hole and confirming its placement through an audible “click”. The cable plug is still able to freely rotate after proper placement in the universal hole.


Insert a cable through the eyelet of the cable plug, tighten the cable, and crimp the sleeve. Cut the cable near the crimped sleeve.

For complete instructions, refer to the Dall-Miles cabling system operative technique (DALLM-SP-1_21060).

Ref #	Instrumentation
541400S	3.5mm cable plug
662202S	5.0mm cable plug
6704-0-520	Ø2.0mm Vitallium beaded cable / sleeve set
6704-9-320	Single sided tensioner
6704-9-150	Crimp tool
6704-9-420	Cable cutter

Cable plug removal

If a cable plug has to be removed, simply cut or remove the cable and pull out the cable plug using forceps. The cable plug can be re-seated up to 3 times intraoperatively.



CAUTION

Do not exceed plate/cable plug relative angulation of 35° from the screw hole axis for universal holes and 25° from the screw hole axis for hybrid holes.



Final fluoroscopic check

After final fixation with all screws, the femur should be internally and externally rotated under continuous fluoroscopy; the true distance of all screws should be inspected in the AP and lateral views to ensure none of the screws have penetrated the joint or are prominent.



WARNING

After the procedure, check that all implants are positioned correctly using an image intensifier.

Implant removal

Removal of the Pangea Distal Femur Plates is not required in general. The additional surgical trauma and risks associated with additional anesthesia should be individually outweighed against the potential benefits for every patient.

In the case of implant removal, the scar of the previous incision is (partly) re-opened and the screws and plates are successfully removed.

Nail Plate Combination²

Femoral nails e.g., T2 Alpha Femur Antegrade GT/PF or T2 Alpha Femur Retrograde nails may be used in combination with any Pangea, PeriPRO, or AxSOS 3 femur or utility plates.


Note: Care must be taken to ensure all supplementary implants, such as nails or additional plates and screws, do not come into contact with the Pangea plate construct.



Note: In the extreme event of broken or stripped screws, the Stryker implant extraction set (literature number IES-ST-1) includes a variety of broken screw removal instruments.

Final construct examples

Distal Lateral Femur Plate




Predetermined trajectory




Screws in the predetermined screw trajectory

Variable angle trajectory



Screw trajectories using variable angle locking to obtain the widest allowable screw trajectory

5.0mm Distal Medial Femur Plate




Predetermined trajectory



Screws in the predetermined screw trajectory



Variable angle trajectory




Screw trajectories using variable angle locking to obtain the widest allowable screw trajectory

Final construct examples

3.5 Distal Medial Femur Plate

**Predetermined trajectory**

Screws in the predetermined screw trajectory

Variable angle trajectory

Screw trajectories using variable angle locking to obtain the widest allowable screw trajectory


Section 04

System Components

System components

Pangea Distal Lateral Femur Plate Tray

The information in this section is not intended to be used for sales and/or promotional purposes. This information is solely intended to be used as a reference for clinical usage.



Ref #	Description
542228	Distal Lateral Femur Plate Tray

Distal Lateral Femur Plates

Left Ref #	Right Ref #	Plate Length	Shaft Holes
540144	540174	160mm	4
540146	540176	196mm	6
540148	540178	231mm	8
540150	540180	267mm	10
540152	540182	302mm	12
540154	540184	337mm	14
540156S*	540186S*	371mm	16
540158S*	540188S*	406mm	18
540160S*	540190S*	441mm	20




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Distal Medial Femur Plate Tray

The information in this section is not intended to be used for sales and/or promotional purposes. This information is solely intended to be used as a reference for clinical usage.



Ref #	Description
542230	Distal Medial Femur Plate Tray
542231	Silicone Mat, Distal Medial Femur Tray

3.5 Distal Medial Femur Plates

Left Ref #	Right Ref #	Plate Length	Shaft Holes
540202	540222	96	2
540203	540223	110	3
540204	540224	124	4
540206	540226	152	6



5.0 Distal Medial Femur Plates

Left Ref #	Right Ref #	Plate Length	Shaft Holes
540244	540264	151	4
540246	540266	189	6
540248S*	540268S*	226	8
540250S*	540270S*	263	10
540252S*	540272S*	301	12
540254S*	540274S*	337	14




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray


Top level consisting of the instruments listed below:



Ref #	Description
542201	Aluminum lid - universal
542203	Small fragment core tray base
542204	Small fragment ORIF instrument insert
542254	Small fragment std. torque handle insert

Pangea Small Fragment Core Tray - ORIF instrument insert

Ref #	Description
542005	Fixed angle sleeve, T8
542006	Drill sleeve insert, T8, 60mm, ø2.0mm
542008	Variable angle drill guide, T8, ø2.0mm
542009	Compression drill guide, T8 / T15, ø2.0mm / ø2.5mm
542010	Lag screw drill guide, T8, ø2.0mm / ø2.7mm
542011	Screwdriver bit, AO, T8, 93mm
542012	Screwdriver bit, AO, T8, 180mm
542016	Depth gauge, T8, 0-80mm
542015	Screw capture sleeve, T8
542017	Depth gauge, T15, 0-120mm
542025	Fixed angle sleeve, T15
542026	Drill sleeve insert, T15, 60mm, ø2.5mm
542028	Variable angle drill guide, T15, ø2.5mm
542030	Lag screw drill guide, T15, ø2.5mm / ø3.5mm
542031	Screwdriver bit, AO, T15, 93mm
542032	Screwdriver bit, AO, T15, 180mm
542035	Soft tissue elevator, T15
45-80040	Countersink for screws 2.7/3.5mm, AO fitting
542103	Threaded guide post, T15
542027	Screw capture sleeve, T15
542098	Delta handle, small, AO, T8 / T15
542099	Delta handle, large, AO, T20
542000	Drill bit, AO, Ø2.0mm x 135mm
542001	Drill bit, AO, Ø2.0mm x 175mm



System components

Pangea Small Fragment Core Tray

Top level consisting of the instruments listed below:

Pangea Small Fragment Core Tray - ORIF instrument insert

Ref #	Description
542002	Drill bit, AO, Ø2.7mm x 125mm
542020	Drill bit, AO, Ø2.5mm x 135mm
542021	Drill bit, AO, Ø2.5mm x 215mm
542022	Drill bit, AO, Ø3.5mm x 135mm
542003	Tap, locking, AO, Ø2.7mm x 125mm
702801	Tap, AO, Ø2.7mm x 125mm
542023	Tap, locking, AO, Ø3.5mm x 125mm
702802	Tap, AO, Ø3.5mm x 125mm
702803	Tap, cancellous, Ø4mm x 125mm
705002	K-wire drill tip, Ø2.0mm x 234mm
390157	K-wire Ø1.25 x 150mm
390164	K-wire Ø1.6 x 150mm
390192	K-wire Ø2.0 x 150mm
542036	K-wire, olive tip, Ø1.6mm / 100mm

Torque limiter insert instruments


Ref #	Description
542014	Torque limiting Delta handle, AO, T8 1.7Nm
542034	Torque limiting T-handle, AO, T15 4.0Nm

System components

Pangea Small Fragment Core Tray - Instruments

Second level consisting of the instruments listed below:

This level can be removed and replaced with the optional inserts: Core tray ankle plate insert, core tray small fragment utility plate insert, or core tray Asnis III 4.0mm cannulated screw insert.




Ref #	Description
542207	Core tray auxiliary insert
542208	Silicone mat, auxiliary insert

System components

Pangea Small Fragment Core Tray - Instruments

Third level consisting of the instruments listed below:



Ref #	Description
542209	Core tray small reduction insert


Pangea Small Fragment Core Tray - Core tray small reduction insert

Ref #	Description
700151	Hook
700153	Ball spike
700664	Hohmann retractor 6mm
700667	Hohmann retractor 15mm
705294	Periosteal elevator, round edge 6mm
705295	Periosteal elevator, flat blade 13mm
705297	Straight reduction clamp, broad
702932	Repositioning forceps, L143mm (lobster claw)
703938	Bending iron
705019	Temporary plate fixator, AO

System components

Pangea Small Fragment Core Tray - Instruments

Third level consisting of the instruments listed below:



Ref #	Description
542209	Core tray small reduction insert

Pangea Small Fragment Core Tray - Core tray small reduction insert

Ref #	Description
700151	Hook
700153	Ball spike
700664	Hohmann retractor 6mm
700667	Hohmann retractor 15mm
705294	Periosteal elevator, round edge 6mm
705295	Periosteal elevator, flat blade 13mm
705297	Straight reduction clamp, broad
702932	Repositioning forceps, L143mm (lobster claw)
703938	Bending iron
705019	Temporary plate fixator, AO

System components

Pangea Small Fragment Core Tray - Optional ankle insert



Ref #	Description
542248	Core tray ankle plate insert

Implants

Ref #	Description	
540644	Distal Lateral Fibula Plate, L	2.7/3.5mm, 81mm/4 holes
540664	Distal Lateral Fibula Plate, R	2.7/3.5mm, 81mm/4 holes
540645	Distal Lateral Fibula Plate, L	2.7/3.5mm, 95mm/5 holes
540665	Distal Lateral Fibula Plate, R	2.7/3.5mm, 95mm/5 holes
540646	Distal Lateral Fibula Plate, L	2.7/3.5mm, 109mm/6 holes
540666	Distal Lateral Fibula Plate, R	2.7/3.5mm, 109mm/6 holes
540647	Distal Lateral Fibula Plate, L	2.7/3.5mm, 123mm/7 holes
540667	Distal Lateral Fibula Plate, R	2.7/3.5mm, 123mm/7 holes
540648	Distal Lateral Fibula Plate, L	2.7/3.5mm, 137mm/8 holes
540668	Distal Lateral Fibula Plate, R	2.7/3.5mm, 137mm/8 holes
541342	3.5 1/3 tubular plate	3.5mm, 26mm/2 holes
541343	3.5 1/3 tubular plate	3.5mm, 38mm/3 holes
541344	3.5 1/3 tubular plate	3.5mm, 50mm/4 holes
541345	3.5 1/3 tubular plate	3.5mm, 62mm/5 holes
541346	3.5 1/3 tubular plate	3.5mm, 74mm/6 holes
541347	3.5 1/3 tubular plate	3.5mm, 86mm/7 holes
541348	3.5 1/3 tubular plate	3.5mm, 98mm/8 holes
541349	3.5 1/3 tubular plate	3.5mm, 110mm/9 holes
541303	2.7 hook plate	2.7mm, 37mm/3 holes
541305	2.7 hook plate	2.7mm, 57mm/5 holes
541321	2.7 T-plate	2.7mm, 3Tx5
541331	3.5 T-plate	3.5mm, 3Tx5

For the full offering of lengths please refer to the fibula and small fragment utility plate optechs.

*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

Pangea Small Fragment Core Tray - Optional Asnis III 4.0mm cannulated screw insert

Ref #	Description
542245	Core tray Asnis III 4.0mm cannulated screw insert
940236	Screw rack Asnis III 4.0mm

Instruments

Ref #	Description
702465	Double drill guide, \varnothing 1.4mm/2.7mm
702446	Cannulated drill, AO, \varnothing 4.0mm
702449	Cannulated drill, AO, \varnothing 2.7mm
702459	Threaded guide wire, \varnothing 1.4x150mm
702454	Cannulated tap, AO, \varnothing 4.0mm
702473	Cannulated countersink, AO, \varnothing 4.0mm screws
702499	Direct measuring gauge for wires \varnothing 1.4/2.0x150mm
702482	Cannulated screwdriver, AO, hex 2.5mm
702485	Solid screwdriver, AO, hex 2.5mm
702489	Holding sleeve for screwdrivers, for screwheads \varnothing 5.0mm
702492	Cleaning stylet \varnothing 1.4mm
702496	Extractor, AO, \varnothing 4.0mm screws

Implants

Ref #	Description
604624	4.0 X 24mm TI cannulated screw
604626	4.0 X 26mm TI cannulated screw
604628	4.0 X 28mm TI cannulated screw
604630	4.0 X 30mm TI cannulated screw
604632	4.0 X 32mm TI cannulated screw
604634	4.0 X 34mm TI cannulated screw
604636	4.0 X 36mm TI cannulated screw
604638	4.0 X 38mm TI cannulated screw
604640	4.0 X 40mm TI cannulated screw
604642	4.0 X 42mm TI cannulated screw

Pangea Small Fragment Core Tray - Optional Asnis III 4.0mm cannulated screw insert**Implants**

Ref #	Description
604644	4.0 X 44mm TI cannulated screw
604646	4.0 X 46mm TI cannulated screw
604648	4.0 X 48mm TI cannulated screw
604650	4.0 X 50mm TI cannulated screw
604655	4.0 X 55mm TI cannulated screw
604660	4.0 X 60mm TI cannulated screw
604665	4.0 X 65mm TI cannulated screw
604670	4.0 X 70mm TI cannulated screw
604724	4.0 X 24mm TI cannulated screw, full thread
604726	4.0 X 26mm TI cannulated screw, full thread
604728	4.0 X 28mm TI cannulated screw, full thread
604730	4.0 X 30mm TI cannulated screw, full thread
604732	4.0 X 32mm TI cannulated screw, full thread
604734	4.0 X 34mm TI cannulated screw, full thread
604736	4.0 X 36mm TI cannulated screw, full thread
604738	4.0 X 38mm TI cannulated screw, full thread
604740	4.0 X 40mm TI cannulated screw, full thread
604742	4.0 X 42mm TI cannulated screw, full thread
604744	4.0 X 44mm TI cannulated screw, full thread
604746	4.0 X 46mm TI cannulated screw, full thread
604748	4.0 X 48mm TI cannulated screw, full thread
604750	4.0 X 50mm TI cannulated screw, full thread
619905	Washer

*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

Pangea Small Fragment Core Tray - Optional small fragment utility plate insert

Ref #	Description
542247	Core tray small utility plate insert

Implants

Ref #	Description	
541342	3.5 1/3 tubular plate	3.5mm, 26mm/2 holes
541343	3.5 1/3 tubular plate	3.5mm, 38mm/3 holes
541344	3.5 1/3 tubular plate	3.5mm, 50mm/4 holes
541345	3.5 1/3 tubular plate	3.5mm, 62mm/5 holes
541346	3.5 1/3 tubular plate	3.5mm, 74mm/6 holes
541347	3.5 1/3 tubular plate	3.5mm, 86mm/7 holes
541348	3.5 1/3 tubular plate	3.5mm, 98mm/8 holes
541349	3.5 1/3 tubular plate	3.5mm, 110mm/9 holes
541303	2.7 hook plate	2.7mm, 37mm/3 holes
541305	2.7 hook plate	2.7mm, 57mm/5 holes
541313	3.5 hook plate	3.5mm, 49mm/3 holes
541315	3.5 hook plate	3.5mm, 79mm/5 holes
541321	2.7 T-plate	2.7mm, 3TX5
541322	2.7 T-plate	2.7mm, 3TX10

Ref #	Description	
541331	3.5 T-plate	3.5mm, 3Tx5
541332	3.5 T-plate	3.5mm, 3Tx10
541036	3.5 straight narrow plate	3.5mm, 81mm/6 holes
541037	3.5 straight narrow plate	3.5mm, 104mm/7 holes
541038	3.5 straight narrow plate	3.5mm, 115mm/8 holes
541039	3.5 straight narrow plate	3.5mm, 138mm/9 holes
541040	3.5 straight narrow plate	3.5mm, 149mm/10 holes
541042	3.5 straight narrow plate	3.5mm, 183mm/12 holes
541004	2.7 straight narrow plate	2.7mm, 31mm/4 holes
541006	2.7 straight narrow plate	2.7mm, 46mm/6 holes
541008	2.7 straight narrow plate	2.7mm, 61mm/8 holes
541010	2.7 straight narrow plate	2.7mm, 76mm/10 holes
541020	2.7 straight narrow plate	2.7mm, 152mm/20 holes

For the full offering of lengths please refer to the fibula and small fragment utility plate optechs.


*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



Ref #	Description
542205	Small fragment screw rack (with lid)
542206	Small fragment screw rack lid

2.7mm locking screw self-tapping, T8 drive

Ref #	Length (mm)	Ref #	Length (mm)
541408	8	541436	36
541410	10	541438	38
541412	12	541440	40
541414	14	541442	42
541416	16	541444	44
541418	18	541446	46
541420	20	541448	48
541422	22	541450	50
541424	24	541455	55
541426	26	541460	60
541428	28	541465	65
541430	30	541470	70
541432	32	541475S*	75
541434	34	541480S*	80




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



2.7mm cortex screw self-tapping, T8 drive

Ref #	Length (mm)	Ref #	Length (mm)
541708	8	541736	36
541710	10	541738	38
541712	12	541740	40
541714	14	541742	42
541716	16	541744	44
541718	18	541746	46
541720	20	541748	48
541722	22	541750	50
541724	24	541755	55
541726	26	541760	60
541728	28	541765	65
541730	30	541770	70
541732	32	541775S*	75
541734	34	541780S*	80




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



3.5mm locking screw self-tapping, T15 drive

Ref #	Length (mm)	Ref #	Length (mm)
541510	10	541546	46
541512	12	541548	48
541514	14	541550	50
541516	16	541555	55
541518	18	541560	60
541520	20	541565	65
541522	22	541570	70
541524	24	541575	75
541526	26	541580	80
541528	28	541585	85
541530	30	541590	90
541532	32	541595	95
541534	34	541600S*	100
541536	36	541605S*	105
541538	38	541610S*	110
541540	40	541615S*	115
541542	42	541620S*	120
541544	44		




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



3.5mm cortex screw self-tapping, T15 drive

Ref #	Length (mm)	Ref #	Length (mm)
661410	10	661446	46
661412	12	661448	48
661414	14	661450	50
661416	16	661455	55
661418	18	661460	60
661420	20	661465	65
661422	22	661470	70
661424	24	661475	75
661426	26	661480	80
661428	28	661485	85
661430	30	661490	90
661432	32	661495	95
661434	34	661500S*	100
661436	36	661505S*	105
661438	38	661510S*	110
661440	40	661515S*	115
661442	42	661520S*	120
661444	44		




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



4.0mm cancellous screw full thread self-tapping, T15 drive

Ref #	Length (mm)	Ref #	Length (mm)
607310	10	607342	42
607312	12	607344	44
607314	14	607346	46
607316	16	607348	48
607318	18	607350	50
607320	20	607355	55
607322	22	607360	60
607324	24	607365	65
607326	26	607370	70
607328	28	607375	75
607330	30	607380	80
607332	32	607385	85
607334	34	607390	90
607336	36	607395	95
607338	38	607400S*	100
607340	40		




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



4.0mm cancellous screw partial thread self-tapping, T15 drive

Ref #	Length (mm)	Thread length (mm)	Ref #	Length (mm)	Thread length (mm)
607410	10	5	607442	42	15
607412	12	5	607444	44	15
607414	14	5	607446	46	15
607416	16	6	607448	48	15
607418	18	7	607450	50	15
607420	20	8	607455	55	16
607422	22	9	607460	60	16
607424	24	10	607465	65	16
607426	26	12	607470	70	16
607428	28	14	607475	75	16
607430	30	14	607480	80	16
607432	32	14	607485	85	16
607434	34	14	607490	90	16
607436	36	14	607495	95	16
607438	38	14	607500S*	100	16
607440	40	14			



Cable plug and washers

Ref #	Description
541400S*	3.5mm cable plug
40-30900	Washer, T8
663001	Washer, T15




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Large Fragment Core Tray


Top level consisting of the instruments listed below:



Ref #	Description
542201	Aluminum lid - universal
542211	Large fragment core tray base
542212	Large fragment ORIF instrument insert
542255	Large fragment std. torque handle insert

Pangea Large Fragment Core Tray - ORIF instrument insert

Ref #	Description
542058	Fixed angle sleeve, T20
542059	Drill sleeve insert, T20, 60mm, ø3.2mm
542060	Drill sleeve insert, T20, 60mm, ø4.3mm
542061	K-wire sleeve insert, T20, ø2.0mm
542062	Variable angle drill guide, T20, ø3.2mm
542063	Variable angle drill guide, T20, ø4.3mm
542064	Compression drill guide, T20, ø3.2mm
542065	Lag screw drill guide, T20, ø4.5mm / ø3.2mm
542066	Screwdriver bit, AO, T20, 93mm
542067	Screwdriver bit, AO, T20, 180mm
542099	Delta handle, large, AO, T20
702812	Countersink, AO, ø4.5/6.5mm
705014	Depth gauge, T20, 0-120mm
542105	Threaded guide post, T20
706416	Soft tissue elevator, T20
705019	Temporary plate fixator, AO
542069	Screw capture sleeve, T20
542050	Drill bit, AO, Ø3.2mm x 145mm
542051	Drill bit, AO, Ø3.2mm x 215mm
542052	Drill bit, AO, Ø4.3mm x 145mm
542053	Drill bit, AO, Ø4.3mm x 215mm
542054	Drill bit, AO, Ø4.5mm x 135mm
542024	Tap, locking, AO, Ø4.0mm x 145mm
542057	Tap, locking, AO, Ø5.0mm x 145mm
705054	Tap, Cancellous, Ø6mm x 180mm
702808	Tap, AO, Ø4.5 x 145mm
390192	K-wire, Ø2.0 x 150mm
705002	K-wire drill tip, Ø2.0 x 234mm



System components

Pangea Small Fragment Core Tray

Top level consisting of the instruments listed below:

Torque limiter insert instruments


Ref #	Description
542068	Torque limiting T-handle, AO, T20 6Nm

System components

Pangea Small Fragment Core Tray - Instruments

Third level consisting of the instruments listed below:

This level can be removed and replaced with the optional inserts: Core tray large fragment reduction insert or core tray Asnis III 4.0mm cannulated screw insert.



Ref #	Description
542207	Core tray auxiliary insert
542208	Silicone mat, auxiliary insert

Pangea Small Fragment Core Tray - Optional large fragment reduction insert

Ref #	Description
542249	Core tray large reduction insert


Instruments

Ref #	Description
700151	Hook
700153	Ball spike
705294	Periosteal elevator, round edge 6mm
705295	Periosteal elevator, flat blade 13mm
702927	Repositioning forceps, L205mm
702940	Reduction forceps with serrated jaws
390084	Reduction pin, AO, Ø5.0 x 180mm
700367	T-handle, AO quick-chuck

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



Ref #	Description
542213	Large fragment screw rack (with lid)
542214	Large fragment screw rack lid

4.0mm locking screw self-tapping, T20 drive


Ref #	Length (mm)	Ref #	Length (mm)
662214	14	662242	42
662216	16	662244	44
662218	18	662246	46
662220	20	662248	48
662222	22	662250	50
662224	24	662255	55
662226	26	662260	60
662228	28	662265	65
662230	30	662270	70
662232	32	662275	75
662234	34	662280	80
662236	36	662285	85
662238	38	662290	90
662240	40	662295	95



System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



5.0mm locking screw self-tapping, T20 drive

Ref #	Length (mm)	Ref #	Length (mm)
662314	14	662348	48
662316	16	662350	50
662318	18	662355	55
662320	20	662360	60
662322	22	662365	65
662324	24	662370	70
662326	26	662375	75
662328	28	662380	80
662330	30	662385	85
662332	32	662390	90
662334	34	662395	95
662336	36	662400S*	100
662338	38	662405S*	105
662340	40	662410S*	110
662342	42	662415S*	115
662344	44	662420S*	120
662346	46		



5.0mm periprosthetic locking screw self-tapping, T20 drive

Ref #	Length (mm)	Ref #	Length (mm)
661210	10	661216	16
661212	12	661218	18
661214	14	661220	20




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



4.5mm cortex screw self-tapping, T20 drive

Ref #	Length (mm)	Ref #	Length (mm)
661714	14	661760	60
661716	16	661765	65
661718	18	661770	70
661720	20	661775	75
661722	22	661780	80
661724	24	661785	85
661726	26	661790	90
661728	28	661795	95
661730	30	661800S*	100
661732	32	661805S*	105
661734	34	661810S*	110
661736	36	661815S*	115
661738	38	661820S*	120
661740	40	661825S*	125
661742	42	661830S*	130
661744	44	661835S*	135
661746	46	661840S*	140
661748	48	661845S*	145
661750	50	661850S*	150
661755	55		



Plate attachments

Ref #	Description
662202S*	5.0mm cable plug
663201	Washer, T20




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



6.0mm cancellous screw full thread, T20 drive

Ref #	Length (mm)	Ref #	Length (mm)
608020	20	608090	90
608025	25	608095	95
608030	30	608100S*	100
608035	35	608105S*	105
608040	40	608110S*	110
608045	45	608115S*	115
608050	50	608120S*	120
608055	55	608125S*	125
608060	60	608130S*	130
608065	65	608135S*	135
608070	70	608140S*	140
608075	75	608145S*	145
608080	80	608150S*	150
608085	85		




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



6.0mm cancellous screw, 16mm thread T20 drive

Ref #	Length (mm)	Ref #	Length (mm)
608230	30	608295	95
608235	35	608300S*	100
608240	40	608305S*	105
608245	45	608310S*	110
608250	50	608315S*	115
608255	55	608320S*	120
608260	60	608325S*	125
608265	65	608330S*	130
608270	70	608335S*	135
608275	75	608340S*	140
608280	80	608345S*	145
608285	85	608350S*	150
608290	90		




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



6.0mm cancellous screw, 32mm thread, T20 drive

Ref #	Length (mm)	Ref #	Length (mm)
608445	45	608500S*	100
608450	50	608505S*	105
608455	55	608510S*	110
608460	60	608515S*	115
608465	65	608520S*	120
608470	70	608525S*	125
608475	75	608530S*	130
608480	80	608535S*	135
608485	85	608540S*	140
608490	90	608545S*	145
608495	95	608550S*	150



*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

1. Schmidt, W et al. "Stryker Orthopaedic Modeling and Analytics (SOMA): A Review." *Surgical Technology International*, vol. 32 (2018): 315-324.
2. Liporace, F, and Yoon, R. "Nail Plate Combination Technique for Native and Periprosthetic Distal Femur Fractures." *Journal of Orthopaedic Trauma*, vol. 33,2 (2019): 64-68. doi:10.1097/BOT.0000000000001332

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The instructions for use, operative techniques, cleaning instructions, patient information leaflets and other associated labeling may be requested online at ifu.stryker.com or stryker.com. If saving the instructions for use, operative techniques, cleaning instructions from the above mentioned websites, please make sure you always have the most up to date version prior to use.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, AxSOS 3, Dall-Miles, and Pangea. All other trademarks are trademarks of their respective owners or holders.



Manufacturer:
Stryker GmbH
Bohnackerweg 1
2545 Selzach, Switzerland
www.stryker.com