

PangeaTM

Fibula Plating System



Operative Technique

- Distal Lateral Fibula
- Distal Posterolateral Fibula
- Distal Posterior Fibula

Contents

| | |
|--|-----------|
| Indications and contraindications | 3 |
| Pangea introduction | 4 |
| Section 01: Pangea Overview | 5 |
| Trays | 6 |
| Color coding | 8 |
| Small and large fragment trays | 9 |
| Instrumentation | 11 |
| Section 02: System Overview Pangea Fibula Plating | |
| System Distal Fibula | 18 |
| Required trays | 19 |
| Distal Fibula Plate offering | 20 |
| Screw offering | 22 |
| Plate details | 24 |
| Section 03: Surgical Protocol Pangea Fibula Plating | |
| System Distal Fibula | 27 |
| Preoperative planning | 28 |
| Surgical approaches and instrumentation | 30 |
| Reduction and plate contouring..... | 31 |
| Plate insertion | 32 |
| Plate positioning..... | 33 |
| Provisional plate fixation | 35 |
| Non-locking screw insertion..... | 36 |
| Locking screw insertion..... | 37 |
| Lag screw technique..... | 38 |
| Syndesmosis application | 39 |
| Final constructs | 40 |
| Section 04: System Components | 43 |

This publication sets forth detailed recommended procedures for using Stryker devices and instruments. It offers guidance that you should heed, but, as with any such technical guide, each surgeon must consider the particular needs of each patient and make appropriate adjustments when and as required.

This document is applicable to US and Canada.

The surgeon must advise patients of surgical risks and make them aware of adverse effects and alternative treatments.

WARNING

The patient should be advised that the device cannot and does not replicate a normal healthy bone, that the device can break or become damaged as a result of strenuous activity or trauma and that the device has a finite expected service life.

Removal or revision of the device may be required sometime in the future due to medical reasons.



Pangea Platform indications

The Pangea Platform is indicated for the internal fixation and stabilization of bone fractures, osteotomies, and arthrodesis in normal and osteopenic bone, including:

- Diaphyseal, metaphyseal, epiphyseal, extra- and intra-articular fractures
- Non-unions, malunions, and deformities
- Periprosthetic fractures

The Pangea Platform is also indicated for children (2-12 years) and adolescents (12 – 21 years) for the internal fixation and stabilization of bone fractures of the diaphysis and metaphysis in which growth plates have fused or in which growth plates will not be crossed by implants.

Pangea Fibula Plating System indications

The Pangea Fibula Plating System is indicated for the internal fixation and stabilization of fibula bone fractures and osteotomies in normal and osteopenic bone, including:

- Diaphyseal, metaphyseal, epiphyseal, extra- and intra-articular fractures
- Non-unions, malunions, and deformities
- Periprosthetic fractures

Compatibility with other systems

Components from the Pangea Fibula Plating System may be used with the following systems:

- AxSOS 3
- Stryker Plating System (SPS)
- Dall-Miles cable system
- Pangea Platform
- Gravity Synchfix

Please remember that the compatibility of different product systems has not been tested unless specified otherwise in the product labeling. Consult instructions for use (www.ifu.stryker.com)

for a complete list of potential adverse effects, contraindications, warnings and precautions.

Contraindications

The physician's education, training, and professional judgement must be relied upon to choose the most appropriate device and treatment.

Conditions presenting an increased risk of failure include:

- Any active or suspected latent infection or marked local inflammation in or about the affected area
- Compromised vascularity that would inhibit adequate blood supply to the fracture or the operative site
- Bone stock compromised by disease, infection or prior implantation that cannot provide adequate support and / or fixation of the devices
- Material sensitivity, documented or suspected
- Patients having inadequate tissue coverage over the operative site
- Implant utilization that would interfere with anatomical structures or physiological performance
- Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care
- Other medical or surgical conditions which would preclude the potential benefit of surgery

WARNING

Follow the instructions provided in our cleaning and sterilization guide (OT-RG-1). All non-sterile devices must be cleaned and sterilized before use.


Multicomponent instruments must be disassembled for cleaning. Please refer to the corresponding assembly / disassembly instructions.

WARNING

Stryker bone screws are not approved or intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

Pangea overview

Pangea systems provide small and large fragment plating solutions for fracture treatment. Pangea's plate offerings include both utility and anatomical plates to address various fracture patterns and anatomy. The plates and non-locking screws are produced from titanium alloy (Ti6Al4V ELI), whereas the locking screws are produced from cobalt-chrome alloy (CoCr).




CoCr locking screw and Ti6Al4V ELI plate hole

Variable angle locking technology

Pangea's variable angle locking technology uses a CoCr locking screw, which is harder than the Ti6Al4V ELI plate, allowing for the screwhead's threads to form a definitive locking position in the plate's locking hole by engaging the softer, Ti6Al4V ELI material.

This technology allows the user to aim and lock the screw into the plate within a true 30° cone of the predetermined hole trajectory. The variable angle drill guide provided with the system offers guidance with respect to the limit of the 30° cone. The locking mechanism remains functional for up to three attempts at locking screw insertion.




Universal holes offer 30° cone of angulation

Hybrid LC Holes (locking/compression)

Hybrid LC Holes allow for either active compression with the use of a non-locking screw in the compression section of the hole or variable angle locking with the use of a locking screw in the universal section of the hole. If locking is not desired, the universal section of the hole also accepts non-locking screws. Each Hybrid LC Hole is designed to provide up to 2mm of compression.

Note: Hybrid LC Holes are not available with every plate type. Refer to the "Plate details" page for additional information.



A: Universal: For locking or non-locking screws
B: Compression: For non-locking screws only

Hybrid LC Hole

Section 01


Pangea Overview

Pangea overview


The Pangea systems require the use of a small or large fragment core tray, which contain the necessary instruments and screws to be used for every case. Plates are contained in separate anatomic plate trays, utility plate trays, or optional tray inserts. Some plates are offered sterile packaged only.

Anatomic and utility plate trays contain plates that will correspond to each color code within the small and large fragment core trays. To complete a case, a surgeon will often need an implant tray and its corresponding core tray.

| | Small fragment | | Large fragment | |
|------------------|----------------|-------------|-------------------|-------|
| Color coding | Purple | Yellow | Orange | Blue |
| Screw diameters | 2.7mm | 3.5 / 4.0mm | 4.0 / 4.5 / 6.0mm | 5.0mm |
| Screwdriver type | T8 | T15 | T20 | T20 |



Pangea Small Fragment Core Tray



Pangea Large Fragment Core Tray



Auxiliary insert and optional inserts


The Pangea small and large fragment core trays are designed to offer users modularity in their set configurations. Each core tray's standard configuration includes an auxiliary insert containing a silicone mat for storage of miscellaneous instrumentation. Listed below are optional inserts that can be ordered separately and placed into the core tray to meet the user's needs.

When using an optional insert with the Pangea small or large fragment core tray, the auxiliary insert may be replaced with one of the optional inserts.


| | Auxiliary insert with silicone mat | Small fragment reduction insert | Large fragment reduction insert | Small fragment standard plate insert | Small fragment ankle plate insert | Asnis III 4.0mm cannulated screw insert |
|--------------------------|------------------------------------|---------------------------------|---------------------------------|--------------------------------------|-----------------------------------|---|
| Small fragment core tray | ✓ | ✓ | | ✓ | ✓ | ✓ |
| Large fragment core tray | ✓ | | ✓ | | | ✓ |

Optional insert configurations


Small fragment utility plate insert



Small fragment ankle plate insert



Asnis III 4.0mm cannulated screw insert
















Note: Optional inserts may not be available in all markets. Check with your local Stryker sales representative for availability.

Color coding





The Pangea systems are color-coded to allow the user to easily identify the proper instrumentation for a particular plate type or screw type. Each color represents the proper drill, drill guide, or screwdriver for a particular screw diameter. The small fragment is color-coded purple for 2.7mm screws and yellow for 3.5mm and 4.0mm screws. The large fragment is color-coded orange for 4.0/4.5/6.0mm screws and blue for 5.0mm screws.

The small and large fragment core trays offer short and long drill bits to account for various depths of bone stock. These ORIF pilot drills are calibrated for the surgeon to measure depth by referencing the associated drill guide. ORIF pilot drill bits are identified by one colored stripe, while lag screw overdrills have one colored stripe and one black stripe. The drill bit's diameter can be found on the AO quick connect.

| |  Small fragment | |  Large fragment | |
|--------------------------------------|---|---|---|---|
| Screw type | 2.7mm cortex 2.7mm locking | 3.5mm cortex 3.5mm locking 4.0mm cancellous | 4.0mm locking 4.5mm cortex 6.0mm cancellous | 5.0mm locking |
| Color code | Purple | Yellow | Orange | Blue |
| ORIF pilot drills |  |  |  |  |
| Lag screw overdrills |  |  |  | N/A |
| ORIF pilot drill lengths | Short: 135mm Long: 175mm | Short: 135mm Long: 215mm | Short: 145mm Long: 215mm | Short: 145mm Long: 215mm |
| ORIF pilot drill calibrations | Short: 0-40mm Long: 0-80mm | Short: 0-40mm Long: 0-120mm | Short: 0-50mm Long: 0-120mm | Short: 0-50mm Long: 0-120mm |
| Drill guides |  |  |  |  |


Small and large fragment trays

Small fragment core tray set content – 2.7 / 3.5 / 4.0mm screws

| T8 torx interface for all 2.7mm screws | T15 torx interface for all 3.5mm screws and 4.0mm cancellous screws |
|---|--|
|  |  |
|  <p>Ø2.7mm locking Ø2.7mm cortex</p>  |  <p>Ø3.5mm locking Ø3.5mm cortex Ø4.0mm cancellous Partially threaded (Variable thread length) Ø4.0mm cancellous Fully threaded</p>  |
| 2.7mm lag screw drill guide | 3.5mm lag screw drill guide |
| 2.0mm pilot drill for 2.7mm screws 2.7mm overdrill for 2.7mm lag screws | 2.5mm pilot drills for 3.5 / 4.0mm screws 3.5mm overdrill for 3.5mm lag screws |
| Torque limiting handle, 1.7Nm for 2.7mm locking screws | Torque limiting handle, 4.0Nm for 3.5mm locking screws |
| Depth gauge for T8 2.7mm screws | Depth gauge for T15 3.5 / 4.0mm screws |
| Variable angle drill guides | |
| Compression drill guide | |
| Fixed angle guides with drill sleeve inserts | |
| Small and large Delta screwdriver handles with AO quick connect | |
| Calibrated soft tissue elevator with distal suture hole for alternative plate insertion techniques | |

Large fragment core tray set content – 4.0 / 4.5 / 5.0 / 6.0mm screws

T20 torx interface for all large fragment screws



4.5mm lag screw drill guide

3.2mm pilot drills for 4.0 / 4.5 / 6.0mm screws
4.3mm pilot drills for 5.0mm screws
4.5mm overdrill for 4.5mm lag screws

Torque limiting handle, 6.0Nm for 4.0mm and 5.0mm locking screws

Depth gauge for T20 4.0/4.5/5.0/6.0mm screws

Variable angle drill guides

Compression drill guide

Fixed angle guides with drill sleeve inserts

Large Delta screwdriver handles with AO quick connect

Calibrated, curved soft tissue elevator with distal suture hole for alternative plate insertion techniques

Fixed angle drill sleeve

First, the appropriate drill sleeve insert is loosely inserted into the fixed angle sleeve. Next, the entire assembly is inserted into the desired plate hole, which can be confirmed for proper placement with tactile feedback upon insertion. Lastly, turn the drill sleeve insert clockwise to tighten. This fastens the assembly securely to the plate (Fig. 1).

The drill sleeve insert is meant to be hand tightened only. The hex interface on the T8, T15, and T20 screwdriver bits facilitate removal of each respective drill sleeve insert (Fig. 3).




Fig. 1: Fixed angle drill sleeve and insert assembly

Threaded guide post

The T15 and T20 threaded guide posts are intended only for the threaded monoaxial hole in some Pangea plates and serve as the attachment point for MIS targeting capabilities (Fig. 2). Additionally, the threaded guide posts may be utilized as a plate insertion handle, joystick, K-wire sleeve, and drill guide for plates with a threaded monoaxial hole. The threaded guide post is meant to be hand tightened only (Fig. 2).

Note: Threaded monoaxial holes are not available with every plate type. Refer to the “Plate details” page for additional information.

The hex interface on the T20 screwdriver bit facilitates removal of the T20 threaded guide post. The T15 threaded guide post does not have a hex interface and may be removed using the torx portion of the T15 screwdriver bit (Fig. 3).




Fig. 2: Threaded guide post placed in threaded monoaxial hole




Fig. 3: Hex interface between screwdriver bit, fixed angle drill insert, and T20 threaded guide post

Variable angle drill guide

The ball-and-cone variable angle drill guides are used in combination with their respective drill bit to determine screw trajectory. The color-coded guide restricts the degree of screw angulation to 15° in any direction resulting in a 30° cone of the predetermined hole trajectory.


When using the ball end of the guide, gently press the guide into the plate's universal or Hybrid LC Hole. The ball end of the drill guide can be gently rotated in these holes using the handle while maintaining 15° of angulation. For small fragment indications, to ensure a precise 15° angulation, use the cone end of the variable angle drill guide by engaging the cone end of the guide into the plate hole.



Variable angle drill guide

Lag screw drill guide


Lag screw drill guides are offered to accommodate 2.7mm, 3.5mm, or 4.5mm lag screws. The color-coded guide features an overdrill sleeve with serrated tip for near cortex drilling and a self-centering pilot drill sleeve to complete the lag screw pilot hole. Lag screws may be placed through a plate or independently.



Lag screw drill guide

Compression drill guide







The color-coded compression drill guide is an eccentric drill guide used to apply compression in plates equipped with Hybrid LC Holes. Arrows are inscribed on the guide and must be aimed towards the fracture site when engaged into the plate's Hybrid LC Hole. This allows the user to properly drill the Hybrid LC Hole and achieve up to 2mm of compression upon insertion of a non-locking screw.







Compression drill guide

Handles

Pangea offers multiple handle options. These options include small and large Delta handles, and three torque limiting handles for locking screws. All handles are equipped with a small AO-coupling.

| Screwdriver type | Locking screw Ø | Torque | Torque limiting handles |
|--|-----------------|--------|--|
|  | 2.7mm | 1.7Nm |  |
|  | 3.5mm | 4.0Nm |  |
|  | 4.0mm 5.0mm | 6.0Nm |  |

| Screwdriver type | Non-locking screw Ø | Small and large Delta handles |
|---|---|---|
|  | 2.7mm 3.5mm 4.0mm 4.5mm 6.0mm |  |
|  | 2.7mm 3.5mm 4.0mm 4.5mm 6.0mm |  |

Depth measuring instruments

Pangea systems provide two options for determining screw length.

Drill calibrations

All pilot drills are calibrated to allow the surgeon to determine the appropriate screw length when drilling through either the fixed angle sleeve assembly, threaded guide post, or the ball end of the variable angle drill guide (Fig. 1). The cone end of the variable angle drill guide is not compatible with the calibrated drill bits to determine screw length. The calibrations when read against the selected drill guide measure the distance to the tip of the drill bit.



Fig. 1: Measuring screw length using calibrated drill bit

Depth gauge

A depth gauge can be used independently or through a plate hole. Depth gauges correspond with the screw head size i.e., T8, T15, or T20 (Fig. 2).



Fig. 2: T8, T15, and T20 depth gauges

Screw selection

To verify a screw's length, use the screw measurement scale found on the screw rack. The screw tip is placed with its tip against the back stop and its length can be read off the scale (Fig. 3).



WARNING

Always check the correct position and length of the inserted screws by fluoroscopy.




Fig. 3: Screw rack with screw measurement scale

Note: The T8 and T15 depth gauges are color-coded purple and yellow, respectively. The T20 depth gauge is not color-coded and is intended to be used for all T20 screws.

Screw capture sleeve

Pangea offers T8, T15, and T20 screw capture sleeves. The screw capture sleeves are optional devices available to offer efficient screw pick-up, insertion, and removal.

First, the T8, T15, or T20 long screwdriver shaft is inserted into its respective screw capture sleeve. Next, the screw is loaded onto the screwdriver shaft and the knob on the screw capture sleeve is depressed to capture the screw head. Once the screw is captured by the inner sleeve, the knob may be released and the screw is successfully captured. The screw may now be inserted into the pilot hole. Lastly, the screw is released by depressing the button on the end of the device prior to final tightening (Fig. 1).




Fig. 1: Screw capture sleeve assembly

Countersinks and Taps

Countersinks are available for all screws sizes to reduce screw head prominence when the screw is used independently from a plate (Fig. 2).

Taps are available for all screw types and diameters. All screws are self-tapping; however, when inserting a screw in hard bone, it is recommended to use the appropriate tap prior to screw insertion (Fig. 3).



Fig. 2: Countersink



Fig. 3: Tap

Temporary plate fixator

Temporary plate fixators may be used as a provisional plate fixation device and should be limited to the shaft holes of the plate. The temporary plate fixator functions by pushing the shaft of the plate to the bone. The temporary plate fixator is designed with a self-drilling, self-tapping tip for quick insertion into cortical bone (Fig. 4).




Fig. 4: Temporary plate fixator

Small fragment reduction instruments and retractors

These instruments facilitate fracture reduction and soft tissue management.



Periosteal elevator, round edge (705294)



Periosteal elevator, flat edge (705295)



Sharp hook (700151)



Straight reduction clamp (705297) Lobster claw (702932)



Ball spike (700153)



Hohmann retractor 6mm (700664)



Hohmann retractor 15mm (700667)

Large fragment reduction instruments and retractors

These instruments facilitate fracture reduction and soft tissue management.



T-handle large, AO coupling (700367)



Repositioning forceps, 205mm (702927)

Reduction forceps with serrated jaws (702940)




Reduction pin, 5.0mm (390083/390084)




Sharp hook (700151)



Ball spike (700153)



Periosteal elevator, round edge (705294)



Periosteal elevator, flat edge (705295)

Note: Large fragment reduction instruments are contained in an optional large fragment reduction insert tray and may not be available in all markets. Check with your local Stryker sales representative.

Bending instruments and plate cutters

These instruments are available to aid in plate bending and plate cutting. Please refer to the table for the bending and cutting specifics of each instrument.

⚠ WARNING

Repetitively bending or rebending the plates may lead to early fatigue failure or inability to lock at a screw hole. Bending the plate at the screw hole risks compromising the locking mechanism.

⚠ WARNING

In most cases the precontoured plate will fit without the need for further bending. However, should additional out of plane bending of the plate be required, refer to the following table for bending devices and respective plate sizes.

| Device Type | 2.7 plates | 3.5 plates | 5.0 plates |
|--------------------|------------|------------|------------|
| Plate bender | ✓ | | |
| Bending irons | ✓ | ✓ * | |
| Plate cutter | ✓ | | |
| Table plate bender | | ✓ | ✓ |

*Bending irons cannot be used to bend the following 3.5 plates: straight broad, extra articular distal humerus, proximal lateral tibia, or distal medial femur




Fig. 3: Table Plate Bender




Fig. 2: Plate Cutter



Fig. 1: Bending Iron



Fig. 4: Plate Bender

Section 02

System Overview


Pangea Fibula Plating System

Distal Fibula

Required trays

ORIF distal fibula

This page details the trays required to use the Pangea Distal Fibula Plates.



Pangea Small Fragment Core Tray



Pangea Distal Fibula Plate Tray



**Distal Lateral
Fibula Plate**



**Distal
Posterolateral
Fibula Plate**



**Distal
Posterior
Fibula Plate**



Hook Plate

Distal fibula plate offering

Implants - Plates

Distal Lateral Fibula Plates



| 81mm | 95mm | 109mm | 123mm | 137mm | 151mm | 165mm | 179mm | 207mm | 235mm |
|---------|---------|---------|---------|---------|---------|----------|-----------|-----------|-----------|
| 4 holes | 5 holes | 6 holes | 7 holes | 8 holes | 9 holes | 10 holes | 11 holes* | 13 holes* | 15 holes* |



Distal Posterolateral Fibula Plates



| 81mm | 95mm | 109mm | 123mm | 137mm | 151mm | 165mm | 179mm | 207mm | 235mm |
|---------|---------|---------|---------|---------|---------|----------|-----------|-----------|-----------|
| 4 holes | 5 holes | 6 holes | 7 holes | 8 holes | 9 holes | 10 holes | 11 holes* | 13 holes* | 15 holes* |



* Available sterile packed only

Distal fibula plate offering

Implants - Plates

Distal Posterior Fibula Plates

T8



79mm
6 holes



94mm
8 holes



109mm
10 holes



125mm
12 holes



147mm
15 holes



170mm
18 holes*



201mm
22 holes*



231mm
26 holes*




* Available sterile packed only

Implants - Screws

The following screws are suitable for use with the Pangea Fibula Plating System:



| Description | Ø2.7mm cortex | Ø2.7mm locking |
|-----------------------|--|---|
| Length (increments) | 8-80mm 8-50mm (2mm) 50- 80mm (5mm) | 8-80mm 8-50mm (2mm) 50- 80mm (5mm) |
| Material | Ti6Al4V ELI | Cobalt-chrome (CoCr) alloy |
| Locking torque | N/A | 1.7Nm |
| Washer | Washer T8 (40-30900) | N/A |
| Screwdriver interface |  Screwdriver bit, AO, T8 short, 93mm (542011) Screwdriver bit, AO, T8 long, 180mm (542012) | |
| Angulation | +/-15° (30° cone) | |
| Pilot drill Ø | Ø2.0mm | |
| Drill bits | Ø2mm x 135mm (542000) Ø2mm x 175mm (542001) | |
| Taps | Tap, locking, Ø2.7mm x 125mm (542003) Tap, cortex, Ø2.7mm x 125mm (702801) | |

Note: Purple represents the color code for the T8 screws. Ensure the drill guides, drills, screwdrivers, and depth gauges correspond to the purple color code.

Screw offering

T15 Implants - Screws

The following screws are suitable for use with the Pangea Fibula Plating System:



| Description | Ø3.5mm cortex | Ø4.0mm cancellous Partially threaded (Variable thread length) | Ø4.0mm cancellous Fully threaded | Ø3.5mm locking |
|------------------------|--|---|---|---|
| Length (increments) | 10-120mm 10-50mm (2mm) 50-120mm (5mm) | 10-100mm 10-50mm (2mm) 50-100mm (5mm) | 10-100mm 10-50mm (2mm) 50-100mm (5mm) | 10-120mm 10-50mm (2mm) 50-120mm (5mm) |
| Material | Titanium alloy (Ti6Al4V ELI) | | | Cobalt-chrome alloy (CoCr) |
| Locking torque | N/A | | | 4Nm |
| Washer | Washer T15 (663001) | | | N/A |
| Screwdriver interface |  Screwdriver bit, AO, T15 short, 93mm (542031) Screwdriver bit, AO, T15 long, 180mm (542032) | | | |
| Angulation | +/-15° (30° cone) | | | |
| Pilot drill Ø | Ø2.5mm | | | |
| Drill bits | Ø2.5mm x 135mm (542020) Ø2.5mm x 215mm (542021) | | | |
| Taps | Tap, locking, Ø3.5mm x 125mm (542023) Tap, cortex, Ø3.5mm x 125mm (702802) Tap, cancellous, Ø4mm x 125mm (702803) | | | |

Note: Yellow represents the color code for the T15 screws. Ensure the drill guides, drills, screwdrivers, and depth gauges correspond to the yellow color code.

Distal Lateral Fibula Plate

- 7 distal T8 universal holes
- 3 T15 syndesmosis screw holes compatible with non-locking screws or flexible fixation i.e., Gravity Synchfix button; allows for syndesmosis fixation when indicated
- Shaft holes are all T15 universal holes
- 1.6mm K-wire holes
- Rounded and tapered ends designed for insertion under soft tissue

Screw options

- T15** T15 – 3.5mm locking and cortex screws and 4.0mm cancellous screws.
- T8** T8 - 2.7mm locking and cortex screws.


Variable angle locking technology

Universal holes that accept both locking and non-locking screws within a 30° cone.

SOMA designed

These plates are designed using Stryker Orthopedics Modeling and Analytics (SOMA) software.¹

- Evidence-based screw hole placement designed to offer a wide range of trajectory options for particular anatomy¹
- Anatomical fit, creating contoured plates with left and right specific options¹



Distal Posterolateral Fibula Plate

- 5 distal T8 universal holes
- Low profile metaphyseal cluster to minimize perineal tendon irritation
- 3 T15 syndesmosis screw holes compatible with non-locking screws or flexible fixation i.e., Gravity Synchfix button; allows for syndesmosis fixation when indicated
- Shaft holes are all T15 universal holes
- 1.6mm K-wire holes
- Designed for buttressing oblique or Supination External Rotation ("SER") ankle fractures
- Rounded and tapered ends designed for insertion under soft tissue

Screw options

T15 T15 – 3.5mm locking and cortex screws and 4.0mm cancellous screws.

T8 T8 - 2.7mm locking and cortex screws.


Variable angle locking technology

Universal holes that accept both locking and non-locking screws within a 30° cone.

SOMA designed

These plates are designed using Stryker Orthopedics Modeling and Analytics (SOMA) software.¹

- Evidence-based screw hole placement designed to offer a wide range of trajectory options for particular anatomy¹
- Anatomical fit, creating contoured plates with left and right specific options¹



Distal Posterior Fibula Plate

- 6 distal T8 universal holes
- Low profile metaphyseal cluster to minimize perineal tendon irritation
- Shaft holes are all T8 universal holes
- 1.6mm K-wire holes
- Designed for direct posterior placement through the posterolateral approach
- Designed to allow independent syndesmotic fixation when necessary
- Rounded and tapered ends designed for insertion under soft tissue

Screw options

-  T8 - 2.7mm locking and cortex screws.


Variable angle locking technology

Universal holes that accept both locking and non-locking screws within a 30° cone.

SOMA designed

These plates are designed using Stryker Orthopedics Modeling and Analytics (SOMA) software.¹

- Evidence-based screw hole placement designed to offer a wide range of trajectory options for particular anatomy¹
- Anatomical fit, creating contoured plates with left and right specific options¹



Section 03

Surgical Protocol


Pangea Fibula Plating System

Distal Fibula

Plate sizes and laser markings

Selecting the appropriate plate length involves considerations including bone quality, fracture configuration and location, and the type of bone healing expected.


The plate's reference number and left/right orientation is written on the distal end of the plate. Additionally, the number of holes, reference number, Left/Right orientation, and plate length is written on the bottom of the implant tray.





Preoperative planning

Digital templates

The use of digital templates in association with adequate X-ray/ fluoroscopy or CT scans may assist in the selection of an appropriately sized implant. Pangea digital templates are licensed to template service providers e.g., PeekMed, Sectra AB, Lexi, and Merge Healthcare. Template service providers provide software tools and access to 3D models in order to allow surgeons to perform preoperative planning.




Preoperative planning using 3D planning software




Scale marker

Details: Scale, marker, and product description



Rotational marker

Example of 2D femur templates using ".dxf" files



3D scanned bone model

".stl" plate model

".igs" axis data

".stl" screw model

Example of 3D femur templates using ".stl" and ".igs" files

NOTICE If digital X-ray images are used, correct magnification has to be verified prior to use.

Note: Digital templates may also be accessed by contacting your local Stryker sales representative.

Surgical approaches

The approaches described below are common approaches for the insertion of the selected distal fibula plate. These plates can be inserted through a variety of approaches based on surgeon preference, soft tissue, and fracture characteristics.

Distal fibula

Surgeons may use a standard lateral or posterolateral approach to the distal fibula (Fig. 1, 2).

Note: Depending on the surgical approach, care must be taken not to damage the superficial peroneal nerve, sural nerve, or the superior peroneal retinaculum.

Periosteal elevators

The soft tissue elevators can be utilized to separate or lift the soft tissues from the bone (Fig. 3). The curved soft tissue elevator is designed to create a pathway for the plate and is calibrated allowing the surgeon to monitor depth or measure plate length (Fig. 4).

| Part # | Instrumentation |
|--------|---|
| 542035 | Curved soft tissue elevator, T15 |
| 705294 | Periosteal elevator, round edge 6mm |
| 705295 | Periosteal elevator, straight edge 13mm |




Fig. 1: Lateral approach to the distal fibula




Fig. 2: Posterolateral approach to the distal fibula




Fig. 3: Periosteal elevators




Fig. 4: Curved soft tissue elevator, T15

Reduction

Indirect fracture fragment reduction can be achieved through traction, manipulation, external fixation, or the surgeon's preferred technique. Direct fragment reduction can be achieved using elevators, clamps, or K-wires. Anatomic reduction of the fracture can be achieved through direct visualization and use of clamps.

The use of K-wires may be helpful in maintaining provisional reduction (Fig. 1). Independent lag screws may definitively stabilize articular reduction prior to plate insertion.




Fig. 1: Reduction of distal fibula fracture

Plate contouring

Should bending of the plate be required, the bending irons (ref 703938) or the table plate bender (ref 702900) may be used. The bending irons are designed to be used as a pair. The slots allow the device to slide over the shaft of the plate for ease of bending (Fig. 2).



Fig. 2: Bending irons

Plate insertion

When a plate insertion handle is desired, the surgeon may use either the T8 or T15 fixed angle sleeve assembly. The T8 or T15 fixed angle sleeve assembly may be attached to any of the plate's respective universal holes (Fig. 1).

Using the T8 or T15 fixed angle sleeve assembly, the surgeon may insert the selected distal fibula plate submuscularly (Fig. 2)

Once the plate has been inserted, the surgeon may use either plate insertion handle as a joystick to adjust plate positioning (Fig. 3).

Assembly instructions of the fixed angle sleeve and drill sleeve insert can be found on page 11.

CAUTION

Avoid plate insertion through the muscle to prevent intramuscular vessel disruption. Minimize periosteal disruption while inserting the plate to help preserve bone blood supply.



Fig. 1: T8 or T15 fixed angle sleeve assembly used as a plate insertion handle or joystick



Fig 2: Plate insertion using the fixed angle sleeve as a plate insertion handle





Fig. 3: Fixed angle sleeve assembly used as a joystick for final plate adjustment

Distal Lateral Fibula Plate Positioning

When appropriate reduction of the fragments is achieved and verified by fluoroscopy, the distal portion of the plate should be placed so the 2.7mm cluster is centered over the lateral aspect of the distal fibula. Direct visualization or lateral fluoroscopic imaging should be utilized to confirm the plate is not placed over the edge of the distal fibula. When positioned correctly, the distal cluster should cup the lateral malleolus.


The surgeon should confirm that the proximal aspect of the plate is centered on the lateral aspect of the fibular shaft through either direct visualization or lateral fluoroscopic imaging.



Distal Posterolateral Fibula Plate Positioning

When appropriate reduction of the fragments is achieved and verified by fluoroscopy, the distal portion of the plate should be placed so the 2.7mm cluster is centered over the fibular ridge, just posterior to the lateral surface. Direct visualization or lateral fluoroscopic imaging should be utilized to confirm there is no overhang over the edge of the distal fibula. When positioned correctly, the distal aspect of the plate should not contact the peroneal tendons.


The surgeon should confirm either through direct visualization or lateral fluoroscopic imaging that the proximal aspect of the plate is centered on the lateral aspect of the shaft of the fibula.



Distal Posterior Fibula Plate Positioning

When appropriate reduction of the fragments is achieved and verified by fluoroscopy, the distal portion of the plate should be placed so the 2.7mm cluster is centered over the posterior metaphyseal fibula, just proximal to the peroneal groove. Direct visualization or lateral fluoroscopic imaging should be utilized to confirm there is no overhang over the edge of the distal fibula. When positioned correctly, the distal aspect of the plate should not contact the peroneal tendons.

The surgeon should confirm either through direct visualization or lateral fluoroscopic imaging that the proximal aspect of the plate is centered on the posterior aspect of the shaft of the fibula.



Provisional plate fixation

Provisional plate fixation may be achieved in various ways:

- Ø1.6mm Olive K-wire may be placed in any of the universal holes
- Ø1.6mm K-wire placed through the K-wire holes
- Ø2.0mm K-wire placed through the T8 or T15 fixed angle sleeve and drill sleeve insert

| Part # | Instrumentation |
|------------------|---|
| 542036 | Olive K-wire Ø1.6mm, T8/T15 |
| 390164 | 1.6mm K-wire |
| 390192 | 2.0mm x 150mm K-wire |
| 705002 | 2.0mm x 234mm K-wire, drill tip |
| 542005 542025 | Fixed angle sleeve, T8 Fixed angle sleeve, T15 |
| 542006 542026 | Drill sleeve insert, T8, 60mm, Ø2.0mm Drill sleeve insert, T15, 60mm, Ø2.5mm |

CAUTION

Be careful when using sharp instruments such as drills, taps, K-wires and temporary plate fixators or when inserting screws to avoid damage to the soft tissue or vessels by going too far past the far cortex.



Fig. 1: Provisional fixation of Distal Lateral Fibula Plate

Non-locking screw insertion

Insertion of a non-locking screw is started with the insertion of the appropriate drill guide for the screw hole.

To achieve the predetermined trajectory of the universal holes, use the T8 or T15 fixed angle sleeve with its respective drill sleeve insert (Fig. 1). Assembly instructions of the fixed angle sleeve and drill sleeve insert can be found on page 11. If a variable angle trajectory is desired, use the T8 or T15 variable angle drill guide (Fig. 2).

Next, using the 2.0mm drill bit for 2.7mm cortex screws or 2.5mm drill bit for 3.5mm cortex and 4.0mm cancellous screws, create a pilot hole by drilling through the selected drill guide (Fig. 3). The drill trajectory may be verified under fluoroscopy if required.

The depth may be measured utilizing the associated T8 or T15 depth gauge or the drill bit calibrations.

The selected screw is then inserted into the pilot hole using the T8 or T15 screwdriver bit (Fig. 4).

The T8 or T15 screw capture sleeve may be used to aid in retention between the screw and screwdriver shaft during screw insertion.

CAUTION

Use bi-cortical fixation when possible.

CAUTION

If excessive resistance is felt during insertion, or if the bone is dense, it is recommended to use a tap.

CAUTION

It is recommended to insert the screws by hand. If power tools are used, use those at low speed to avoid improper alignment.



Fig. 1: T8 and T15 fixed angle sleeve assembly for predetermined screw trajectory




Fig. 2: Variable angle drill guide for variable angle drilling



Fig. 3: Drill pilot hole with 2.0mm or 2.5mm drill bit



Fig. 4: Screw placement with appropriate screwdriver

Locking screw insertion

Locking screws can be placed within a 30° cone in any universal screw hole.

Uni-cortical fixation is recommended when bi-cortical fixation cannot be achieved safely. For example, when bi-cortical fixation will cause screw penetration into the articular surface. In situations of uni-cortical screws, an increased screw count may be needed to obtain sufficient fixation.

To achieve the predetermined trajectory of the universal holes, use the T8 or T15 fixed angle sleeve with its respective drill sleeve insert. Assembly instructions of the fixed angle sleeve and drill sleeve insert can be found on page 11. If a variable angle trajectory is desired, use the T8 or T15 variable angle drill guide (Fig. 2).

Next, using the 2.0mm drill bit for 2.7mm locking screws or 2.5mm drill bit for 3.5mm locking screws, create a pilot hole by drilling through the selected drill guide (Fig. 3). The drill trajectory may be verified with fluoroscopy if required.

The depth may be measured utilizing the associated T8 or T15 depth gauge or the drill bit calibrations.

The selected screw is then inserted into the pilot hole using the T8 or T15 screwdriver bit.

A T8 or T15 screw capture sleeve may be used to aid in retention between the screw and screwdriver shaft during screw insertion.

Use the appropriate torque limiting handle and respective T8 or T15 screwdriver bit to ensure proper seating of the locking screw. The torque limiter will produce an audible “click” when the required torque is achieved (Fig. 4).

⚠ WARNING

Always perform final tightening by hand using the appropriate torque limiter and screwdriver bit as final tightening with the power tool can cause over-torquing or damage to the screw-plate interface, which can lead to breaking or stripping screw heads.



Fig. 1: Fixed angle sleeve for predetermined trajectory drilling




Fig. 2: Variable angle drill guide for variable angle drilling



Fig. 3: Drill pilot hole with 2.0mm or 2.5mm drill bit



Fig. 4: Screw placement with appropriate torque limiter

Lag screw technique

Independent lag screw

To insert an independent lag screw, select the appropriate lag screw drill guide for the desired lag screw size. Then use the serrated end of the guide and the appropriate overdrill to drill the near cortex (Step 1a).

Next, insert the self-centering end of the lag screw drill guide into the gliding hole and use the appropriate pilot drill to drill through the far cortex (Step 2a).


Measure the depth of the hole using the pilot drill or appropriate depth gauge and insert the selected cortex screw. If required, countersink the gliding hole or apply a washer to the screw. Upon screw insertion, this technique will serve to lag the far cortex towards the near cortex, thus applying compression (Step 3a).

Lag screw through plate

To insert a lag screw through a plate hole, use the serrated end of the appropriate lag screw drill guide and its respective overdrill to drill the near cortex through the plate hole (Step 1b).

Next, drill the far cortex with the appropriate size pilot drill by placing either the variable angle drill guide into the plate hole or the self-centering end of the lag screw drill guide into the gliding hole through the plate. When drilling the far cortex using a variable angle drill guide, ensure the trajectory of the pilot drill is co-linear with the gliding hole (Step 2b).

Measure the depth of the hole using the pilot drill or appropriate depth gauge and insert the selected cortex screw. Upon screw insertion, this technique will serve to lag the far cortex towards the near cortex, thus applying compression through the plate. Screw holes in the plate may be populated to complete the construct (Step 3b).



Syndesmosis fixation

After completion of all other fixation, an intra-operative stress test is performed under fluoroscopy to assess syndesmotic stability.

If syndesmotic fixation is necessary, the Pangea Distal Lateral and Posterolateral Fibula Plates are designed with three non-locking syndesmotic holes to facilitate anatomic reduction of the syndesmosis through flexible fixation or non-locking screw fixation (Fig. 1). To achieve optimal anatomic reduction of the syndesmosis, direct visualization of the reduction is recommended. Screw trajectory should be orientated parallel to the tibial plafond in the coronal plane. Doing so will minimize the likelihood of mal-reduction (Fig. 2).




Fig. 1 T15 Non-locking syndesmotic holes in Distal Lateral Fibula Plate and Distal Posterolateral Fibula Plate




Fig. 2: Suggested syndesmotic screw trajectory

Flexible fixation

Syndesmotic holes may accept flexible fixation. Please refer to the surgical technique for the Gravity Synchfix for more information. Content ID 015664A_02-Feb-2017



GRAVITY™ Syndesmosis LP

| Ref # | Product | Contents |
|----------|--|---|
| 86SYN005 | GRAVITY SYNCHFIX Syndesmotic Fixation Device Sterile Packed Implant/ Instrument Kit | 2.8mm Guide Wire Countersink, Lateral Knot Capsule #5 Suture Construct Blunt Guide Wire Medical Button Holder Medial Button (Qty. 2) |

Screw fixation

A clamp or direct pressure may be used to obtain an anatomic syndesmotic reduction. The reduction may be maintained with K-wires when necessary. Reduction should be confirmed fluoroscopically on the A/P, mortise, and lateral view.

Using the T15 variable angle drill guide and long 2.5mm drill bit, create a pilot hole parallel to the tibio-talar articulation and the trans-malleolar axis. Care must be taken to avoid excessive angulation through the plate beyond the 30° cone. One or more tri-cortical or tetra-cortical T15 3.5mm cortex or T15 4.0mm cancellous screws can be used to achieve syndesmotic fixation based on patient factors, injury factors, and surgeon preference.



Final fluoroscopic check

After final fixation of the plate with all screws, the ankle should be internally and externally rotated under continuous fluoroscopy; the true distance of all screws should be inspected in the mortise, AP, and lateral views to ensure none of the screws have penetrated the joint or are prominent.



After the procedure, check that all implants are positioned correctly using an image intensifier.

Implant removal

Removal of the Pangea Distal Fibula Plates is not required in general. The additional surgical trauma and the risks associated with additional anesthesia should be individually outweighed against the potential benefits for every patient.

In the case of implant removal, the scar of the previous incision is (partly) re-opened and the screws and the plate are successively removed.



The syndesmosis screw is recommended to be taken out in a timely manner.

NOTICE

In the extreme event of broken or stripped screws, the Stryker implant extraction set (literature number IES-ST-1) includes a variety of broken screw removal instruments.

Final construct examples

Distal Lateral Fibula Plate




Predetermined trajectory



Screws in the predetermined screw trajectory

Variable angle trajectory




Screw trajectories using variable angle locking to obtain the widest allowable screw trajectory

Distal Posterolateral Fibula Plate




Predetermined trajectory



Screws in the predetermined screw trajectory

Variable angle trajectory




Screw trajectories using variable angle locking to obtain the widest allowable screw trajectory

Final construct examples

Distal Posterior Fibula Plate


Predetermined trajectory


Screws in the predetermined screw trajectory

Variable angle trajectory


Screw trajectories using variable angle locking to obtain the widest allowable screw trajectory


Section 04

System Components

System components

Pangea Distal Fibula Plate Tray

The information in this section is not intended to be used for sales and/or promotional purposes. This information is solely intended to be used as a reference for clinical usage.



| Ref # | Description |
|--------|--------------------------|
| 542223 | Distal Fibula Plate Tray |

Distal Lateral Fibula Plates

| Left Ref # | Right Ref # | Plate Length | Shaft Holes |
|------------|-------------|--------------|-------------|
| 540644 | 540664 | 81 mm | 4 |
| 540645 | 540665 | 95 mm | 5 |
| 540646 | 540666 | 109 mm | 6 |
| 540647 | 540667 | 123 mm | 7 |
| 540648 | 540668 | 137 mm | 8 |
| 540649 | 540669 | 151 mm | 9 |
| 540650 | 540670 | 165 mm | 10 |
| 540651S* | 540671S* | 179 mm | 11 |
| 540653S* | 540673S* | 207 mm | 13 |
| 540655S* | 540675S* | 235 mm | 15 |



Distal Posterolateral Fibula Plates

| Left Ref # | Right Ref # | Plate Length | Shaft Holes |
|------------|-------------|--------------|-------------|
| 540704 | 540724 | 81 mm | 4 |
| 540705 | 540725 | 95 mm | 5 |
| 540706 | 540726 | 109 mm | 6 |
| 540707 | 540727 | 123 mm | 7 |
| 540708 | 540728 | 137 mm | 8 |
| 540709 | 540729 | 151 mm | 9 |
| 540710 | 540730 | 165 mm | 10 |
| 540711S* | 540731S* | 179 mm | 11 |
| 540713S* | 540733S* | 207 mm | 13 |
| 540715S* | 540735S* | 235 mm | 15 |



*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Distal Posterior Fibula Plates

| Left Ref # | Right Ref # | Plate Length | Shaft Holes |
|------------|-------------|--------------|-------------|
| 540746 | 540776 | 79 mm | 6 |
| 540748 | 540778 | 94 mm | 8 |
| 540750 | 540780 | 109 mm | 10 |
| 540752 | 540782 | 125 mm | 12 |
| 540755 | 540785 | 147 mm | 15 |
| 540758S* | 540788S* | 170 mm | 18 |
| 540762S* | 540792S* | 201 mm | 22 |
| 540766S* | 540796S* | 231 mm | 26 |



2.7 Hook Plate

| Ref # | Length (mm) | Holes |
|--------|-------------|-------|
| 541303 | 37 | 3 |
| 541305 | 57 | 5 |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray


Top level consisting of the instruments listed below:



| Ref # | Description |
|--------|--|
| 542201 | Aluminum lid - universal |
| 542203 | Small fragment core tray base |
| 542204 | Small fragment ORIF instrument insert |
| 542254 | Small fragment std. torque handle insert |

Pangea Small Fragment Core Tray - ORIF instrument insert

| Ref # | Description |
|----------|--|
| 542005 | Fixed angle sleeve, T8 |
| 542006 | Drill sleeve insert, T8, 60mm, ø2.0mm |
| 542008 | Variable angle drill guide, T8, ø2.0mm |
| 542009 | Compression drill guide, T8 / T15, ø2.0mm / ø2.5mm |
| 542010 | Lag screw drill guide, T8, ø2.0mm / ø2.7mm |
| 542011 | Screwdriver bit, AO, T8, 93mm |
| 542012 | Screwdriver bit, AO, T8, 180mm |
| 542016 | Depth gauge, T8, 0-80mm |
| 542015 | Screw capture sleeve, T8 |
| 542017 | Depth gauge, T15, 0-120mm |
| 542025 | Fixed angle sleeve, T15 |
| 542026 | Drill sleeve insert, T15, 60mm, ø2.5mm |
| 542028 | Variable angle drill guide, T15, ø2.5mm |
| 542030 | Lag screw drill guide, T15, ø2.5mm / ø3.5mm |
| 542031 | Screwdriver bit, AO, T15, 93mm |
| 542032 | Screwdriver bit, AO, T15, 180mm |
| 542035 | Soft tissue elevator, T15 |
| 45-80040 | Countersink for screws 2.7/3.5mm, AO fitting |
| 542103 | Threaded guide post, T15 |
| 542027 | Screw capture sleeve, T15 |
| 542098 | Delta handle, small, AO, T8 / T15 |
| 542099 | Delta handle, large, AO, T20 |
| 542000 | Drill bit, AO, Ø2.0mm x 135mm |
| 542001 | Drill bit, AO, Ø2.0mm x 175mm |



System components

Pangea Small Fragment Core Tray

Top level consisting of the instruments listed below:

Pangea Small Fragment Core Tray - ORIF instrument insert

| Ref # | Description |
|--------|-----------------------------------|
| 542002 | Drill bit, AO, Ø2.7mm x 125mm |
| 542020 | Drill bit, AO, Ø2.5mm x 135mm |
| 542021 | Drill bit, AO, Ø2.5mm x 215mm |
| 542022 | Drill bit, AO, Ø3.5mm x 135mm |
| 542003 | Tap, locking, AO, Ø2.7mm x 125mm |
| 702801 | Tap, AO, Ø2.7mm x 125mm |
| 542023 | Tap, locking, AO, Ø3.5mm x 125mm |
| 702802 | Tap, AO, Ø3.5mm x 125mm |
| 702803 | Tap, cancellous, Ø4mm x 125mm |
| 705002 | K-wire drill tip, Ø2.0mm x 234mm |
| 390157 | K-wire Ø1.25 x 150mm |
| 390164 | K-wire Ø1.6 x 150mm |
| 390192 | K-wire Ø2.0 x 150mm |
| 542036 | K-wire, olive tip, Ø1.6mm / 100mm |

Torque limiter insert instruments


| Ref # | Description |
|--------|--|
| 542014 | Torque limiting Delta handle, AO, T8 1.7Nm |
| 542034 | Torque limiting T-handle, AO, T15 4.0Nm |

System components

Pangea Small Fragment Core Tray - Instruments

Second level consisting of the instruments listed below:

This level can be removed and replaced with the optional inserts: Core tray ankle plate insert, core tray small fragment utility plate insert, or core tray Asnis III 4.0mm cannulated screw insert.




| Ref # | Description |
|--------|--------------------------------|
| 542207 | Core tray auxiliary insert |
| 542208 | Silicone mat, auxiliary insert |

System components

Pangea Small Fragment Core Tray - Instruments

Third level consisting of the instruments listed below:



| Ref # | Description |
|--------|----------------------------------|
| 542209 | Core tray small reduction insert |

Pangea small fragment core tray - Core tray small reduction insert

| Ref # | Description |
|--------|--|
| 700151 | Hook |
| 700153 | Ball spike |
| 700664 | Hohmann retractor 6mm |
| 700667 | Hohmann retractor 15mm |
| 705294 | Periosteal elevator, round edge 6mm |
| 705295 | Periosteal elevator, flat blade 13mm |
| 705297 | Straight reduction clamp, broad |
| 702932 | Repositioning forceps, L143mm (lobster claw) |
| 703938 | Bending iron |
| 705019 | Temporary plate fixator, AO |

System components

Pangea Small Fragment Core Tray - Optional ankle insert



| Ref # | Description |
|--------|------------------------------|
| 542248 | Core tray ankle plate insert |

Implants

| Ref # | Description | |
|--------|--------------------------------|--------------------------|
| 540644 | Distal lateral fibula plate, L | 2.7/3.5mm, 81mm/4 holes |
| 540664 | Distal lateral fibula plate, R | 2.7/3.5mm, 81mm/4 holes |
| 540645 | Distal lateral fibula plate, L | 2.7/3.5mm, 95mm/5 holes |
| 540665 | Distal lateral fibula plate, R | 2.7/3.5mm, 95mm/5 holes |
| 540646 | Distal lateral fibula plate, L | 2.7/3.5mm, 109mm/6 holes |
| 540666 | Distal lateral fibula plate, R | 2.7/3.5mm, 109mm/6 holes |
| 540647 | Distal lateral fibula plate, L | 2.7/3.5mm, 123mm/7 holes |
| 540667 | Distal lateral fibula plate, R | 2.7/3.5mm, 123mm/7 holes |
| 540648 | Distal lateral fibula plate, L | 2.7/3.5mm, 137mm/8 holes |
| 540668 | Distal lateral fibula plate, R | 2.7/3.5mm, 137mm/8 holes |
| 541342 | 3.5 1/3 tubular plate | 3.5mm, 26mm/2 holes |
| 541343 | 3.5 1/3 tubular plate | 3.5mm, 38mm/3 holes |
| 541344 | 3.5 1/3 tubular plate | 3.5mm, 50mm/4 holes |
| 541345 | 3.5 1/3 tubular plate | 3.5mm, 62mm/5 holes |
| 541346 | 3.5 1/3 tubular plate | 3.5mm, 74mm/6 holes |
| 541347 | 3.5 1/3 tubular plate | 3.5mm, 86mm/7 holes |
| 541348 | 3.5 1/3 tubular plate | 3.5mm, 98mm/8 holes |
| 541349 | 3.5 1/3 tubular plate | 3.5mm, 110mm/9 holes |
| 541303 | 2.7 hook plate | 2.7mm, 37mm/3 holes |
| 541305 | 2.7 hook plate | 2.7mm, 57mm/5 holes |
| 541321 | 2.7 T-plate | 2.7mm, 3Tx5 |
| 541331 | 3.5 T-plate | 3.5mm, 3Tx5 |

For the full offering of lengths please refer to the fibula and small fragment utility plate optechs.

*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

Pangea Small Fragment Core Tray - Optional Asnis III 4.0mm cannulated screw insert

| Ref # | Description |
|--------|---|
| 542245 | Core tray Asnis III 4.0mm cannulated screw insert |
| 940236 | Screw rack Asnis III 4.0mm |

Instruments

| Ref # | Description |
|--------|--|
| 702465 | Double drill guide, $\varnothing 1.4\text{mm}/2.7\text{mm}$ |
| 702446 | Cannulated drill, AO, $\varnothing 4.0\text{mm}$ |
| 702449 | Cannulated drill, AO, $\varnothing 2.7\text{mm}$ |
| 702459 | Threaded guide wire, $\varnothing 1.4 \times 150\text{mm}$ |
| 702454 | Cannulated tap, AO, $\varnothing 4.0\text{mm}$ |
| 702473 | Cannulated countersink, AO, $\varnothing 4.0\text{mm}$ screws |
| 702499 | Direct measuring gauge for wires $\varnothing 1.4/2.0 \times 150\text{mm}$ |
| 702482 | Cannulated screwdriver, AO, hex 2.5mm |
| 702485 | Solid screwdriver, AO, hex 2.5mm |
| 702489 | Holding sleeve for screwdrivers, for screwheads $\varnothing 5.0\text{mm}$ |
| 702492 | Cleaning stylet $\varnothing 1.4\text{mm}$ |
| 702496 | Extractor, AO, $\varnothing 4.0\text{mm}$ screws |

Implants

| Ref # | Description |
|--------|--------------------------------|
| 604624 | 4.0 X 24mm TI cannulated screw |
| 604626 | 4.0 X 26mm TI cannulated screw |
| 604628 | 4.0 X 28mm TI cannulated screw |
| 604630 | 4.0 X 30mm TI cannulated screw |
| 604632 | 4.0 X 32mm TI cannulated screw |
| 604634 | 4.0 X 34mm TI cannulated screw |
| 604636 | 4.0 X 36mm TI cannulated screw |
| 604638 | 4.0 X 38mm TI cannulated screw |
| 604640 | 4.0 X 40mm TI cannulated screw |
| 604642 | 4.0 X 42mm TI cannulated screw |

Pangea Small Fragment Core Tray - Optional Asnis III 4.0mm cannulated screw insert**Implants**

| Ref # | Description |
|--------|---|
| 604644 | 4.0 X 44mm TI cannulated screw |
| 604646 | 4.0 X 46mm TI cannulated screw |
| 604648 | 4.0 X 48mm TI cannulated screw |
| 604650 | 4.0 X 50mm TI cannulated screw |
| 604655 | 4.0 X 55mm TI cannulated screw |
| 604660 | 4.0 X 60mm TI cannulated screw |
| 604665 | 4.0 X 65mm TI cannulated screw |
| 604670 | 4.0 X 70mm TI cannulated screw |
| 604724 | 4.0 X 24mm TI cannulated screw, full thread |
| 604726 | 4.0 X 26mm TI cannulated screw, full thread |
| 604728 | 4.0 X 28mm TI cannulated screw, full thread |
| 604730 | 4.0 X 30mm TI cannulated screw, full thread |
| 604732 | 4.0 X 32mm TI cannulated screw, full thread |
| 604734 | 4.0 X 34mm TI cannulated screw, full thread |
| 604736 | 4.0 X 36mm TI cannulated screw, full thread |
| 604738 | 4.0 X 38mm TI cannulated screw, full thread |
| 604740 | 4.0 X 40mm TI cannulated screw, full thread |
| 604742 | 4.0 X 42mm TI cannulated screw, full thread |
| 604744 | 4.0 X 44mm TI cannulated screw, full thread |
| 604746 | 4.0 X 46mm TI cannulated screw, full thread |
| 604748 | 4.0 X 48mm TI cannulated screw, full thread |
| 604750 | 4.0 X 50mm TI cannulated screw, full thread |
| 619905 | Washer |

Pangea Small Fragment Core Tray - Optional small fragment utility plate insert

| Ref # | Description |
|--------|--------------------------------------|
| 542247 | Core tray small utility plate insert |

Implants

| Ref # | Description | |
|--------|-----------------------|----------------------|
| 541342 | 3.5 1/3 tubular plate | 3.5mm, 26mm/2 holes |
| 541343 | 3.5 1/3 tubular plate | 3.5mm, 38mm/3 holes |
| 541344 | 3.5 1/3 tubular plate | 3.5mm, 50mm/4 holes |
| 541345 | 3.5 1/3 tubular plate | 3.5mm, 62mm/5 holes |
| 541346 | 3.5 1/3 tubular plate | 3.5mm, 74mm/6 holes |
| 541347 | 3.5 1/3 tubular plate | 3.5mm, 86mm/7 holes |
| 541348 | 3.5 1/3 tubular plate | 3.5mm, 98mm/8 holes |
| 541349 | 3.5 1/3 tubular plate | 3.5mm, 110mm/9 holes |
| 541303 | 2.7 hook plate | 2.7mm, 37mm/3 holes |
| 541305 | 2.7 hook plate | 2.7mm, 57mm/5 holes |
| 541313 | 3.5 hook plate | 3.5mm, 49mm/3 holes |
| 541315 | 3.5 hook plate | 3.5mm, 79mm/5 holes |
| 541321 | 2.7 T-plate | 2.7mm, 3TX5 |
| 541322 | 2.7 T-plate | 2.7mm, 3TX10 |

| Ref # | Description | |
|--------|---------------------------|-----------------------|
| 541331 | 3.5 T-plate | 3.5mm, 3Tx5 |
| 541332 | 3.5 T-plate | 3.5mm, 3Tx10 |
| 541036 | 3.5 straight narrow plate | 3.5mm, 81mm/6 holes |
| 541037 | 3.5 straight narrow plate | 3.5mm, 104mm/7 holes |
| 541038 | 3.5 straight narrow plate | 3.5mm, 115mm/8 holes |
| 541039 | 3.5 straight narrow plate | 3.5mm, 138mm/9 holes |
| 541040 | 3.5 straight narrow plate | 3.5mm, 149mm/10 holes |
| 541042 | 3.5 straight narrow plate | 3.5mm, 183mm/12 holes |
| 541004 | 2.7 straight narrow plate | 2.7mm, 31mm/4 holes |
| 541006 | 2.7 straight narrow plate | 2.7mm, 46mm/6 holes |
| 541008 | 2.7 straight narrow plate | 2.7mm, 61mm/8 holes |
| 541010 | 2.7 straight narrow plate | 2.7mm, 76mm/10 holes |
| 541020 | 2.7 straight narrow plate | 2.7mm, 152mm/20 holes |

For the full offering of lengths please refer to the fibula and small fragment utility plate optechs.


*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



| Ref # | Description |
|--------|--------------------------------------|
| 542205 | Small fragment screw rack (with lid) |
| 542206 | Small fragment screw rack lid |

2.7mm locking screw self-tapping, T8 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 541408 | 8 | 541436 | 36 |
| 541410 | 10 | 541438 | 38 |
| 541412 | 12 | 541440 | 40 |
| 541414 | 14 | 541442 | 42 |
| 541416 | 16 | 541444 | 44 |
| 541418 | 18 | 541446 | 46 |
| 541420 | 20 | 541448 | 48 |
| 541422 | 22 | 541450 | 50 |
| 541424 | 24 | 541455 | 55 |
| 541426 | 26 | 541460 | 60 |
| 541428 | 28 | 541465 | 65 |
| 541430 | 30 | 541470 | 70 |
| 541432 | 32 | 541475S* | 75 |
| 541434 | 34 | 541480S* | 80 |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



2.7mm cortex screw self-tapping, T8 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 541708 | 8 | 541736 | 36 |
| 541710 | 10 | 541738 | 38 |
| 541712 | 12 | 541740 | 40 |
| 541714 | 14 | 541742 | 42 |
| 541716 | 16 | 541744 | 44 |
| 541718 | 18 | 541746 | 46 |
| 541720 | 20 | 541748 | 48 |
| 541722 | 22 | 541750 | 50 |
| 541724 | 24 | 541755 | 55 |
| 541726 | 26 | 541760 | 60 |
| 541728 | 28 | 541765 | 65 |
| 541730 | 30 | 541770 | 70 |
| 541732 | 32 | 541775S* | 75 |
| 541734 | 34 | 541780S* | 80 |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



3.5mm locking screw self-tapping, T15 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 541510 | 10 | 541546 | 46 |
| 541512 | 12 | 541548 | 48 |
| 541514 | 14 | 541550 | 50 |
| 541516 | 16 | 541555 | 55 |
| 541518 | 18 | 541560 | 60 |
| 541520 | 20 | 541565 | 65 |
| 541522 | 22 | 541570 | 70 |
| 541524 | 24 | 541575 | 75 |
| 541526 | 26 | 541580 | 80 |
| 541528 | 28 | 541585 | 85 |
| 541530 | 30 | 541590 | 90 |
| 541532 | 32 | 541595 | 95 |
| 541534 | 34 | 541600S* | 100 |
| 541536 | 36 | 541605S* | 105 |
| 541538 | 38 | 541610S* | 110 |
| 541540 | 40 | 541615S* | 115 |
| 541542 | 42 | 541620S* | 120 |
| 541544 | 44 | | |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



3.5mm cortex screw self-tapping, T15 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 661410 | 10 | 661446 | 46 |
| 661412 | 12 | 661448 | 48 |
| 661414 | 14 | 661450 | 50 |
| 661416 | 16 | 661455 | 55 |
| 661418 | 18 | 661460 | 60 |
| 661420 | 20 | 661465 | 65 |
| 661422 | 22 | 661470 | 70 |
| 661424 | 24 | 661475 | 75 |
| 661426 | 26 | 661480 | 80 |
| 661428 | 28 | 661485 | 85 |
| 661430 | 30 | 661490 | 90 |
| 661432 | 32 | 661495 | 95 |
| 661434 | 34 | 661500S* | 100 |
| 661436 | 36 | 661505S* | 105 |
| 661438 | 38 | 661510S* | 110 |
| 661440 | 40 | 661515S* | 115 |
| 661442 | 42 | 661520S* | 120 |
| 661444 | 44 | | |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



4.0mm cancellous screw full thread self-tapping, T15 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 607310 | 10 | 607342 | 42 |
| 607312 | 12 | 607344 | 44 |
| 607314 | 14 | 607346 | 46 |
| 607316 | 16 | 607348 | 48 |
| 607318 | 18 | 607350 | 50 |
| 607320 | 20 | 607355 | 55 |
| 607322 | 22 | 607360 | 60 |
| 607324 | 24 | 607365 | 65 |
| 607326 | 26 | 607370 | 70 |
| 607328 | 28 | 607375 | 75 |
| 607330 | 30 | 607380 | 80 |
| 607332 | 32 | 607385 | 85 |
| 607334 | 34 | 607390 | 90 |
| 607336 | 36 | 607395 | 95 |
| 607338 | 38 | 607400S* | 100 |
| 607340 | 40 | | |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Small Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



4.0mm cancellous screw partial thread self-tapping, T15 drive

| Ref # | Length (mm) | Thread length (mm) | Ref # | Length (mm) | Thread length (mm) |
|--------|-------------|--------------------|----------|-------------|--------------------|
| 607410 | 10 | 5 | 607442 | 42 | 15 |
| 607412 | 12 | 5 | 607444 | 44 | 15 |
| 607414 | 14 | 5 | 607446 | 46 | 15 |
| 607416 | 16 | 6 | 607448 | 48 | 15 |
| 607418 | 18 | 7 | 607450 | 50 | 15 |
| 607420 | 20 | 8 | 607455 | 55 | 16 |
| 607422 | 22 | 9 | 607460 | 60 | 16 |
| 607424 | 24 | 10 | 607465 | 65 | 16 |
| 607426 | 26 | 12 | 607470 | 70 | 16 |
| 607428 | 28 | 14 | 607475 | 75 | 16 |
| 607430 | 30 | 14 | 607480 | 80 | 16 |
| 607432 | 32 | 14 | 607485 | 85 | 16 |
| 607434 | 34 | 14 | 607490 | 90 | 16 |
| 607436 | 36 | 14 | 607495 | 95 | 16 |
| 607438 | 38 | 14 | 607500S* | 100 | 16 |
| 607440 | 40 | 14 | | | |



Cable plug and washers

| Ref # | Description |
|----------|------------------|
| 541400S* | 3.5mm cable plug |
| 40-30900 | Washer, T8 |
| 663001 | Washer, T15 |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Large Fragment Core Tray


Top level consisting of the instruments listed below:



| Ref # | Description |
|--------|--|
| 542201 | Aluminum lid - universal |
| 542211 | Large fragment core tray base |
| 542212 | Large fragment ORIF instrument insert |
| 542255 | Large fragment std. torque handle insert |

Pangea Large Fragment Core Tray - ORIF instrument insert

| Ref # | Description |
|--------|---|
| 542058 | Fixed angle sleeve, T20 |
| 542059 | Drill sleeve insert, T20, 60mm, ø3.2mm |
| 542060 | Drill sleeve insert, T20, 60mm, ø4.3mm |
| 542061 | K-wire sleeve insert, T20, ø2.0mm |
| 542062 | Variable angle drill guide, T20, ø3.2mm |
| 542063 | Variable angle drill guide, T20, ø4.3mm |
| 542064 | Compression drill guide, T20, ø3.2mm |
| 542065 | Lag screw drill guide, T20, ø4.5mm / ø3.2mm |
| 542066 | Screwdriver bit, AO, T20, 93mm |
| 542067 | Screwdriver bit, AO, T20, 180mm |
| 542099 | Delta handle, large, AO, T20 |
| 702812 | Countersink, AO, ø4.5/6.5mm |
| 705014 | Depth gauge, T20, 0-120mm |
| 542105 | Threaded guide post, T20 |
| 706416 | Soft tissue elevator, T20 |
| 705019 | Temporary plate fixator, AO |
| 542069 | Screw capture sleeve, T20 |
| 542050 | Drill bit, AO, Ø3.2mm x 145mm |
| 542051 | Drill bit, AO, Ø3.2mm x 215mm |
| 542052 | Drill bit, AO, Ø4.3mm x 145mm |
| 542053 | Drill bit, AO, Ø4.3 x 215mm |
| 542054 | Drill bit, AO, Ø4.5mm x 135mm |
| 542024 | Tap, locking, AO, Ø4.0mm x 145mm |
| 542057 | Tap, locking, AO, Ø5.0mm x 145mm |
| 705054 | Tap, Cancellous, Ø6mm x 180mm |
| 702808 | Tap, AO, Ø4.5 x 145mm |
| 390192 | K-wire, Ø2.0 x 150mm |
| 705002 | K-wire drill tip, Ø2.0 x 234mm |



System components

Pangea Large Fragment Core Tray

Top level consisting of the instruments listed below:

Torque limiter insert instruments


| Ref # | Description |
|--------|---------------------------------------|
| 542068 | Torque limiting T-handle, AO, T20 6Nm |

System components


Pangea Large Fragment Core Tray - Instruments

Third level consisting of the instruments listed below:

This level can be removed and replaced with the optional inserts: Core tray large fragment reduction insert or core tray Asnis III 4.0mm cannulated screw insert.



| Ref # | Description |
|--------|--------------------------------|
| 542207 | Core tray auxiliary insert |
| 542208 | Silicone mat, auxiliary insert |

Pangea Large Fragment Core Tray - Optional large fragment reduction insert

| Ref # | Description |
|--------|----------------------------------|
| 542249 | Core tray large reduction insert |


Instruments

| Ref # | Description |
|--------|--------------------------------------|
| 700151 | Hook |
| 700153 | Ball spike |
| 705294 | Periosteal elevator, round edge 6mm |
| 705295 | Periosteal elevator, flat blade 13mm |
| 702927 | Repositioning forceps, L205mm |
| 702940 | Reduction forceps with serrated jaws |
| 390084 | Reduction pin, AO, Ø5.0 x 180mm |
| 700367 | T-handle, AO quick-chuck |

System components

Pangea Large Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



| Ref # | Description |
|--------|--------------------------------------|
| 542213 | Large fragment screw rack (with lid) |
| 542214 | Large fragment screw rack lid |

4.0mm locking screw self-tapping, T20 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|--------|-------------|
| 662214 | 14 | 662242 | 42 |
| 662216 | 16 | 662244 | 44 |
| 662218 | 18 | 662246 | 46 |
| 662220 | 20 | 662248 | 48 |
| 662222 | 22 | 662250 | 50 |
| 662224 | 24 | 662255 | 55 |
| 662226 | 26 | 662260 | 60 |
| 662228 | 28 | 662265 | 65 |
| 662230 | 30 | 662270 | 70 |
| 662232 | 32 | 662275 | 75 |
| 662234 | 34 | 662280 | 80 |
| 662236 | 36 | 662285 | 85 |
| 662238 | 38 | 662290 | 90 |
| 662240 | 40 | 662295 | 95 |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Large Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



5.0mm locking screw self-tapping, T20 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 662314 | 14 | 662348 | 48 |
| 662316 | 16 | 662350 | 50 |
| 662318 | 18 | 662355 | 55 |
| 662320 | 20 | 662360 | 60 |
| 662322 | 22 | 662365 | 65 |
| 662324 | 24 | 662370 | 70 |
| 662326 | 26 | 662375 | 75 |
| 662328 | 28 | 662380 | 80 |
| 662330 | 30 | 662385 | 85 |
| 662332 | 32 | 662390 | 90 |
| 662334 | 34 | 662395 | 95 |
| 662336 | 36 | 662400S* | 100 |
| 662338 | 38 | 662405S* | 105 |
| 662340 | 40 | 662410S* | 110 |
| 662342 | 42 | 662415S* | 115 |
| 662344 | 44 | 662420S* | 120 |
| 662346 | 46 | | |



5.0mm periprosthetic locking screw self-tapping, T20 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|--------|-------------|
| 661210 | 10 | 661216 | 16 |
| 661212 | 12 | 661218 | 18 |
| 661214 | 14 | 661220 | 20 |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Large Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



4.5mm cortex screw self-tapping, T20 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 661714 | 14 | 661760 | 60 |
| 661716 | 16 | 661765 | 65 |
| 661718 | 18 | 661770 | 70 |
| 661720 | 20 | 661775 | 75 |
| 661722 | 22 | 661780 | 80 |
| 661724 | 24 | 661785 | 85 |
| 661726 | 26 | 661790 | 90 |
| 661728 | 28 | 661795 | 95 |
| 661730 | 30 | 661800S* | 100 |
| 661732 | 32 | 661805S* | 105 |
| 661734 | 34 | 661810S* | 110 |
| 661736 | 36 | 661815S* | 115 |
| 661738 | 38 | 661820S* | 120 |
| 661740 | 40 | 661825S* | 125 |
| 661742 | 42 | 661830S* | 130 |
| 661744 | 44 | 661835S* | 135 |
| 661746 | 46 | 661840S* | 140 |
| 661748 | 48 | 661845S* | 145 |
| 661750 | 50 | 661850S* | 150 |
| 661755 | 55 | | |



Plate attachments

| Ref # | Description |
|----------|------------------|
| 662202S* | 5.0mm cable plug |
| 663201 | Washer, T20 |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Large Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



6.0mm cancellous screw full thread, T20 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 608020 | 20 | 608090 | 90 |
| 608025 | 25 | 608095 | 95 |
| 608030 | 30 | 608100S* | 100 |
| 608035 | 35 | 608105S* | 105 |
| 608040 | 40 | 608110S* | 110 |
| 608045 | 45 | 608115S* | 115 |
| 608050 | 50 | 608120S* | 120 |
| 608055 | 55 | 608125S* | 125 |
| 608060 | 60 | 608130S* | 130 |
| 608065 | 65 | 608135S* | 135 |
| 608070 | 70 | 608140S* | 140 |
| 608075 | 75 | 608145S* | 145 |
| 608080 | 80 | 608150S* | 150 |
| 608085 | 85 | | |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Large Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



6.0mm cancellous screw, 16mm thread T20 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 608230 | 30 | 608295 | 95 |
| 608235 | 35 | 608300S* | 100 |
| 608240 | 40 | 608305S* | 105 |
| 608245 | 45 | 608310S* | 110 |
| 608250 | 50 | 608315S* | 115 |
| 608255 | 55 | 608320S* | 120 |
| 608260 | 60 | 608325S* | 125 |
| 608265 | 65 | 608330S* | 130 |
| 608270 | 70 | 608335S* | 135 |
| 608275 | 75 | 608340S* | 140 |
| 608280 | 80 | 608345S* | 145 |
| 608285 | 85 | 608350S* | 150 |
| 608290 | 90 | | |




*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

System components

Pangea Large Fragment Core Tray - Screws

Screw rack consisting of the implants listed below:



6.0mm cancellous screw, 32mm thread, T20 drive

| Ref # | Length (mm) | Ref # | Length (mm) |
|--------|-------------|----------|-------------|
| 608445 | 45 | 608500S* | 100 |
| 608450 | 50 | 608505S* | 105 |
| 608455 | 55 | 608510S* | 110 |
| 608460 | 60 | 608515S* | 115 |
| 608465 | 65 | 608520S* | 120 |
| 608470 | 70 | 608525S* | 125 |
| 608475 | 75 | 608530S* | 130 |
| 608480 | 80 | 608535S* | 135 |
| 608485 | 85 | 608540S* | 140 |
| 608490 | 90 | 608545S* | 145 |
| 608495 | 95 | 608550S* | 150 |



*Sterile packed only

Note: All non-sterile plates and screws may be ordered sterile by placing an "S" at the end of the ref number. Sterile option availability varies by market

1. Schmidt, W et al. "Stryker Orthopaedic Modeling and Analytics (SOMA): A Review." *Surgical Technology International*, vol. 32 (2018): 315-324.

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The instructions for use, operative techniques, cleaning instructions, patient information leaflets and other associated labeling may be requested online at ifu.stryker.com or stryker.com. If saving the instructions for use, operative techniques, cleaning instructions from the above mentioned websites, please make sure you always have the most up to date version prior to use.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, AxSOS 3, Dall-Miles, Gravity, Synchfix, and Pangea. All other trademarks are trademarks of their respective owners or holders.



Manufacturer:
Stryker GmbH
Bohnackerweg 1
2545 Selzach, Switzerland
www.stryker.com