

Sonopet iQ[®]

Ultrasonic Aspirator System

Bone cutting tips

12cm iQ Micro Claw

5500-25B-315
L 4.8in (121mm) **CL** 2.2mm **CW** 2mm **WL** 12.25mm



12cm iQ Claw

5500-25B-318
L 4.8in (121mm) **CL** 1.7mm **CW** 2.8mm **WL** 12.25mm



12cm iQ 120 Claw

5500-25B-317
L 4.8in (121mm) **CL** 2.58mm **CW** 3.15mm **WL** 12.25mm



12cm iQ Apex 360

5500-25B-312
L 4.8in (121mm) **CL** 2.58mm **CW** 3.15mm **WL** 12.25mm



20cm iQ Apex 360

5500-25B-302
L 7.9in (200mm) **CL** 2.5mm **CW** 3.15mm **WL** 11.05mm



18cm iQ Micro Claw

5500-25B-303
L 7.2in (182.8mm) **CL** 1.86mm **CW** 2mm **WL** 12mm



18cm iQ Apex Micro Claw

5500-25B-319
L 7.2in (182.8mm) **CL** 1.86mm **CW** 2mm **WL** 12mm



11cm iQ Knife

5500-25B-107
L 4.2in (106mm) **CL** 11.6mm **CW** 0.83mm **WL** 15.25mm



11cm iQ Apex Knife

5500-25B-114
L 4.2in (106mm) **CL** 12mm **CW** 0.83mm **WL** 15.6mm



10cm iQ Broad Knife

5500-25B-110
L 3.98in (101mm) **CL** 8.2mm **CW** 1.4mm **WL** 25mm



4/box

L = Length **CL** = Cutting length **CW** = Cutting width **WL** = Working length

Removal of the frangible sleeve adds 7.5mm to working lengths. On 5500-25B-114 and 5500-25B-107 there are 2 frangible sleeves meaning their working length can be increased by 7.5mm or 15mm. For 5500-25B-110 10cm iQ Broad Knife and 5500-25S-601 37cm iQ Large, there is no frangible section on the sleeve.

12cm iQ Micro

5500-25S-309

L 4.7in (120mm) **OD** 1.6mm **ID** 1.20mm **WL** 4.25mm



12cm iQ Standard

5500-25S-307

L 4.5in (115mm) **OD** 1.92mm **ID** 1.5mm **WL** 4.25mm



12cm iQ Large

5500-25S-308

L 4.4in (112.1mm) **OD** 2.49mm **ID** 2.0mm **WL** 4.25mm



20cm iQ Standard

5500-25S-301

L 8.1in (205mm) **OD** 1.92mm **ID** 1.5mm **WL** 3.8mm



12cm iQ Barracuda

5500-25S-310

L 4.5in (115mm) **OD** 1.92mm **ID** 1.5mm **WL** 4.25mm



12cm iQ Serrated Large

5500-25S-316

L 4.6in (115.7mm) **OD** 2.51mm **ID** 2.0mm **WL** 4.2mm



37cm iQ Large

5500-25S-601

L 14.6in (370.5mm) **OD** 2.52mm **ID** 1.8mm **WL** 322.8mm



4/box

L = Length **CL** = Cutting length **CW** = Cutting width **WL** = Working length

Removal of the frangible sleeve adds 7.5mm to working lengths. On 5500-25B-114 and 5500-25B-107 there are 2 frangible sleeves meaning their working length can be increased by 7.5mm or 15mm. For 5500-25B-110 10cm iQ Broad Knife and 5500-25S-601 37cm iQ Large, there is no frangible section on the sleeve.

Neurosurgical

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Apex, Barracuda, Sonopet iQ and Stryker. All other trademarks are trademarks of their respective owners or holders.

D0000011657 Rev. AC
G55/PS
Copyright © 2023 Stryker
Printed in USA

Stryker
1941 Stryker Way
Portage, MI 49002 USA
t: 269 323 7700 f: 269 323 2742
toll free: 800 253 3210

neurosurgical.stryker.com