

Stryker ENT navigation system powered by
Scopis ENT software with TGS 3.6.1

Next-level software functionality to help drive next-level outcomes

Real-time feedback

Notification and assistance for field generator alignment designed to empower users for optimal set-up.

CT auto-adjustment

Optimizes scan imports by automatically removing unwanted artifacts.

Virtual endoscopy

Designed to provide more visual information with a view from the tip of an instrument with built-in-tracking technology for complex anatomy.



3.6.1 software at a glance

Software features	3.6.1	3.4
Virtual endoscopy	✓	
Real-time feedback	✓	
Sterile user interaction: interact with the system within the sterile field	✓	
Field generator alignment for set-up	✓	
Supports the use of multiple instrument clamps to track rigid instrumentation	✓	
CT merging	✓	✓
CT auto adjustment	✓	
Procedure planning	✓	✓
Display integration	✓	✓
Augmented reality	✓	✓
Customization of screen layouts and configurations	✓	
Improved cyber security protections	✓	

Ordering information

Product	Ordering no.	Description
Electromagnetic navigation unit	8000-010-003	The electromagnetic measuring technology comes together with the Field Generator (FG). Together with the FG the unit builds the platform for the ENT Navigation System and provides the electromagnetic measuring technology
Field Generator	8000-010-004	The electromagnetic measuring technology
Scopis ENT software with TGS 3.6.1	8000-021-002	Premium target guided software with advanced planning features, virtual endoscopy and augmented reality technology


To order

Customer Service

P 800 253 3210

F 877 790 7772

inst.stryker.cs@stryker.com



Building Blocks software

Access a complimentary 1-year license of Stryker's Building Blocks software.

Stryker Building Blocks software is an educational software only



Navigation in-service video

Watch the Stryker ENT navigation system powered by Scopis software in-service video.



Instruments brochure

View our navigation instruments brochure.

ENT

See Instructions for Use - visit: <https://ent.stryker.com/physician-resources/instructions-for-use>

Please see products' Instructions For Use for a complete listing of warnings, precautions, and adverse events.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Building Blocks, Scopis, Stryker and TGS. All other trademarks are trademarks of their respective owners or holders.

ENT-NAV-BROC-971925_REV-1 | 02/2024 | Copyright © 2024 Stryker

Stryker ENT
3600 Holly Lane North,
Suite 40
Plymouth, MN 55447
O 866 620 7615
F 763 463 1599

We're on social media
Follow @StrykerENT



ent.stryker.com