

Eye Surgery Stretcher and Surgistool® Chair



Eye Surgery Stretcher

Accessible. Simple. Functional.

Stryker's Eye Surgery Stretcher is designed to offer maximum head end clearance, giving surgeons ample leg room to maneuver while performing delicate eye procedures. The work surface and mobile design combine to help with surgical throughput.



Enhanced head end clearance designed to maximize the surgeon's access for either a temporal or a superior approach to surgery.



Optional contoured or flat headpiece comfortably supports patient during procedures.

Retractable fifth wheel helps improve mobility.

Metal tooth and groove braking system helps prevent swivel and roll.

Option of two- or three-position fold-down siderails.

Multiple bumper color options allow you to customize the stretcher for identification purposes.



Navy



Gray



Red

Surgistool Chair

Customizable surgeon's stool

Surgistool Chair offers customizable seat configurations with fully adjustable backrest and armrests, which provide support to meet the surgeon's specific positioning preferences. Enhanced pedal access and a stable, yet mobile, design helps to facilitate clinical performance.



Three spindle base provides control and pedal access.



Optional padded armrests allows versatile positioning and stability.

Various seat configuration for personalized comfort.



Teardrop seat



Wedge seat

Eye Surgery Stretcher and Surgistool® Chair



Standard features: Eye Surgery Stretcher

- 3" (8 cm) Enhanced Comfort Mattress
- 8" (20 cm) Omni Surface casters with wheel covers
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side hydraulic controls
- Fold-down siderails
- Four IV receptacles
- Central-locking brakes with foot-end brake and steer pedal
- Hands-free Trendelenburg
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Rubber bumpers
- Stationary foot-end litter
- Uni-lower pedal
- Dual-articulation, enhanced clearance headpiece

Standard features: Surgistool Chair

- Adjustable contoured backrest
- Dual-density foam seat cushions
- Foot-controlled, hydraulic, three-spindle base
- Front-end locking system
- Thermoformed plastic hood
- Teardrop or wedge seat

Specifications

	Eye Surgery Stretcher	Surgistool Chair
Model number	1089	830
Overall length	90" (228.6 cm)	
Overall width	31.5" (80 cm)	
Weight capacity	500 lb. (225 kg)	300 lb. (136 kg)
Height range		
High	34" (86.4 cm) (to litter top)	28.5" (72 cm) / 30.5" (77 cm) (with riser kit)
Low	22.25" (56.5 cm) (to litter top)	21" (53 cm) / 23" (58 cm) (with riser kit)
Seat cushion style		Teardrop or wedge
Litter positioning		
Backrest	0°– 90°	
Knee gatch	0°– 30°	
Trend./Reverse Trend.	±18°	
Patient surface	26" x 87" (66 x 221 cm)	
Siderails	13" x 55" (33 x 139.5 cm)	
Caster diameter	8" (20 cm)	

Stryker reserves the right to change specifications without notice.

Optional features: Eye Surgery Stretcher

- 4" (10 cm) Ultra Comfort Mattress
- Pre/post-op shoulder pads assembly (crank backrest only)
- Air delivery system/ drape support
- Crank backrest
- Crank knee gatch
- Defibrillator tray/ foot extender
- Temporal or superior wristrest
- Footboard/chartholder
- Three-sided hydraulic controls
- Three-sided brake and steer pedals
- Integrated patient transfer/ arm board
- IV caddy
- Restraint straps
- Serving/instrument tray
- Siderail pads
- Surgery accessory rail
- Three-position siderails
- Upright oxygen bottle holder
- Oxygen bottle retaining ring
- Serving/instrument tray holder
- Headpiece - Ultra Comfort
 - 3" (8 cm), flat
 - 3" (8 cm), concave
- IV poles
 - Permanent folding two- or three-stage at foot end
 - Removable
- Extended warranty available

Optional features: Surgistool Chair

- 2" (5 cm) riser kit
- Multi-position armrest

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or Instructions for Use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Surgistool and Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

