

# ProCeed™ Hospital Bed



Designed for **durability.**  
Designed for **care.**  
Designed for **ease.**  
Designed for **you.**

Your hospital beds aren't just furniture—they play a key role on a patient's journey back to health.

Meet ProCeed—elevating experiences for nurses, staff and patients alike. Designed with input from clinicians around the world for durability and cleanability, ProCeed is designed to help you deliver efficient, quality care.



## Durability

### Cleanable

With an IPX6 rating, ProCeed's frame and electronic panel are resilient against powerful jets of water, allowing you to confidently clean your hospital bed. The support surface features welded seams to prevent water ingress.

### Durable design

With a 260 kg safe working load, the bed has been rigorously strength tested. The entire frame has been comprehensively tested for impact resistance.<sup>1</sup>



## Safety

### Full coverage folding siderails

ProCeed's full coverage siderails, with integrated status indicators, help keep patients in bed when they need to be, and fold down and tuck away when they don't.

### Low height

ProCeed's low bed height of 34 cm not only offers an ergonomic position for nurses and staff to assist patients, but also supports patient mobility and reduces their risk of injury from falls.



## Ease of use

### Radiolucent backrest

With ProCeed's radiolucent backrest, staff can take x-rays without transferring the patient to another surface, saving valuable time and energy in the process.

### Headboard options

Choose between a headboard that stays stationary when the frame tilts, or a mobile headboard that tilts when the frame tilts.

# ProCeed

## Standard features

- Moveable or stationary headboard
- Battery backup
- Manual CPR on both sides of bed
- Electric lift, electric leg rest, electric backrest
- Full coverage lockable and tuck-away siderails
- Dual siderail controls
- Removable litter covers
- Four corner brake, neutral, steer pedals
- Anti-static castor to ground frame
- Crash bumpers
- Easy adjust knee and foot elevation articulation
- Angle indicator for backrest, Trendelenburg, Reverse Trendelenburg
- Indicators for battery, articulation lockout, bed height and power
- IPX6 protection from water jets
- Ten (10) years of life expectancy
- Dual-wheel castors
- Underbed light
- Electro-static painted metal frame
- Removable headboard and footboard

## Optional features

- Brake indicator
- 5th wheel
- Nurse pendant
- Patient pendant
- X-ray translucent backrest and cassette holder
- Bed extender with linen tray

## Specifications

	ProCeed MedSurg Bed
<b>Model number</b>	8500
<b>Safe working load*</b>	260 kg
<b>Maximum patient weight</b>	215 kg
<b>Overall product size</b>	
Length	220 cm (± 1 cm)
Length (w/ bed extender - option)	251 cm (± 1 cm)
Width	99 cm (± 1 cm)
<b>Product height (w/o mattress)</b>	
Low	34 cm (+1.5/-2.5cm)
High (patient controls)	48.8 cm (± 1 cm)
High (operator controls)	76.5 cm (± 1 cm)
<b>Under product clearance</b>	15.5 cm
<b>Caster size</b>	Ø15 cm
<b>Product angle indicator</b>	0° – 15°
<b>Backrest angle indicator</b>	0° – 90°
<b>Backrest angle</b>	0° – 65°
<b>Trendelenburg/Reverse</b>	-12° – 12°
<b>Gatch angle</b>	0° – 30°
<b>Electronics</b>	100-240 VAC, 50/60 Hz

## Packages

- ProCeed Basic, includes all standard features
- ProCeed Plus, includes all standard features plus:
  - Brake indicator
  - 5th wheel
  - Nurse pendant
  - Patient pendant
  - X-ray backrest
  - Bed extender with linen tray
- ProCeed Custom\*\*
  - Fully customizable

\*Safe working load indicates the sum of the patient, mattress and accessory weight.  
 \*\* Contact your sales representative for customizable options.

### Reference: 1. Test data on file

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, ProCeed. All other trademarks are trademarks of their respective owners or holders. ProCeed : class I. The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

