

ProForm[®] Support Surface



ProForm support surface is used to assist in the prevention and treatment of pressure injuries. ProForm's foam design helps to promote immersion and envelopment, helping to enhance pressure redistribution.

Specifications

Model number	2815
Overall length	
Standard	203.2 cm (80") / 213.4 cm (84")
Overall width	88.9 cm (35")
Thickness	15.2 cm (6")
Safe working load	226.8 kg (500 lb)

Warranty

Three-year, non-prorated warranty on mattress and cover

Features

- Non-powered design
- Fire barrier*
- No unintended latex

Optional features

- Fire barrier* (standard only)

*Fire retardant inner cover meets flammability standard 16CFR1632, 16CFR1633, CGSB CAN 2-4.2 Method 27.7-M77, CAL TB129, BFDIX-11, BS EN 597-1, BS EN 597-2, UNI 9175, BS 6807;2006 Clause 9 Source 5
Stryker reserves the right to change specifications without notice.
Pro Form : Class I

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it. The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, ProForm. All other trademarks are trademarks of their respective owners or holders. The products depicted are CE marked in accordance with applicable EU Regulations and Directives. This material is not intended for distribution outside the EU and EFTA.

ProForm is a non-powered, pressure redistributing support surface through its multi-layered foam.

360 watershed prevents fluid intrusion

High density siderail foam provides a supportive crib to assist with ingress and egress

Covered foam crib increases patients surface area.