

stryker

LIFEPAK[®] 35

monitor/defibrillator



**Intuitive.
Clinically advanced.
Proven.**



When you respond to emergencies,
you need a monitor/defibrillator that's
clinically advanced and tough.^{1,2}



Built on a legacy of trust.

Ready for the future.

For decades, we've developed technology and devices that first responders like you reach for during critical events. Our hard-earned reputation has helped set the standard for innovative emergency care. It's a legacy and standard that includes exceptional attention to detail and performance that goes above and beyond.

Developing ground-breaking products takes time – but the results are worth it. Every step of our process was inspired by your commitment as we addressed your needs with the goal of surpassing your expectations. This LIFEPAK 35 may be different, both inside and out, but the focus remains the same: to help save lives.

Introducing the **LIFEPAK 35**

The LIFEPAK 35 is a clinically advanced monitor/defibrillator with proprietary tools and technology built on an intuitive,³ modern platform for advanced patient care. It's a future-ready device designed to promote confident cardiac care¹ and enable clinical excellence in today's modern healthcare environment.



Intuitive

LIFEPAK 35's large, easy-to-use touchscreen³ and advanced clinical decision support tools provide a customizable clinical experience that helps reduce cognitive burden.³



Clinically advanced

When every second matters, you need patient care that moves as fast as you do. With advanced connected capabilities providing insights and guidance, LIFEPAK 35 is reliably on your side – when time is not.



Proven

Built on a legacy of life-saving products and reimaged for today's modern healthcare professional, the LIFEPAK 35 platform carries a foundation of trust and toughness that will advance patient care.

Advancing patient care now and **into the future**

Working in partnership with customers like you, we've designed a device to help give you more control during critical events. The LIFEPAK 35 emergency response monitor/defibrillator balances advanced clinical technologies and ease-of-use^{1,2} in a device that you can rely on in challenging situations.³

SunVue™ mode allows screen to be easily viewed in different light settings.³

Built-in documentation tools can help reduce time to treatment and calculation errors.¹

cprINSIGHT® helps improve CPR performance by reducing pauses* during chest compressions.⁵

Create a custom events list of up to 80 medications and/or treatments^{2,7}

Quickly and easily streamline patient care by setting custom-timed reminders for medications and therapies.¹

Front-facing right-angle cable ports enable easy access to monitoring accessories.³



* Pre-shock pauses were drastically reduced to an average of 8 seconds vs. an average of 22 seconds with the conventional AED. CCF with cprINSIGHT was 86 percent vs. 80 percent in the conventional AED group.

LIFEPAK Next Gen: 10.4" XGA, 1024x768, 52.9in²

iPad

LIFEPAK 15:

8.4" VGA, 640x480, 33.8in²

Competition:

6.5" VGA, 640x480, 20.3in²

A big, tough touchscreen

The high-definition, 10.4" touchscreen enables a customizable display of our next-generation capabilities. The durable, chemically reinforced glass² is double glove-friendly and is designed to withstand both direct impacts and drops.² You can view up to 11 monitoring parameters at once while also sorting and renaming leads.³

Enhanced therapy options, including pediatric AED mode, help treat patients with escalating biphasic energy from 1J to 360J.^{2,6}

Streamline workflow by transmitting data for quick analysis via built-in WiFi or Bluetooth.¹

Code/Real-time clock allows a caregiver to switch between real-time and elapsed time during an event to help document and enhance patient care.¹

The intuitive user interface is easy to use and requires minimal training.³

Live view 12/15-lead with STJ Insight™ provides a graphical representation of the ECG to help diagnose myocardial injury.¹




FLEX lithium-ion dual battery system allows for nine hours of monitoring.⁴

The integrated kickstand can be easily positioned for optimal viewing.

Advanced technology promotes **clarity and confidence**

Focus on patient care with our clinically advanced therapy tools




Post-sync decision menu

Promotes patient safety by requiring confirmation to determine the next step of patient care.³



Digital calipers

Eliminate the need to print and help streamline the process of accurately calculating irregular cardiac rhythms.¹



STJ Insight

Presents a new, enhanced view of ST segment measurements to aid in assessing the location and magnitude of ischemic ST segment changes while also providing an automated alert of clinically significant changes to help clinicians diagnose myocardial injury.¹

Enhance **your response**

From the street to the emergency room, our focus is on creating seamless workflows that empower you to provide rapid, high-quality care under pressure. The LIFEPAK 35 allows you to customize events and medication lists,^{1,7} track CPR time,¹ and create on-screen reminders.^{1,3} By streamlining processes, we enable you to focus on what truly matters—saving lives and improving patient outcomes.

Includes cloud-based LIFENET® solutions⁸



Stay ready with LIFENET Asset

Remotely manage your devices, receive status alerts and view static/dynamic device location.



Respond effectively with LIFENET Care

Drive better care by transmitting, sharing and capturing patient data to enhance team communication and coordination.



Review quickly with CODE-STAT™

Automatically access or share data, review performance and identify improvement opportunities.



LIFEPAK 35 TOUGH™.

It's been rigorously tested and has an IP55 rating. So it's ready for wherever your next call takes you.

A legacy of **trust**

For six decades, we've partnered with our customers to design solutions and develop technology that help save lives. Our devices have been carried to the top of Mount Everest and in orbit on the International Space Station. You'll find more than half a million units in use today on fire rescue rigs, ambulances, and hospital crash carts worldwide. When you partner with us you'll also get industry-leading service and support. Our uncompromising commitment to quality produces devices you can count on in critical situations – so you can focus on what you do best: saving lives in the community, the field and the hospital.



For more information about the LIFEPAK 35 monitor/defibrillator, contact your local Stryker sales representative or visit [stryker/emergencycare.com](https://stryker.com/emergencycare.com).

References:

1. Data on file, Stryker OI 3337705. 2. Data on file, Stryker SHVR 3338410. 3. Data on file, Stryker HFVR 3338373. 4. Data on file, Stryker 3352191 Rev A. 5. Panchal A, Bartos J, Cabañas J, et al 2020 American Heart Association Guidelines for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care. Part 3: Adult basic and advanced life support. *Circulation*. 2020;142(suppl 2):S366-S468. 6. Data on file, Stryker UND-3331509 section 6.5. 7. Data on file, Stryker Technical Manual 3337706. 8. Data on file, Stryker LCCR 3349947 and 3347969.

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate Stryker's product offerings. A healthcare professional must always refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CODE-STAT, cprINSIGHT, LIFEENET, LIFEPAK, LIFEPAK TOUGH, STJ Insight, SunVue and Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

EC-GSNPS-BROC-1105204_REV-0_en_us © 2024 Stryker



Manufactured by:

Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
stryker.com/emergencycare