

stryker

M-1

**Roll-in ambulance
stretcher.**



M-1 Roll-in ambulance stretcher.

Reliable technology

When transporting patients, you already have enough responsibility. Instead of having to think about your materials, rely on the M-1 roll-in ambulance stretcher. Reliable technology resulting from years of experience.

Designed in close consultation with ambulance services from around the world, we have developed a stretcher that gives equal consideration to the ergonomic needs of both patients and operators.

Seen from a "safety first" perspective, this stretcher is distinguished by its outstanding lateral stability despite its 360° mobility on any terrain in the highest position.

For the optimal positioning of patients, the M-1 roll-in stretcher allows multiple settings, i.e., the shock position with knee gatch and infinite positions for the backrest. Cleaning the stretcher is also well thought out and simplifies the work involved; thanks to the removable litter and the lubrication-free system, it can be completely cleaned using a power washer.

The M-1 roll-in stretcher – the standard for first-class flexibility.



Reliable technology in every detail

The M-1 ambulance stretcher allows operation by a single person. Ergonomic lift-and-grip design with colour-coded controls provide unprecedented ease of use.

Maximum flexibility

The stretcher is easy to use and was designed with ergonomics in mind. With no need for lubrication, it can be completely cleaned using a power washer. The M-1 roll-in ambulance stretcher offers first-class flexibility.



Removable litter

Flexibility when transporting patients, litter can also be reversed to allow patients to face forward in the vehicle.

Pneumatic backrest

Infinitely adjustable.

Four corner lift handles

Retractable as needed.

Optional head-end retractable lift handles

Retractable as needed. Integrated release lever.



Optional features



Defibrillator platform



Removable O2 bottle



Knee gatch



IV Pole (right or left side)

Parallel fold-down side rails

Full functionality in tight spaces.

Multiple positioning options

Shock position, knee gatch or flat legs.

Other options

- One-hand release, fold-down side rails
- Head-end retractable lifthandles (base)



Head extension with pillow

Push/pull handle

Comfortable position, maximum maneuverability.

Foot-end controls

Patented* steer-lock system

The wheels can swivel or lock into a straight line position.

Oversized rugged wheels

Designed for performance on or off-road.





Customer support services

Technical support

Stryker technical support comprises a team of professionals available to help with your InTouch needs. Contact via phone at 1 800 667 558 or email at ssptechservices@stryker.com.

ProCare service from Stryker

Product reliability is just the beginning of the Stryker commitment to lifetime customer satisfaction. Stryker utilises advanced metrics to ensure equipment uptime and assist customers in achieving their patient and caregiver goals.

ProCare service solutions are customisable and scalable to any budget. All ProCare offerings are supported by our talented factory-trained Stryker technicians, ensuring all equipment and documentation is maintained to the highest standards.

All ProCare agreements provide:

- Stryker-authorised service representative
- Stryker-direct factory parts
- Two-hour callback response time
- Fixed service costs up front
- Increased uptime
- 24-72 hour equipment turnaround time*

*Based on the provisions of the Service Agreement and the location of the product.

Flex financial program

Our financial programs provide a range of smart alternatives designed to fit your organisation's needs. We offer flexibility beyond a cash purchase with payment structures that can be customised to meet budgetary needs and help to build long-term financial stability. Contact your account manager for more information.

Reference

1 Patent US6, 735, 794

Stryker reserves the right to change specifications without notice.

The information here is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Stryker Australia Pty Ltd

8 Herbert Street
St Leonards NSW 2065

61 2 9467 1000
www.stryker.com.au

Stryker New Zealand Ltd

515 Mt Wellington Highway
Mt Wellington, Auckland 1060

64 9 573 1890
www.stryker.com