

M-1

Roll-in ambulance Stretcher.



Key features

- One-person operation
- Roll-in loading with two-stage load release safety mechanism
- Oversized wheels with sealed castor and wheel bearings
- Colour-coded controls
- High visibility powder-coated frame
- Light-weight, rugged aluminium construction
- Shock, knee gatch or flat leg positioning
- Head-end steer lock system
- Ergonomic lift-and-grip design
- Seven height positions
- Foot-end wheel lock
- Reflective labeling
- Release lever guards
- Removable litter
- Extendable one-hand release lift handles
- Integrated bumper system
- One-hand release, infinite positioning, lift assist, pneumatic backrest
- One four-point shoulder harness restraint
- Two lap-belt restraints
- Lubrication-free maintenance
- Pressure-washable

Optional features

- Premium sealed bolster mattress that is antimicrobial
- Sealed bolster mattress
- Flat mattress
- Head extension with pillow
- Pillow (premium mattress only)
- One-hand release, fold-down siderails
- Heavy-duty, three-stage IV pole (patient right or patient left)
- DIN lock-down pins
- Head-end retractable lift handles (base)
- Foot-end stretcher retaining loop
- Foot-end stretcher retaining post
- Removable oxygen bottle holder

Warranty

- Two year

For more information contact your local sales representative

Stryker Australia Pty Ltd
8 Herbert Street
St Leonards NSW 2065

Stryker New Zealand Ltd
515 Mt Wellington Highway
Mt Wellington, Auckland 1060

61 2 9467 1000
www.stryker.com.au

64 9 573 1890
www.stryker.com

Specifications

Model number	6100	
Height¹		
Position 1	14"	(37 cm)
Position 2	24"	(62 cm)
Position 3	30"	(76 cm)
Position 4	35"	(88 cm)
Position 5	38"	(96 cm)
Position 6	39"	(100 cm)
Position 7	13"	(33 cm)
Length	77"	(197 cm)
Width	22"	(56 cm)
Weight²	74 lbs	(33 kg)
Wheels		
Diameter	6"	(15 cm)
Width	2"	(5 cm)
Maximum weight capacity	500 lbs	(228 kg)
Recommended loading height	27" - 31"	(69 - 79 cm)
Minimum operators required		
Occupied stretcher	1	
Unoccupied stretcher	1	

Litter specifications*

	Flat foot		Knee gatch	
Length	75"	(190 cm)	75"	(190 cm)
Width	22"	(56 cm)	22"	(56 cm)
Height	7"	(18 cm)	7"	(18 cm)
Weight²	43.5 lbs	(20 kg)	49 lbs	(22 kg)
Maximum weight capacity	500 lbs	(228 kg)	500 lbs	(228 kg)
Handle extension	8"	(20 cm)	8"	(20 cm)
Articulation				
Backrest	0° - 75°		0° - 75°	
Shock position	N/A		+17	
Knee gatch	N/A		30°	

¹Height measured from bottom of mattress at seat section to ground level

²Stretcher is weighed without mattress and restraints

Stryker reserves the right to change specifications without notice.

The information here is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.