

LIFENET® System

Selecting the appropriate gateway

Gateways provide a connection to your LIFEPAK® 15 or 12 monitor/defibrillator so it can connect to the LIFENET System. Stryker offers several gateway options so you can select the one which best meets your needs.



4G cellular gateway



Wi-Fi gateway



Wi-Fi/4G gateway



LIFENET PC gateway

	4G cellular gateway	Wi-Fi gateway	Wi-Fi/4G gateway	LIFENET PC gateway
Connection type	<ul style="list-style-type: none"> • 4G cellular 	<ul style="list-style-type: none"> • Wireless network 	<ul style="list-style-type: none"> • Wireless network/4G cellular 	<ul style="list-style-type: none"> • Windows® PC/tablet connection
Use environment need	<ul style="list-style-type: none"> • Transmit data anywhere you have a cellular signal • Transmission from patient's side 	<ul style="list-style-type: none"> • Use existing wireless network (ambulance, hospital, station) 	<ul style="list-style-type: none"> • Audio recording (optional) • Use existing wireless network (ambulance, hospital, station) or transmit anywhere you have a cellular signal 	<ul style="list-style-type: none"> • Use ePCR tablets that have a connection as a gateway
Benefits	<ul style="list-style-type: none"> • Always on and available when LIFEPAK device is powered on • Can transmit from patient's side; no need to be in ambulance and connected to wireless network to connect • No pairing or additional batteries needed to manage 	<ul style="list-style-type: none"> • Easy connection to existing wireless network • Supports multiple wireless networks 	<ul style="list-style-type: none"> • Automatic audio recording once LIFEPAK device turned on • Easy connection to existing wireless network • Supports multiple wireless networks • Can transmit 12-lead ECG reports and vital signs from patient's side 	<ul style="list-style-type: none"> • Use existing computer with connection (i.e., ePCR tablet) as a gateway • Use PC Gateway interface to quickly add additional information to 12-lead ECG (i.e., name, chief complaint)
Type of data transmitted	<ul style="list-style-type: none"> • 12-lead ECG reports and vital signs • Code data 	<ul style="list-style-type: none"> • 12-lead ECG reports and vital signs • Code data 	<ul style="list-style-type: none"> • 12-lead ECG reports and vital signs • Code data • Audio data 	<ul style="list-style-type: none"> • 12-lead ECG reports and vital signs • Code data
Connection to LIFEPAK device	<ul style="list-style-type: none"> • Cable (included as part of gateway) 	<ul style="list-style-type: none"> • Cable (included as part of gateway) 	<ul style="list-style-type: none"> • Cable (included as part of gateway) 	<ul style="list-style-type: none"> • Cable or Bluetooth®

	4G cellular gateway	Wi-Fi gateway	Wi-Fi/4G gateway	LIFENET PC gateway
Important details	<ul style="list-style-type: none"> • If sending large amounts of data, such as entire code files, use an unlimited data plan (a code file can be up to 2MB) • You can add gateways to your existing data plan contracts. To do so, carrier contact information is required and you must authorize the new gateways on your contract. • Determine which carrier and type of data plan before you order so that you receive the appropriate equipment • Supports LIFEPAK 15 v4+ live streaming 	<ul style="list-style-type: none"> • Requires existing wireless connection 	<p>For Cellular:</p> <ul style="list-style-type: none"> • Stryker offers annual, unlimited data plans for either Verizon® or AT&T® which allow you to get started quickly <p>For Audio:</p> <ul style="list-style-type: none"> • Requires existing wireless connection • Requires LIFEPAK 15 monitor/defibrillator with software version 3306808* or later) • Requires CODE-STAT™ 9 or higher software with LIFENET Connect for download and for playback • High-speed wireless network required for audio transmission (i.e., at station) 	<ul style="list-style-type: none"> • Software installed on Windows tablet computer • Tablet computer must have existing data connection
What to order (U.S. plans)	<p>Gateways for use with Stryker data plans</p> <ul style="list-style-type: none"> • 11996-000471 Verizon 4G Gateway and • 99428-000305 Verizon 1 year data plan or • 11996-000475 AT&T 4G Gateway and • 99428-000304 AT&T 1 year data plan <p>Gateways for use with customer data plans</p> <ul style="list-style-type: none"> • 11996-000474 Verizon 4G • 11996-000476 AT&T 4G 	<ul style="list-style-type: none"> • 21996-000109 Titan III Wireless Gateway 	<p>Gateways for use with Stryker data plans Wireless/4G/Audio (audio on/off configurable by customer)</p> <ul style="list-style-type: none"> • 11996-000479 Verizon • 11996-000484 AT&T <p>Gateways for use with customers data plans Wireless/4G/Audio (audio on/off configurable by customer)</p> <ul style="list-style-type: none"> • 11996-000480 Verizon, or AT&T 	<ul style="list-style-type: none"> • 21340-000706 LIFENET PC Gateway <p>For cable connection:</p> <ul style="list-style-type: none"> • 11230-000369 LIFEPAK Monitor to PC Cable (USB)
What to order (Canada plans)	<p>Gateway for streaming</p> <ul style="list-style-type: none"> • 11996-000503 Assembly, Gateway, 4G, Multitech, no SIM 	<ul style="list-style-type: none"> • 21996-000109 TITAN III -Wifi Gateway 	<ul style="list-style-type: none"> • 11996-000480 4G TITAN III Trio: WiFi/Cellular/Audio – no SIM 	<ul style="list-style-type: none"> • 21340-000706 LIFENET PC Gateway release 5.0

For further information, please contact Stryker at 800 442 1142 (U.S.), 800 668 8323 (Canada) or visit our website at [stryker.com](https://www.stryker.com)

Emergency Care

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CODE-STAT, LIFENET, LIFEPAK, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

M0000010338 REV AA
Copyright © 2023 Stryker

Distributed by:

Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
[stryker.com](https://www.stryker.com)

Distributed in Canada by:

Stryker Canada
2 Medicorum Place
Waterdown, Ontario
L8B 1W2
Canada
Toll free 800 668 8323