

# CLAW™ II


Polyaxial Compression Plating System  
ORTHOLOC™ 3DSi Locking Screw Technology

## The Next Generation of Compression Plates

- **Polyaxial Locking Screws** – Universal 2.7mm and 3.5mm ORTHOLOC™ 3DSi screws provide secure 30° polyaxial locking capability.
- **Expanded Plate Offerings** – Multiple styles provide versatility and specificity of plate selection.
- **Dynamic Plating System** – Surgeon-controlled compression achieved through a locked plate construct.



ORTHOLOC™ 3DSi  
2.7 & 3.5 mm Screws



### More Information:

- ❑ CLAW™ II Surgical Technique 008982
- ❑ CLAW™ II Animation FA809-1011

### Alternative Products:


- ❑ CHARLOTTE™ CLAW™
- ❑ CHARLOTTE™ Compression Staples
- ❑ DARCO™ Locked Plating System

### Complementary Products:

- ❑ CHARLOTTE™ Multi-Use Compression Screws
- ❑ DARCO™ Headed and Headless Compression Screws
- ❑ ALLO-PURE™ Bone Wedges
- ❑ BIOFOAM™ Wedges



CHARLOTTE™ Compression Screws



BIOFOAM™ Wedge System

Stryker South Pacific  
8 Herbert Street  
St. Leonards, NSW  
Australia  
Ph. 02 9467 1000

™ and ® denote Trademarks and Registered Trademarks of Wright Medical Group N.V. or its affiliates.  
©2018 Wright Medical Group N.V. or its affiliates. All Rights Reserved. 014535A\_08-Aug-2016

 **WRIGHT™**  
FOCUSED EXCELLENCE