


Phalinx™

Hammertoe Fixation






Simple and effective hammertoe fixation


Hammertoe Fixation

System Advantages

- Versatile system addressing hammertoe correction and MTP joint instability
- Titanium implants eliminate need for freezer storage
- Implants available in straight cannulated and 10° angled
- Color-coded implants and instruments increase system efficiency



Implant Description/Size	Length (mm)			Diameter (mm)		
	A Overall	B Proximal	C Distal	D Proximal	E Distal	F Thickness
Extra Small Cannulated/Angled	11.2	6	4.2	2.5	3.2	1.8
Small Cannulated/Angled	13.0	7	5.5	2.8	4.0	1.9
Medium Cannulated/Angled	14.1	7	6.0	3.0	4.5	2.0
Large Cannulated/Angled	15.2	8	6.2	3.3	5.0	2.0



Foot & Ankle

Visit footankle.stryker.com or email us at footandankle@stryker.com

The information contained in this document is intended for healthcare professionals only. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Phalinx, Stryker. All other trademarks are trademarks of their respective owners or holders.