

Sage

PrimoFit+™

External Urine Management for the Male Anatomy



Easily track usage

Convenient time/date label encourages compliance to protocol


Stays in place

Adhesive pad and base adhesive keeps device in place and allows for easy application and removal

Customizable fit

Eliminates the need for sizing charts and adjusts to fit most patients' anatomy, from adolescent to adult

Back



Easy to see

Blue liners are easy to see after removal to reduce the risk of misplacement

Reduces leakage

Semi-closed top reduces the risk of leakage while maintaining airflow²

Helps keep skin dry

Ultra-soft fabric diverts urine away from the skin

1-step connection

Tiered suction tubing adapter allows for easy connection to hospital suction tubing



Can be left in place for up to **24 hours**

A breakthrough in male urine management

PrimoFit+ provides an easily customizable fit that stays in place and uses continuous suction to effectively manage urinary incontinence while keeping patients comfortable. PrimoFit+ serves as an alternative to external collection devices or, when appropriate, indwelling catheters.

PrimoFit+ diverts urine away from the skin. Urine is then diverted through the device and suctioned into a canister, helping keep the skin dry and addressing a risk factor of incontinence-associated dermatitis (IAD).¹

Approximately

75%

of hospital-acquired UTIs in the US are **associated with an indwelling urinary catheter.**³



(1) PrimoFit+, adolescent - adult size/package

40 devices/case
Reorder #5405-X

References: 1. Gray M, et al., Incontinence-associated Dermatitis: A Consensus, Journal of Wound Ostomy & Continence Nursing, 2007;(34(1)):45-54. 2. Data on file. Sage Products LLC. P23-0041R and P22-0040R. 3. Centers for Disease Control and Prevention, Catheter-associated Urinary Tract Infections (CAUTIs). Available at: https://www.cdc.gov/hai/ca_uti/uti.html. Accessed August 24, 2017.

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: PrimoFit+, Sage, and Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

PrimoFit+ is a Class I medical device CE.

This document is not intended for distribution outside EU and EFTA.

SAGE-UM-BROC-865240_REV-0_en_us © 2024 Stryker



Sage Products LLC
3909 Three Oaks Road
Cary, IL 60013
800 323 2220

stryker.com