

Sage
Prevalon

TAP LC™

Lift Compatible Patient Repositioning System



3-in-1 mobility

1 Boost

Help address
patient migration

2 Turn

Help address risk factors
of pressure injuries

3 Lift

Safely transfer patient
to the bedside chair

Promoting **safety and mobility**

stryker

We understand the importance of early mobility. To achieve your mobility protocols and comply with guidelines, you need a readily available solution that stays under your patient and fits seamlessly into your workflow. That's why we combined our innovative turning and positioning technology with lift compatibility to promote safety and mobility.

30° Body and Anchor Wedge System

- Helps sustain a 30° side-lying position and reduces pressure by offloading sacrum
- Minimize manual activity as wedges slide into a turn, reducing the need for a log roll
- Minimizes patient migration down the bed
- Reduces the need for boosting, minimizing shear and friction

M² Microclimate Body Pad

Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through

Hi-tech fabric

- Promotes a healthy microclimate for the skin
- Easy-to-clean surface

Easy-roll handles


Encourage proper posture for healthcare worker

Black positioning handles

Encourage proper ergonomic position for the healthcare worker when preparing to reposition patient




Intuitive strap system



Enables multiple lifting positions

TAP LC Lift Compatible Patient Repositioning System

41 in x 65 in
Weight limit: 600 lb / 272 kg



- (1) Lift Compatible Glide Sheet
- (1) 30° Body and Anchor Wedge System
- (1) M² Microclimate Body Pad

5 systems/case
Reorder #7281

- (1) Lift Compatible Glide Sheet
- (1) M² Microclimate Body Pad

10 Glide Sheets/case
Reorder #7282



Scan for more information

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Prevalon, TAP LC, Sage and Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

SAGE-PVN-BROC-1020799_REV-0_en_us © 2024 Stryker

Sage Products LLC | 3909 Three Oaks Road | Cary, IL 60013 | stryker.com | 800 323 2220