

# XBraid suture kits quick guide

XBraid TT and NiceLoop offerings together create comprehensive suture kits to address your orthopedic surgical needs.



**XBraid TT** is a smooth, non-absorbable low profile, 100% ultra-high molecular weight polyethylene suture tape designed to lie flat against the tissue and increase tendon-to-bone contact in soft tissue repairs.<sup>1</sup>



**NiceLoop** is a non-absorbable braided suture for reliable fixation. It comes pre-looped, sterile and on a needle with a self-stabilized NiceKnot that prevents slippage.<sup>2</sup>

## Kit contents:



### Small open soft tissue repair suture kit

1.8mm XBraid TT suture tape

XBraid S needle CP-1 Cobra Black

1.4mm XBraid TT suture tape

XBraid S needle CCS Cobra Black

### Large open soft tissue repair suture kit

2.2mm XBraid TT suture tape

XBraid S needle CP-1 Cobra Black

1.4mm XBraid TT suture tape

XBraid S needle CCS Cobra Black

## Suture kit offerings

Product number	Description
3910900145	<p><b>Small open soft tissue repair suture kit</b></p> <p>1.4mm XBraid TT, blue premium, suture tape, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 40" suture</p> <p>1.4mm XBraid TT, white/blue premium, suture tape, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 40" suture</p> <p>1.8mm XBraid TT, black/white premium, suture tape, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 40" suture</p> <p>1.8mm XBraid TT, blue/black premium, suture tape, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 40" suture</p> <p>1.8mm XBraid TT, blue premium, suture tape, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 40" suture</p> <p>1.8mm XBraid TT, white/black premium, suture tape, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 40" suture</p> <p>XBraid S, needle, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 10mm loop</p> <p>XBraid S, needle, CCS round body ribbed (conventional cutting, 1/2c, 48mm, ribbed) Cobra Black, 10mm loop</p>
3910900146	<p><b>Large open soft tissue repair suture kit</b></p> <p>1.4mm XBraid TT, blue premium, suture tape, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 40" suture</p> <p>1.4mm XBraid TT, white/blue premium, suture tape, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 40" suture</p> <p>2.2mm XBraid TT, black premium, suture tape, CCS-1 round body ribbed (conventional cutting, 1/2c, 60mm) Cobra Black, 24" suture loop</p> <p>2.2mm XBraid TT, blue/black premium, suture tape, CCS-1 round body ribbed (conventional cutting, 1/2c, 60mm) Cobra Black, 24" suture loop</p> <p>2.2mm XBraid TT, blue premium, suture tape, CCS-1 round body ribbed (conventional cutting, 1/2c, 60mm) Cobra Black, 24" suture loop</p> <p>2.2mm XBraid TT, white/blue premium, suture tape, CCS-1 round body ribbed (conventional cutting, 1/2c, 60mm) Cobra Black, 24" suture loop</p> <p>2.2mm XBraid TT, white/black premium, suture tape, CCS-1 round body ribbed (conventional cutting, 1/2c, 60mm) Cobra Black, 24" suture loop</p> <p>2.2mm XBraid TT, solid black premium, suture tape, CCS-1 round body ribbed (conventional cutting, 1/2c, 60mm) Cobra Black, 24" suture loop</p> <p>XBraid S, needle, CP-1 (reverse cutting 1/2c 36mm) Cobra Black, 10mm loop</p> <p>XBraid S, needle, CCS-1 round body ribbed (conventional cutting, 1/2c, 60mm) Cobra Black, 10mm loop</p>

### References

1. Stryker DHFD12758
2. SSA-8484 Rev C DCR 180107

## Sports Medicine

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/ or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets.

Please contact your representative if you have questions about the availability of Stryker's products in your area. Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: XBraid, NiceLoop and Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

Stryker  
Sports Medicine  
5670 Greenwood Plaza Blvd.  
Ste. 200  
Greenwood Village, CO 80111

t: 866 596 2022  
www.sportsmedicine.stryker.com